

3050 Spruce Street Saint Louis, Missouri 63103 USA Telephone 800-325-5832 • (314) 771-5765 Fax (314) 286-7828 email: techserv@sial.com sigma-aldrich.com

AURAMINE O

Product Number 85,653-3 Store at Room Temperature Replacement for Product Number A 1396

CAS #: 2465-27-2

Synonyms: 1 C.I. 41000; Basic Yellow 2; Canary yellow; Pyoktaninum aureum; 4,4'-(imidocarbonyl)bis(N,Ndimethylaniline)monohydrochloride

Product Description

$$\begin{array}{c|c} H_3C & & \text{NH} \\ \hline \\ H_3C & & \text{CH}_3 \end{array} \bullet \text{HC}$$

Appearance: Yellow powder with orange-brown cast

Molecular formula: C₁₇H₂₁N₃·HCl

Molecular weight: 303.8

spectral data: λ_{exc} = 438 nm (water); λ_{em} = 505 nm. The absorbance spectrum has two peaks:

432 nm and 370 nm in water.3

Auramine O usually contains dextrin or other diluents that are added to standardize dye content. A9655 has been approved by the BSC (Biological Stain Commission).

Used alone, Auramine O colors many materials a golden yellow. In combination with other dyes, it creates or enhances blues or reds. The quantum yield of this fluorescent dye depends on solvent viscosity; in 95% ethanol or 60% sucrose it is 87 and 10 times higher, respectively, than in water.2

Several procedures using auramine O as a fluorochrome are reported:5

· used with rhodamine B, causing acid-fast organisms to fluoresce.

ProductInformation

- used with acridine orange, showing acid-fast organisms and fungi, mycobacteria. (Auramine O is reportedly specific for mycolic acid)
- · used as a Schiff reagent, causing DNA in nuclei to fluoresce yellow-green.
- used with safranin O chloride to stain spores in bacterial cultures.

A solution of Auramine O in 0.1 M sodium phosphate buffer has no detectable fluorescence. However, the dye forms a complex with horse liver alcohol dehydrogenase; in the presence of very low concentrations of the enzyme, an intense fluorescence at 523 nm is observed (excitation at approx. 440-450 nm). It did not form similar complexes with sixteen other proteins tested.6

This product has been synthesized by heating 4,4'bis(dimethylamino)benzophenone with ammonium chloride. It readily converts to the precursor, with the release of ammonia, when it is subjected to acid or water hydrolysis (~60 °C).3

Preparation Instructions

Sigma tests solubility in ethanol with heating. At 1 mg/mL, a clear bright yellow solution is obtained. At 10 mg/mL, it does not dissolve entirely in ethanol, even with heating and sonication.⁴ Variation in dye content will affect solubility in water or ethanol. This product is reported to be soluble in water at 10 mg/mL.

Solutions are expected to be stable several months at 2-8 °C (light-protected)^{5a}, but fresh dye solutions are recommended, since stored solutions show a rapid quenching on exposure to ultraviolet light.5b

Storage/Stability

When stored at room temperature, this product has a shelf-life of four years.

References

1. H.J. Conn's Biological Stains, 9th ed. (Williams and Wilkins, 1977), R.D. Lillie, ed., p. 238-239.

- 2. Methods in Enzymology, 32, 230 (1974), "Fluorescent Probes".
- 3. Green, Floyd J., *Sigma-Aldrich Handbook of Stains Dyes and Indicators*, p. 103-104.
- 4. Sigma quality control.

- 5. Staining Procedures, 4th ed. (Williams and Wilkins, 1981) G. Clark, ed., a) p. 382-384; b) p. 391; 51-52, 59-61.
- 6. Conrad, R.H. et al., Biochem., 9, 1540-1546 (1970).

feb 04/03