Product Information

Hyaluronidase from bovine testes

Storage Temperature -20 °C

CAS RN 37326-33-3

EC 3.2.1.35

Synonyms: Hyaluronoglucosaminidase, Hyaluronate 4-glycanohydrolase

Product Description

Molecular mass: Molecular mass:¹ 60 kDa

Extinction Coefficient:² $E^{1\%} = 8.0$ (280 nm)

Hyaluronidase from bovine testes is a tetramer of 4 equal subunits with a molecular mass of 14 kDa each.³ This glycoprotein enzyme contains 5% mannose and 2.2% glucosamine.¹ It randomly hydrolyzes 1,4-linkages between N-acetyl- β D glucosamine and D-glucuronate residues in hyaluronate. It also hydrolyzes 1,4- β -D-glycosidic linkages between N-acetyl-galactosamine or N-acetyl-galactosamine sulfate and glucuronic acid in chondroitin, chondroitin 4- and 6-sulfates, and dermatan.⁴ The enzyme's pH optimum is 4.5 to 6.0.⁵

Hyaluronidase is often used in conjunction with collagenase to dissociate the extracellular matrix between cells of animal tissue, in order to release viable cells for use in tissue culture. It may also be used to clarify synovial fluids in order to make cell counts possible.

The following compounds inhibit hyaluronidase: Fe⁺², Fe⁺³, Zn⁺², Cu⁺², heparin, aurintricarboxylic acid, sulfated, nitrated, or acetylated hyaluronic acids, sulfated cellulose esters, sulfated chitin esters, sulfated carboxycellulose, sulfated xylan, bile salts, sulfated steroids, hexylresorcinol, *o*- and *p*- quinones, and sulfated aliphatic alcohols.⁵

Specific Activity: 400-1,000 units/mg solid

Unit definition: One unit will degrade 0.75 μ g of the polysaccharide hyaluronic acid per minute at pH 5.35 at 37 °C (as measured by turbidimetric absorbance (λ : 600 nm) when complexed with BSA after 45 minutes).

Precautions and Disclaimer

For R&D use only. Not for drug, household, or other uses. Please consult the Safety Data Sheet for information regarding hazards and safe handling practices.

Preparation Instructions

This enzyme is soluble in 0.02 M phosphate buffer, pH 7, with 77 mM sodium chloride and 0.01% BSA at 1 mg/mL, yielding a clear solution.

Storage/Stability

It is recommended that fresh solutions be prepared before use. Stability testing of frozen solutions of this product has not been performed in our laboratories. One published reference cites storage of 1.0 mg/mL stock solutions of hyaluronidase in PBS at -20 °C.⁷



References

- Borders, C.L., Jr., and Raftery, M.A., J. Biol. Chem., 243(13), 3756-3762 (1968).
- Worthington Enzyme Manual (Freehold, NJ: 1993), pp. 224.
- Khorlin, A.Y. et al., FEBS Lett., 31(1), 107-110 (1973).
- 4. Enzyme Nomenclature, Academic Press (San Diego, CA), p. 350 (1992).
- The Enzyme Handbook, Schomburg, D., and Salzmann, M., Springer-Verlag (Berlin/Heidelberg) Vol. 4, EC 3.2.1.35, p. 2 (1991).
- De Salegui, M. et al., Arch. Biochem. Biophys., 121(3), 548-554 (1967).
- Brem, G., and Nowshari, M.A., in Methods in Molecular Biology, Volume 254: Germ Cell Protocols (Volume 2; H. Schatten, ed.). Humana Press (Totowa, NJ), p. 218 (2004).

Notice

We provide information and advice to our customers on application technologies and regulatory matters to the best of our knowledge and ability, but without obligation or liability. Existing laws and regulations are to be observed in all cases by our customers. This also applies in respect to any rights of third parties. Our information and advice do not relieve our customers of their own responsibility for checking the suitability of our products for the envisaged purpose.

The information in this document is subject to change without notice and should not be construed as a commitment by the manufacturing or selling entity, or an affiliate. We assume no responsibility for any errors that may appear in this document.

Technical Assistance

Visit the tech service page at <u>SigmaAldrich.com/techservice</u>.

Terms and Conditions of Sale

Warranty, use restrictions, and other conditions of sale may be found at <u>SigmaAldrich.com/terms</u>.

Contact Information.

For the location of the office nearest you, go to <u>SigmaAldrich.com/offices</u>.

The life science business of Merck operates as MilliporeSigma in the U.S. and Canada.

Merck and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

@ 2023 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved. H3506pis Rev 08/23

