For life science research only. Not for use in diagnostic procedures.



# **Chromozym PK**Benzoyl-Pro-Phe-Arg-4-nitranilide acetate

Version: 17
Content Version: May 2021

Cat. No. 10 378 445 001 20 mg

Store the product at +15 to +25°C.

1.	General Information	3
1.1.	Contents	3
1.2.	Storage and Stability	
	Storage Conditions (Product)	
1.3.	Additional Equipment and Reagent required	3
1.4.	Application	3
2.	How to Use this Product	4
2.1.	Before you Begin	
	General Considerations	
	Calculation of plasma kallikrein activity Working Solution	
2.2.	Protocols	
Z.Z.	Preparation of sample material	
	Plasma kallikrein assay	
2.3.	Parameters	6
	Chemical Formula	
	Chemical Name	
	ContaminantsMolecular Weight	
	Purity	
	Working Concentration	
3.	Additional Information on this Product	7
3.1.	Test Principle	7
	Reaction principle	7
4.	Supplementary Information	8
4.1.	Conventions	8
4.2.	Changes to previous version	8
4.3.	Ordering Information	8
4.4.	Trademarks	9
4.5.	License Disclaimer	9
4.6.	Regulatory Disclaimer	9
4.7.	Safety Data Sheet	9
4.8.	Contact and Support	9

## 1. General Information

#### 1.1. Contents

Vial / Bottle	Label	Function / Description	Content
1	Chromozym PK	Powder	1 vial,
	$C_{33}H_{39}N_8O_6COOCH_3 - M_7 702.9$		20 mg

# 1.2. Storage and Stability

## **Storage Conditions (Product)**

When stored at +15 to +25°C, the product is stable through the expiry date printed on the label.

Vial / Bottle	Label	Storage
1	Chromozym PK	Store at +15 to +25°C.

# 1.3. Additional Equipment and Reagent required

#### For preparation of working solutions

- Dextran sulfate
- 1 M HCI
- Tris\* [2-Amino-2-(hydroxymethyl)-1.3-propandiol]
- Imidazole
- NaCl
- Double-distilled water

#### For plasma kallikrein assay

- See section, Working Solution for information on preparing solutions.
- Plastic cuvettes
- Cuvette rack with thermostat
- 0.1 1 M sodium citrate solution

# 1.4. Application

Chromozym PK is used as substrate for the determination of serine proteases, especially of plasma kallikrein, such as in citrated plasma.

# 2. How to Use this Product

# 2.1. Before you Begin

## **General Considerations**

## Calculation of plasma kallikrein activity

U/ml citrated plasma = 
$$\frac{V}{v \times \epsilon \times d} \times F \times \Delta A$$

Abbreviation	Description		
F	Dilution factor = 2		
V	Assay volume (2.022 ml)		
V	Sample volume (0.02 ml)		
ε405 nm	10.4 [mmol <sup>-1</sup> × l × cm <sup>-1</sup> ]  i Absorbance coefficient for 4-nitraniline.		
d	Light path in the cuvette (1 cm)		
$\Delta A/\Delta t = \Delta A/min$	Change in absorbance/minute		
U/ml citrated sample =	× 2 × <u>ΔΑ</u>		

U/ml citrated sample = 
$$\frac{2.02}{0.02 \times 10.4 \times 1} \times 2 \times \underline{\Delta A}$$
 min

# **Working Solution**

Preparation of the solutions for approximately 50 determinations.

Solution No.	Solution	Composition/Preparation	Storage and Stability
1	Dextran sulfate solution	<ul> <li>25 mg/ml dextran sulfate</li> <li>Dissolve 2.5 g dextran sulfate 500, sodium salt, in double-distilled water and fill up to 100 ml.</li> </ul>	Store for at least 4 weeks at +2 to +8°C.  •• Keep at 0°C before use in the assay.
2	HCI	1.0 M HCl solution	Store for several months at +15 to +25°C.
3	Tris/ imidazole buffer	15 mM, pH 7.9; ion strength 0.15  a) 0.25 M Tris, 0.25 M imidazole: Dissolve 3.03 g Tris and 1.70 g imidazole in 50 ml 1 M HCl (Solution 2) and fill up to 100 ml with double-distilled water. b) 0.33 M Tris, 0.33 M imidazole, 0.33 M NaCl: Dissolve 4.04 g Tris, 2.27 g imidazole, and 1.95 g NaCl in approximately 75 ml double-distilled water; fill up to 100 ml with double-distilled water. c) NaCl solution (2 M NaCl): Dissolve 11.63 g NaCl in approximately 75 ml double-distilled water and fill up to 100 ml.	Store for at least 4 weeks at +2 to +8°C.  Avoid contamination with germs.
4	Incubation buffer	<ul> <li>Mix Solutions 3a and 3b to reach a resulting buffer of pH 7.9.</li> <li>Mix 5 ml of this buffer with 5 ml of Solution 3c.</li> <li>Dilute this buffer concentrate 1 + 9 immediately before use.</li> </ul>	⚠ Dilute freshly and warm up to +37°C before use in the assay.

5 Chromozym PK solution PK solution PK solution PK solution ■ 1.3 mM Store for at least 4 weeks at +2 to +8°C. 

■ Warm up to +37°C before use in the assay.

#### 2.2. Protocols

#### **Preparation of sample material**

- Remove blood from an uncongested vein.
- 2 Immediately after removal, mix 9 volume parts of blood with 1 volume part 0.1 1 M sodium citrate solution.
- 3 Centrifuge the sample at approximately  $2,000 \times g$  (approximately 3,000 rpm with conventional lab centrifuges) for 10 minutes within 2 hours after removal.
- *i* Citrated plasma is stable for up to 4 hours stored at +2 to +8 $^{\circ}$ C, or for 6 months stored at -15 to  $-25^{\circ}$ C.

#### Plasma kallikrein assay

Apply the following photospectrometer parameters:

Wavelength	Hg 405 nm
Cuvettes	plastic  i Use a cuvette rack with thermostat.
Light path	1 cm
Temperature	+37°C
Assay volume	2.02 ml
Measurement	against air

- ▲ Equilibrate Solutions 4 and 5 to +37°C before using in the assay. Keep Solution 1 at 0°C.
- A reagent blank is not necessary.
- 3 See section, Working Solution for information on preparing solutions.
- For the Activation, pipette into a plastic cuvette:
  - 0.1 ml plasma
  - 0.1 ml Dextran sulfate solution (Solution 1).
- Mix and incubate for 7 minutes at 0°C (= activated sample).
- 3 For the Kallikrein catalyzed reaction, pipette into a plastic cuvette:
  - 1.52 ml Incubation buffer (Solution 4).
  - 0.52 ml Chromozym PK solution (Solution 5).
  - 0.02 ml activated sample (from Step 1).
- Mix immediately.
  - Read absorption within 30 seconds and simultaneously start the stop watch.
  - Repeat reading after 1, 2, 3, and 4 minutes or follow the reaction by a recorder.
  - For the calculation of the plasma kallikrein activity, see section, General Considerations.

## 2.3. Parameters

## **Chemical Formula**

 $C_{33}H_{39}N_8O_6COOCH_3$ 

#### **Chemical Name**

Benzoyl-prolyl-phenalanyl-arginine-4-nitranilide acetate

## **Contaminants**

<0.5% free 4-nitraniline.

# **Molecular Weight**

702.9 Da

# **Purity**

90% Benzoyl-Pro-Phe-Arg-4-nitranilide acetate (enzymatic).

# **Working Concentration**

Approximately 0.5 mM.

# 3. Additional Information on this Product

# 3.1. Test Principle

## **Reaction principle**

Chromozym PK is cleaved by plasma kallikrein into a residual peptide and free 4-nitraniline which is measured at 405 nm

1 The absorbance difference per minute is used for the determination of the kallikrein activity in U/ml:

Starting material	reacts with	to
Factor XII	dextran sulfate	Factor XIIa
Plasma prokallikrein	factor XIIa	Plasma-kallikrein
Bz-Pro-Phe-Arg-4-NA + H <sub>2</sub> O	plasma-kallikrein	B-Pro-Phe-Arg + 4-nitraniline

# 4. Supplementary Information

## 4.1. Conventions

To make information consistent and easier to read, the following text conventions and symbols are used in this document to highlight important information:

Text convention and symbols		
1 Information Note: Additional information about the current topic or procedure.		
⚠ Important Note: Information critical to the success of the current procedure or use of the product.		
1) 2) 3) etc.	Stages in a process that usually occur in the order listed.	
1 2 3 etc.	Steps in a procedure that must be performed in the order listed.	
* (Asterisk)	The Asterisk denotes a product available from Roche Diagnostics.	

# 4.2. Changes to previous version

Layout changes. Editorial changes.

# 4.3. Ordering Information

Product	Pack Size	Cat. No.	
Reagents, kits	Reagents, kits		
Tris base	1 kg, <i>Not available in US</i>	10 708 976 001	
	1 kg	03 118 142 001	
	5 kg	11 814 273 001	

Version: 17

#### 4.4. Trademarks

All product names and trademarks are the property of their respective owners.

#### 4.5. License Disclaimer

For patent license limitations for individual products please refer to: **List of biochemical reagent products**.

# 4.6. Regulatory Disclaimer

For life science research only. Not for use in diagnostic procedures.

# 4.7. Safety Data Sheet

Please follow the instructions in the Safety Data Sheet (SDS).

# 4.8. Contact and Support

To ask questions, solve problems, suggest enhancements or report new applications, please visit our **Online Technical Support Site**.

To call, write, fax, or email us, visit **sigma-aldrich.com**, and select your home country. Country-specific contact information will be displayed.