

3050 Spruce Street
Saint Louis, Missouri 63103 USA
Telephone 800-325-5832 • (314) 771-5765
Fax (314) 286-7828
email: techserv@sial.com
sigma-aldrich.com

ProductInformation

1-Thioglycerol

Product Number **M 1753** Storage Temperature 2-8 °C

Product Description

Molecular Formula: C₃H₈O₂S Molecular Weight: 108.2 CAS Number: 96-27-5 Density: 1.245 g/ml¹

Molarity: 11.5 M (calculated based on density and

molecular weight)

Boiling Point: 100-101 °C¹

Synonyms: α -monothioglycerol, thioglycerin,

3-mercapto-1,2-propanediol²

1-Thioglycerol is a reagent that is used in analytical chemistry, materials science, and cell culture research. In mass spectrometry (MS), it is notably used as a matrix substrate in FAB (fast atom bombardment) MS methods.^{3,4} Compounds that have been investigated by FAB-MS using 1-thioglycerol as the FAB-MS matrix include lipids, gangliosides, benzo[a]pyrenes, and iron-siderophore complexes.^{5,6,7,8,9}

In materials and biomaterials science, 1-thioglycerol has been used in the preparation of thiol modified gold surfaces for human immunoglobulin binding and of water soluble luminescent CdS quantum dots. 10,11,12 Cell culture studies of embryonic cortical and hippocampal neurons, mouse bone marrow mast cell lines, and human B cell lines have utilized 1-thioglycerol as a component of the culture medium to stimulate proliferation. 13,14,15,16

Precautions and Disclaimer

For Laboratory Use Only. Not for drug, household or other uses.

Preparation Instructions

This product is miscible in ethanol (1 ml/ml, 50%, v/v), yielding a clear, colorless solution. It is also miscible in water (0.1 M).

References

- 1. CRC Handbook of Chemistry and Physics, 82nd ed., Lide, D. R., ed., CRC Press (Boca Raton, FL: 2001), p. **3**-274.
- 2. The Merck Index, 12th ed., Entry# 9471.
- Seifert, W. E., Jr., and Caprioli, R. M., Fast atom bombardment mass spectrometry. Methods Enzymol., 270, 453-486 (1996).
- 4. Dass, C., The role of a liquid matrix in controlling FAB-induced fragmentation. J. Mass Spectrom., **31(1)**, 77-82 (1996).
- Lenz, C., et al., Improved ternary matrices for the analysis of ferri-pyoverdins by fast atom bombardment mass spectrometry. J. Mass Spectrom., 33(10), 984-987 (1998).
- Egge, H., et al., Analysis of gangliosides using fast atom bombardment mass spectrometry. Chem. Phys. Lipids, 37(2), 127-141 (1985).
- 7. Siegel, M. M., et al., Fast atom bombardment mass spectrometry of cisplatin analogs. Biomed. Environ. Mass Spectrom., **13(1)**, 25-32 (1986).
- 8. Riviere, M., et al., Particular matrix for fast atom bombardment mass spectrometric analysis of phenolic glycolipid antigens isolated from pathogen mycobacteria. Biomed. Environ. Mass Spectrom., **16(1-12)**, 275-278 (1988).
- 9. Greaves, J., and Bieri, R. H., Fast-atom bombardment mass spectrometry of conjugated benzo[a]pyrene metabolites. Rapid Commun. Mass Spectrom., **3(9)**, 286-288 (1989).
- Tengvall, P., et al., Complement activation by 3mercapto-1,2-propanediol immobilized on gold surfaces. Biomaterials, 17(10), 1001-1007 (1996).
- Tang, L., et al., Complement activation and inflammation triggered by model biomaterial surfaces. J. Biomed. Mater. Res., 41(2), 333-340 (1998).
- 12. Chen, Y., and Rosenzweig, Z., Luminescent CdS quantum dots as selective ion probes. Anal. Chem., **74(19)**, 5132-5138 (2002).

- Hultner, L., et al., Thiol-sensitive mast cell lines derived from mouse bone marrow respond to a mast cell growth-enhancing activity different from both IL-3 and IL-4. J. Immunol., 142(10), 3440-3446 (1989).
- Grill, R. J., Jr., and Pixley, S. K.,
 2-Mercaptoethanol is a survival factor for olfactory, cortical and hippocampal neurons in short-term dissociated cell culture. Brain Res., 613(1), 168-172 (1993).
- Brielmeier, M., et al., Improving stable transfection efficiency: antioxidants dramatically improve the outgrowth of clones under dominant marker selection. Nucleic Acids Res., 26(9), 2082-2085 (1998).
- Iciek, M., and Wlodek, L., Biosynthesis and Biological Properties of Compounds Containing Highly Reactive, Reduced Sulfane Sulfur. Pol. J. Pharmacol., 53, 215-225 (2001).

GCY/NSB 7/03