

3dGRO® Gastrointestinal Organoids

SCC310-SCC370

FOR RESEARCH USE ONLY

Not for use in diagnostic procedures. Not for human or animal consumption.

Product Overview

Patient derived organoids (PDOs) are novel *in vitro* 3D cell models that preserve original tissue physiology and molecular pathology, thus representing a clinically relevant alternative to traditional 2D cell lines and an effective tool to refine and reduce animal models. PDOs can be derived from adult patient biopsies or resected tissues containing native LGR5+ stem cell populations and cultured within an ECM-rich substrate using specialized organoid media adult tissue derived organoids are phenotypically and genetically stable in long term culture, presenting more mature phenotypes compared to iPSC-derived organoids. Importantly, PDOs have shown to be able to predict patient clinical responses to chemotherapeutics.

We are now offering a comprehensive biobank of highly characterized tissue-derived human gastrointestinal organoids from normal and diseased patients. The intestinal organoid biobank contains over 50 highly characterized intestinal organoids from both normal and diseased patients derived from multiple regions of the digestive systems including small intestine (duodenum, ileum), stomach, rectum and colon.

The 3dGRO® organoids were derived utilizing HUB Organoid Technology. The purchaser of this product shall agree to HUB's [Terms of Use](#), which shall be separately acknowledged and accepted by such purchaser, prior to transfer of this product to purchaser.

ID	Age	Sex	Organ	Catalog Number	Disease Type or Normal
14-881-CR	46	M	Colon (Rectum)	SCC310	Adenocarcinoma
Prep 78-C	52	F	Colon	SCC311	Normal
Prep 80-C	29	F	Colon	SCC312	Normal
Prep 81-C	49	F	Colon	SCC313	Normal
Prep 81-D	49	F	Duodenum	SCC314	Normal
Prep 83-C	45	F	Colon	SCC315	Normal
Prep 83-D	45	F	Duodenum	SCC316	Normal
Prep 84-C	56	F	Colon	SCC317	Normal
Prep 84-D	56	F	Duodenum	SCC318	Normal
Prep 85-C	62	F	Colon	SCC319	Normal
Prep 85-D	62	F	Duodenum	SCC320	Normal
Prep 87-C	21	M	Colon	SCC321	Normal
Prep 87-D	21	M	Duodenum	SCC322	Normal

ID	Age	Sex	Organ	Catalog Number	Disease Type or Normal
Prep 88-C	33	F	Colon	SCC323	Normal
Prep 88-D	33	F	Duodenum	SCC324	Normal
Prep 89-C	55	M	Colon	SCC325	Normal
Prep 89-D	55	M	Duodenum	SCC326	Normal
18-1417-N	45	M	Colon	SCC327	Normal (CRC tissue adjacent)
ht-131-C	58	F	Colon	SCC328	Adenoma
ht-099-C	59	F	Colon	SCC330	Normal
ht-101-C	59	F	Colon	SCC332	Normal
ht-103-C	63	M	Colon	SCC335	Normal
ht-104-C	51	M	Colon	SCC336	Normal
ht-104-I	51	M	Ileum	SCC337	Normal
ht-105-C	19	M	Colon	SCC338	Normal
ht-105-I	19	M	Ileum	SCC339	Normal
ht-118-C	43	M	Colon	SCC340	Normal
ht-119-C	53	F	Colon	SCC341	Normal
ht-127-S	52	F	Stomach	SCC342	Normal
ht-128-S	60	F	Stomach	SCC343	Normal
ht-129-S	60	M	Stomach	SCC344	Normal
ht-130-S	68	M	Stomach	SCC345	Normal
ht-215-CD	59	F	Colon Descending	SCC346	Normal
ht-215-CT	59	F	Colon Transcending	SCC347	Normal
ht-215-D	59	F	Duodenum	SCC348	Normal
ht-222-D	33	M	Duodenum	SCC349	Normal
ht-180-D	38	F	Duodenum	SCC350	Normal
ht-185-D	45	M	Duodenum	SCC351	Normal
ht-188-D	50	F	Duodenum	SCC352	Normal
ht-205-CS	33	F	Colon Sigmoid	SCC353	Normal Margin to Crohn's
ht-205-D	33	F	Duodenum	SCC354	Normal Margin to Crohn's
ht-205-I	33	F	Ileum	SCC355	Normal Margin to Crohn's
ht-213-CA	36	F	Colon Ascending	SCC356	Normal margin to UC
ht-213-CD	36	F	Colon Descending	SCC357	Normal margin to UC
ht-213-CT	36	F	Colon Transcending	SCC358	Normal margin to UC
ht-213-D	36	F	Duodenum	SCC359	Normal margin to UC
ht-214-I	46	M	Ileum	SCC360	Normal margin to UC
ht-214-I-J	46	M	Ileum J-pouch	SCC361	Normal margin to UC
ht-227-D	62	M	Duodenum	SCC362	Normal margin to UC

ID	Age	Sex	Organ	Catalog Number	Disease Type or Normal
ht-230-D	59	F	Duodenum	SCC363	Normal margin to UC
ht-206-CR	24	F	Colon	SCC365	Crohn's
ht-206-I-CR	24	F	Ileum	SCC366	Crohn's
ht-208-CR-CR	23	M	Colon (Rectum)	SCC367	Crohn's
ht-506-I	34	F	Ileum	SCC368	Normal Margin to Crohn's
ht-405-CA-UC	36	F	Colon Ascending	SCC369	Ulcerative Colitis
ht-352-D	55	F	Duodenum	SCC370	Normal

Materials Provided

3dGRO® Gastrointestinal Organoids (SCC310-SCC370): One (1) vial containing ≥ 1500 viable organoids/vial. Store in liquid nitrogen.

Materials Required (Not supplied)

Growth-factor reduced (GFR) Matrigel® (LS356231), TrypLE™ Express Enzyme (1X) (Thermo Fisher, 12605-010), culture ware and media components ([see Protocols](#)).

Storage and Handling

Store in liquid nitrogen. The organoids can be cultured for at least 10 passages after initial thawing without significantly affecting the cell marker expression and functionality.

Quality Control Testing

- Viability: ≥ 1500 viable organoids/vial
- Organoid Growth: Pass
- Cells are tested negative for infectious diseases by a Human Essential CLEAR panel by Charles River Animal Diagnostic Services.
- Mycoplasma Contamination: Negative
- STR Profile: Pass

Representative Data




Figure 1: Various human gastrointestinal organoids grown in Matrigel® domes.




Figure 2: Time course of organoid growth over a 6-day time period shows positive cell proliferation after thaw (SCC318, Duo-84).




Figure 3: Immunocytochemical characterization of human colon organoids (Colon-87, SCC321). Human colon PDOs are positive for colon-specific markers: CA II, CA IV and Mucin5B, posterior hindgut marker: CDX2, stem cell markers: Lgr5 and Sca1 and epithelial markers: TPH-1 and E-Cad.

Duo-85

Colon-87




Figure 4. Adult tissue derived colon (SCC321, Colon-87) and duodenum (SCC320, DUO-85) organoids express SARS-CoV-2 related proteins ACE2 and TMPRSS2.

Protocols

All protocols are performed within a Class II laminar flow biohood and with an aspirator unless otherwise specified. Incubators are humidified and are set to 37 °C and 5% CO₂. PPE should be worn such as gloves, lab coat, and safety glasses. An example protocol for culturing Human Intestinal Organoids is shown below.

Preparing 1X Medium for Human Intestinal Organoids

To prepare 100 mL of 1X Medium for Human Intestinal Organoids prepare the following formulation:

Component	Volume	Final Conc
DMEM/F12 Plus Basal Medium	42.6 mL	1X
L-WRN Conditioned Media Supplement	50 mL	1X
Gibco™ N-2 Supplement (100X)	1 mL	1X
B-27™ Supplement (50X), without vitamin A	2 mL	1X
L-Glutamine Solution (100X), 200 mM	1 mL	1X
HEPES Solution, 1 M in water	1 mL	10 mM
Niacinamide, prepared as 1 M solution in water*	1 mL	10 mM
N-Acetyl-L-cysteine, prepared as 500 mM solution in water	200 µL	1 mM
Human EGF, reconstituted to 200 µg/mL in PBS/0.1% BSA	50 µL	100 ng/mL
[Leu ¹⁵]-Gastrin I, reconstituted to 100 µM in PBS/0.1% BSA	10 µL	10 nM
Prostaglandin E2, reconstituted to 100 µM in DMSO	10 µL	10 nM
A-83-01, reconstituted to 1 mM in DMSO	50 µL	0.5 µM
SB202190, reconstituted to 20 mM in DMSO**	50 µL	10 µM
Penicillin-Streptomycin, 100X solution	1 mL	1X

* SCC327 does not require Niacinamide

** SCC342, SCC343, SCC344, SCC345 requires 10 µM SB431542 instead of SB20290

1. Filter the 1X Medium using a 0.2 µm filter.
2. Use the 1X Medium within 2 weeks and store at 2-8 °C.
3. Y-27632 (ROCK inhibitor) at 10 µM final concentration is recommended to be added to 1X Medium with each fresh media change.
4. For increased viability of human gut organoids during passaging, calculate the total volume of 1X Medium required for passaging and add 2.5 µM CHIR99021 to just that volume. CHIR99021 should only be present for the first two days of each passage. It is not necessary to add CHIR99021 to the 1X Medium during regular media exchanges.
5. Primocin® at 100 µg/mL final concentration may be added to 1X Medium as an additional antimicrobial supplement.

Thawing Organoids

1. Prior to thawing human intestinal organoids, thaw sufficient growth-factor reduced (GFR) Matrigel® for seeding 10 domes of 10 µL per dome. Place thawed Matrigel® on ice.
2. Thaw one vial of human intestinal organoids by submerging 3/5 of vial into a 37 °C water bath until only a sliver of ice remains. Spray the outside of the vial with 70% ethanol or isopropanol.
3. In a sterile tissue culture hood, quickly and gently transfer the vial contents into a 15 mL conical tube containing 4 mL of 1X Medium and 10 µM ROCK inhibitor.
4. Centrifuge the 15 mL conical tube for 5 minutes at 500 x g.

5. Carefully aspirate the media with a pipet, being cautious to avoid disturbing the organoid pellet. Organoids should appear as an opaque layer at the bottom of the tube. If a clear layer of residual Matrigel® is visible overlaying the organoid pellet, remove as much of the Matrigel® as possible without disturbing the organoid pellet.
6. Transfer 100 µL of thawed Matrigel® to the organoid pellet. Quickly and gently resuspend the organoid pellet by pipetting up and down 5 times with a P-100 or P-200 micropipette set to 80 µL, being careful to avoid causing air bubbles.
7. Place the organoid suspension on ice for 5 minutes to cool the Matrigel® plus organoid suspension.
8. Remove the organoid suspension from ice and briefly swirl to mix. Dispense 10 µL domes into one well of the 6-well plate. Work quickly and minimize formation of air bubbles during pipetting.
9. Allow the domes to incubate for 10 minutes in a 37 °C humidified incubator with 5% CO₂.
10. Transfer the 6-well plate containing the organoid domes to a sterile tissue culture hood. Gently add 3 mL of 1X Medium containing 10 µm ROCK Inhibitor and 2.5 µm CHIR99021 into the well containing the organoid domes.
11. Incubate the 6-well plate in a 37 °C humidified incubator with 5% CO₂.
12. After 24-hours, check for recovery by observation by bright field microscope. Live organoids should begin to show rounded or cystic-like morphology. Replace wells with 3 mL fresh 1X Media containing 10 µm ROCK Inhibitor (without CHIR99021) per well.
13. Exchange media every 1-2 days. By day 4, each dome should contain approximately 100-150 organoids.

Passaging Organoids

1. To maintain a continuous culture, human intestinal organoids should be passaged every 5-7 days. Passage of organoids after 7 days is not recommended.
2. Prepare enough GFR Matrigel® for a volume of 10 µL per dome. To maintain approximate equal density, split organoids upon passage at least 1:2 (For example, 10 domes into 20 domes).
3. Prepare ice-cold 1X PBS plus 10 µm ROCK inhibitor.
4. Pipette and transfer the medium in the 6-well plate to a 15 mL conical tube, leaving approximately 800 µL behind in the well.
5. Detach the Matrigel® domes in the well with a P-1000 micropipette. Break organoids by placing the pipette tip perpendicular to the bottom of the well and expel organoids with a scraping motion. Repeat 10 times.
6. Transfer the organoid suspension to a new 15 mL conical tube.
7. Using the same pipette tip, transfer 1 mL of the existing culture medium from the 15 mL tube and rinse and collect any residual organoid suspensions from the well. Transfer to the conical tube containing the organoid suspension.
8. Centrifuge for 5 minutes at 500 x *g*.
9. Carefully aspirate the media, being careful not to disturb the organoid pellet. If a clear layer of Matrigel® is present overlaying the organoid pellet, remove as much Matrigel® as possible. If the Matrigel® layer is not clear but appears opaque, do not attempt to remove.
10. Add 1 mL of ice-cold 1X PBS plus 10 µm ROCK inhibitor to the organoid pellet.
11. Break the organoids with a P-1000 micropipette by placing pipet tip perpendicular to bottom of tube. Repeat 30 times.
12. Add 4 mL of ice-cold 1X PBS plus 10 µm ROCK inhibitor to the organoid suspension.
13. Centrifuge organoid suspension at 4 °C for 5 minutes at 500 x *g*.
14. Carefully aspirate the media with a pipet, removing as much residual Matrigel® as possible without disturbing the opaque layer of organoid pellet.
15. Immediately transfer the appropriate volume of thawed GFR Matrigel® to the organoid pellet with a P-1000 micropipette. Quickly and gently resuspend the organoid pellet about 5 times with a P-1000 tip set to 20 µL below the volume transferred. Avoid causing formation of air bubbles.
16. Place Matrigel® plus organoid suspension on ice for 5 minutes.
17. Briefly swirl the organoid suspension to mix and dispense 10 µL domes into wells of the culture plate. After each 1 minute of dispensing, place suspension on ice for 5 minutes to cool and prevent solidification of the Matrigel®.

18. Allow the domes to incubate for 10 minutes in a 37 °C humidified incubator with 5% CO₂.
19. Remove plate from incubator and add 1X Medium containing 10 µm ROCK inhibitor and 2.5 µm CHIR99021 in sufficient volume to cover the domes (typically 3 mL per well of a 6-well plate, or 750 µL per well of a 12-well plate).
20. Incubate culture in a 37 °C humidified incubator with 5% CO₂.
21. After 24 hours, exchange media with fresh 1X Medium containing 10 µm ROCK Inhibitor (without CHIR99021). Exchange media every 1-2 days.

Single Cell Passaging of Organoids

1. Aspirate medium from each well containing an organoid dome.
2. Add 1 mL of 1X PBS to each well. Pipette up and down 5-10 times with a P-1000 pipet to break up the organoid domes and release the organoids. Do this for each well containing the organoid domes.
3. Combine the organoid suspension from all the wells and transfer to a sterile 15 mL conical tube.
4. Centrifuge at 1100 rpm for 5 minutes at 4 °C.
5. Carefully aspirate the media by connecting a P-200 pipette tip to the end of an aspirating pipette. Be careful not to aspirate the organoid pellet.

Note: A layer of Matrigel® may be visible on top of the organoid pellet. Carefully remove as much of the Matrigel® as possible using a P-200 pipette tip.

6. Add 1 mL of TrypLE™ Express with 10 µm Rock inhibitor to the pellet. Pipette up and down 10 times with a P-1000 pipette tip.
7. Add another 1-2 mL of TrypLE™ Express plus 10 µm ROCK inhibitor and place in a 37 °C incubator for 30-60 minutes. Swirl the tube to mix every 5 minutes.
8. Check a small amount on a hemocytometer at 15 minutes after pipetting up and down a few times to see if organoids are dissociated into single cells. If not, incubate for an additional 10-15 minutes with swirling every 5 minutes and check again. Repeat for up to 60 minutes.
9. Once single cells are observed, remove enough sample for counting on a hemocytometer. Then immediately add 5 mL of plus 10 µm ROCK inhibitor to the rest.
10. Centrifuge at 1100 rpm for 5 minutes at 4 °C. Carefully aspirate the media by connecting a P-200 pipette tip to the end of an aspirating pipette. Be careful not to aspirate the pellet.
11. Wash the pellet with an additional 3-5 mL of 1X Medium plus 10 µm ROCK inhibitor and centrifuge at 1100 rpm for 5 minutes at 4 °C.
12. Carefully aspirate the media by connecting a P-200 pipette tip to the end of an aspirating pipette. Be careful to not aspirate the pellet.
13. Quickly add the appropriate amount of the thawed GFR Matrigel® plus 10 µm ROCK inhibitor to the pellet at the desired cell concentration and resuspend 5-10 times. Avoid generating air bubble during pipetting.
14. Immediately place the conical tube containing the cell suspension on ice for 3-5 minutes to cool down. Remove the tube from ice and spray the outside of the tube with 70% ethanol or isopropanol before putting it back into the TC hood.

Note: Start at 5K cells, 10K cells, 15k cells and 20k cells per 25 µL dome to determine the desired concentration.

Tip: It is important to culture these in 1X Medium plus 10 µm ROCK inhibitor for up to 8-10 days; otherwise, cells will die.

Cryopreservation of Organoids (PDOs)

Important Notes Before Starting:

- We recommend freezing 10 x 10 μ L domes to 1 cryovial. The average organoid density of each dome should be ~90% at the time of freezing. If the density is less than 90%, increase the number of domes to freeze.
 - PDOs may be frozen down in 1X Intestinal Organoid Medium plus 10% DMSO or in 3dGRO® Organoid Freeze Medium.
 - Prepare Thermo Scientific™ Mr. Frosty® Freezing Container and have it ready for storing organoid fragments for freezing down at -80°C .
1. Chill enough 1X PBS plus 10 μ M ROCK inhibitor on ice for at least 30 minutes before starting.
 2. Remove the 6-well plate containing 10 μ L domes from 37°C incubator and place in a TC hood.
 3. Pipette and transfer the medium in the 6-well plate to a 15 mL conical tube, leaving approximately 800 μ L behind in the well.
 4. Detach the Matrigel® domes in the well with a P-1000 micropipette. Break organoids by placing the pipette tip perpendicular to the bottom of the well and expel organoids with a scraping motion. Repeat 10 times.
 5. Transfer the organoid suspension to a new 15 mL conical tube. Using the same pipette tip, transfer 1 mL of the existing culture medium from the 15 mL tube and rinse and collect any residual organoid suspensions from the well. Transfer to the conical tube containing the organoid suspension.
 6. Centrifuge for 5 minutes at $500 \times g$.
 7. Carefully aspirate the media by connecting a P-200 pipette tip to the end of an aspirating pipette. Be careful not to aspirate the organoid pellet. If a clear layer of Matrigel® is present overlaying the organoid pellet, remove as much Matrigel® as possible. If the Matrigel® layer is not clear but appears opaque, do not attempt to remove. This is because there are organoid fragments present in the Matrigel® layer. Repeat centrifugation again at $500 \times g$ for 10 minutes.
 8. Add 1 mL of ice-cold 1X PBS plus 10 μ M ROCK inhibitor to the organoid pellet.
 9. Break the organoids with a P-1000 micropipette by placing pipet tip perpendicular to bottom of tube. Repeat 30 times.
 10. Add 4 mL of ice-cold 1X PBS plus 10 μ M ROCK inhibitor to the organoid suspension.
 11. Centrifuge organoid suspension at 4°C for 5 minutes at $500 \times g$.
 12. Carefully aspirate the media by connecting a P-200 tip to the end of an aspirating pipette, removing as much residual Matrigel® as possible without disturbing the organoid pellet.
 13. Resuspend organoid pellet with freeze medium at 10 μ L domes per mL.
 14. Quickly transfer 1 mL into a cryovial and repeat until all transferred in cryovials.
 15. Quickly transfer cryovials to a Mr. Frosty® Freezing Container and place in a -80°C freezer for 24 hours.
 16. After 24 hours, transfer the cryovial to liquid nitrogen (-135°C) for long term storage.

Whole Mount Immunocytochemistry of Organoids

Important Notes Before Starting:

The following protocol is meant to serve as a guidance for first time users and is based on organoids cultured in 24-well plates. The protocol may be modified and adapted once users are more familiar with the process.

- We recommend using a pair of scissors that have been sterilized with 70% ethanol or isopropanol to cut the ends of P-1000 tips to enlarge the opening. Modified P-1000 tips are used to transfer fixed organoids without shearing them. Do not use serological pipettes as they are too bulky to handle small volumes and organoids may stick to the side of the pipettes.
- During PBS washes, gravity is used to collect organoids. Do not use centrifugation as the centrifugal force will result in mis-sharpened organoids.

- 4% paraformaldehyde performs the dual function of fixing the organoids and to help partially dissolve the Matrigel® and release the organoids. It is important to remove as much of the Matrigel® as possible from the organoids. Matrigel® may result in increased background autofluorescence. The more confluent the organoids are inside the domes at the time of fixing, the more readily the Matrigel® will dissolve.
1. Prepare a 4% paraformaldehyde (PFA) solution by diluting an 8% PFA Solution (Electron Microscopy Sciences, 157-8-100) 1:1 with 1X PBS.
 2. Prepare modified P-1000 and P-200 pipette tips by cutting the ends with a sterilized scissor.
 3. Aspirate the medium from each well containing an organoid dome. Wash each well twice with 1 mL 1X PBS. Aspirate between PBS washes.
 4. Add 1 mL of the 4% PFA solution to each well. Incubate 45-60 minutes at room temperature on a gently rocking or shaking platform. The shaker/rocker will help expedite detaching the Matrigel® domes and the release of the organoids from the Matrigel®.
- Note:** GFR Matrigel® domes will partially dissolve when fixed in PFA. At the end of the incubation period, you will notice that many (but not all) the domes are dislodged and that some of the organoids (but not all) will have been released from the domes.
5. Using the modified P-1000 pipette tips, collect any released organoids along with the fixative solution and transfer the contents to a 50 mL conical tube. Allow the organoids to settle to the bottom of the conical tube by gravity (~10-15 minutes). DO NOT CENTRIFUGE.
 6. In the meantime, add 1 mL 1X PBS per well to the 24-well plate containing the organoid domes. Incubate 10-15 minutes at room temperature. This is done to dilute the PFA in the dome.
 7. Carefully aspirate the fixative from the conical tube containing the released organoids (from step 5) and leave a small amount of liquid behind. This will ensure that the organoid pellet will not be aspirated off.
 8. Using modified P-1000 pipette tips, collect any released organoids along with the PBS solution from each well (from step 6) and transfer the contents to the 50 mL conical tube. Allow the organoids to settle to the bottom of the conical tube by gravity (~10-15 minutes).
 9. Repeat steps 5-8 two more times.
 10. Add 0.8 mL of 1X PBS into each well that contains residual organoid domes.
 11. Carefully aspirate the supernatant from the conical tube containing released organoids. Leave a small amount of liquid behind. This will ensure that the organoid pellet will not be aspirated off.
 12. Add 4.8 mL 1X PBS to the organoid pellet. Swirl the conical tube to resuspend the organoid pellet. Using a modified P-1000 tip, transfer 200 µL of the organoid suspension into each well containing the 0.8 mL volume of residual organoid domes (from step 10).
- Note:** Some organoids may stick to the modified P-1000 tip.
13. If staining will not be performed immediately, seal the 24-well plate containing fixed organoids with Parafilm® and store in the fridge at 2-8 °C for up to 1 month.
 14. When ready to perform ICC, transfer the 24-well plate containing the fixed organoids to a dissecting microscope.
 15. Using modified P-200 tips (from step 2), pipette 1-4 organoids into each well of an 8-well chamber slide. Remove any residual PBS using an unmodified P-200 pipette tip. Avoid accidentally pipetting up the organoids and shearing them through the P-200 tip.
 16. Add 0.4 mL Blocking Buffer (5% horse serum plus 0.5% Triton® X-100 in 1X PBS) to each well of an 8-well chamber slide containing the fixed organoids. Block at 2-8 °C overnight or at room temperature for 2-4 hours.
- Note:** Use the serum from the same species as the host secondary antibody.
17. Using an unmodified P-200 pipet, remove the blocking buffer while tilting the chamber slide. Avoid pipetting the organoids through the P-200 tip.
 18. Prepare primary antibodies or directly conjugated antibodies (300-500 µL) in Blocking Buffer.
 19. Add primary antibodies. Incubate overnight at 2-8 °C on a gently shaking or rotating platform.
 20. Next day, wash 3X with 1X PBS for 10-15 minutes each on the shaking/rotating platform.
- Note:** Do not use centrifugation. Remove the PBS with a P-200 pipet after each wash while tilting the chamber slide.
21. Prepare secondary antibodies (300-500 µL) in Blocking Buffer.

22. Add secondary antibodies. Incubate overnight at 2-8 °C on a gently shaking or rotating platform.
23. Next day, wash with 1X PBS for 10-15 minutes on the shaking/rotating platform.
- Note:** Do not use centrifugation. Remove the PBS with a P-200 pipet after each wash while tilting the chamber slide.
24. Counterstain with DAPI (5 µg/mL in 1X PBS) for 15-20 minutes.
25. Wash 3X with 1X PBS for 10-15 minutes each on the shaking/rotating platform.
- Note:** Do not use centrifugation. Remove the PBS with a P-200 pipet after each wash while tilting the chamber slide.
26. Add 300-400 µL of 1X PBS into each well. Samples are now ready to be imaged on a confocal microscope.

Academic Use Agreement

Subject to local law

THIS PRODUCT MAY ONLY BE USED BY INDIVIDUALS EMPLOYED BY AN ACADEMIC INSTITUTION AND IS INTENDED SOLELY TO BE USED FOR ACADEMIC RESEARCH, WHICH IS FURTHER DEFINED BELOW. BY OPENING THIS PRODUCT, YOU ("PURCHASER") HEREBY REPRESENT THAT YOU HAVE THE RIGHT AND AUTHORITY TO LEGALLY BIND YOURSELF AND/OR YOUR EMPLOYER INSTITUTION, AS APPLICABLE, AND CONSENT TO BE LEGALLY BOUND BY THE TERMS OF THIS ACADEMIC USE AGREEMENT. IF YOU DO NOT AGREE TO COMPLY WITH THESE TERMS, YOU MAY NOT OPEN OR USE THE PRODUCT AND YOU MUST CALL CUSTOMER SERVICE (1-800-645-5476) TO ARRANGE TO RETURN THE PRODUCT FOR A REFUND.

"Product" means 3dGRO® Gastrointestinal Organoids (SCC310-SCC370).

"Academic Research" means any internal in vitro research use by individuals employed by an academic institution. Academic Research specifically excludes the following uses of whatever kind or nature:

- Re-engineering or copying the Product
- Making derivatives, modifications, or functional equivalents of the Product
- Obtaining patents or other intellectual property rights claiming use of the Product
- Using the Product in the development, testing, or manufacture of a Commercial Product
- Using the Product as a component of a Commercial Product
- Reselling or licensing the Product
- Using the Product in clinical or therapeutic applications including producing materials for clinical trials
- Administering the Product to humans
- Using the Product in collaboration with a commercial or non-academic entity

"Commercial Product" means any product intended for: (i) current or future sale; (ii) use in a fee-for-service; or (iii) any diagnostic, clinical, or therapeutic use.

Access to the Product is limited solely to those officers, employees, and students of PURCHASER's academic institution who need access to the Product to perform Academic Research. PURCHASER shall comply with all applicable laws in its use and handling of the Product and shall keep it under reasonably safe and secure conditions to prevent unauthorized use or access.

These use restrictions will remain in effect for as long as PURCHASER possesses the Product. Purchaser may contact HUB at marketing@huborganoids.nl or bd@huborganoids.nl to request a license to use Organoids for other purposes.



3dGRO® organoids were derived utilizing HUB Organoid Technology.

COMMERCIAL OR NON-ACADEMIC ENTITIES INTERESTED IN PURCHASING OR USING THE PRODUCT MUST
CONTACT licensing@emdmillipore.com AND AGREE TO SEPARATE TERMS OF USE PRIOR TO USE OR PURCHASE.

Notice

We provide information and advice to our customers on application technologies and regulatory matters to the best of our knowledge and ability, but without obligation or liability. Existing laws and regulations are to be observed in all cases by our customers. This also applies in respect to any rights of third parties. Our information and advice do not relieve our customers of their own responsibility for checking the suitability of our products for the envisaged purpose.

The information in this document is subject to change without notice and should not be construed as a commitment by the manufacturing or selling entity, or an affiliate. We assume no responsibility for any errors that may appear in this document.

Technical Assistance

Visit the tech service page at SigmaAldrich.com/techservice.

Terms and Conditions of Sale

Warranty, use restrictions, and other conditions of sale may be found at SigmaAldrich.com/terms.

Contact Information

For the location of the office nearest you, go to SigmaAldrich.com/offices.

The life science business of Merck operates
as MilliporeSigma in the U.S. and Canada.

Merck, 3dGRO and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

© 2025 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.

Document Template 20769660 Ver 4.0

00013673 Ver 2.0, Rev 07MAR2025,AV

