

3050 Spruce Street, St. Louis, MO 63103 USA Tel: (800) 521-8956 (314) 771-5765 Fax: (800) 325-5052 (314) 771-5757 email: techservice@sial.com sigma-aldrich.com

Product Information

Trizma[®] base

Catalog Number **T1503** Store at Room Temperature

CAS RN 77-86-1

Synonyms: Tris base, Tris(hydroxyamino)methane, 2-Amino-2-(hydroxymethyl)-1,3-propanediol, Tris(hydroxymethyl)aminomethane, THAM, Trometamol <u>Note</u>: The "Tris" described in this document is **not** the "Tris" used to flame-proof fabric. That compound, Tris(2,3-dibromopropyl)phosphate, has been reported to be a cancer suspect agent.

Molecular Formula: $C_4H_{11}NO_3$ Molecular Weight: 121.14 pK_a (20 °C):¹ 8.3 Melting Point:¹ 171-172 °C

Product Description

Trizma[®] is the registered trademark of Sigma-Aldrich applied to various compounds of Tris(hydroxymethyl)-aminomethane that are prepared by Sigma-Aldrich.

Tris(hydroxyamino)methane, or "Tris" for short, is an established basimetric standard and buffer used in biochemistry and molecular biology.¹ It may be used by itself as a buffer or as a component of mixed buffer formulations.² These different buffer formulations include:

- Tris-EDTA (TE) buffer
- Tris magnesium buffer
- Tris-acetate-EDTA (TAE) buffer
- Tris-borate-EDTA (TBE) buffer
- Tris-buffered saline (TBS)
- Tris-buffered saline with dextrose (TBS-D)
- Tris-glycine buffer
- Tris-phosphate EDTA buffer
- Tris-SDS buffer
- Tris-sucrose
- Tris-Tricine-SDS buffer

When preparing Trizma solutions at a given pH and temperature, it is necessary to choose the proper mixture of Trizma free base and a corresponding Trizma salt to give the desired final pH at the desired temperature.

Trizma has a significant temperature coefficient, which affects the pH of the solution. For a given concentration the following changes are observed:

- From 5 °C to 25 °C, the pH decreases an average of 0.03 pH units per °C.
- From 25 °C to 37 °C, the pH decreases an average of 0.025 pH units per °C.

This product has been used in various studies and application fields, including:

- Gene transfer and expression in cultured cells³
- Investigation of nuclear pores inside cells⁴
- Enzymatic assays^{5,6}
- Tissue sample analysis⁷
- 2-D protein electrophoresis⁸
- Immunohistochemistry⁹
- Fluorescence in situ Hybridization (FISH)¹⁰

Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Safety Data Sheet for information regarding hazards and safe handling practices.

Preparation Instructions

This product is soluble in water (666 mg/mL), yielding a clear, colorless solution.

Storage/Stability

Trizma solutions can be autoclaved.

References

- Gomori, G., Preparation of Buffers for Use in Enzyme Studies. *Methods Enzymol.*, 1, 138-146 (1955).
- Sambrook, J. and Russell, D.W., *Molecular Cloning: A Laboratory Manual* (3rd ed.). Cold Spring Harbor Laboratory Press (Cold Spring Harbor, NY), pp. 5.8, 5.30, 5.43, 5.60, 5.76, 10.25-10.26, 12.75, 12.84, 12.87, 13.52, 16.29-16.31, A1.2-1.3, A1.7-A1.8, A1.17-A1.18, A1.22, A8.42-A8.43 (2001).
- Kampert, S.E. *et al.*, "Gene Transfer and Expression in Animal Cells", in *Handbook of Molecular and Cellular Methods in Biology and Medicine* (L.J. Cseke, P.B. Kaufman, A. Kirakosyan, M.V. Westfall, eds.). CRC Press (Boca Raton, FL), Chapter 32, pp. 557-578 (2011).
- Otsuka, S. et al., "Imaging the Assembly, Structure, and Function of the Nuclear Pore Inside Cells", in *Methods in Cell Biology: Nuclear Pore Complexes* and Nucleocytoplasmic Transport – Methods (V. Doye, ed.). Academic Press / Elsevier (San Diego, CA), Volume 122, pp. 219-238 (2014).
- Butcher, R.H., Am. J. Clin. Pathol., 68(5), 566-569 (1977).

- Cornell, N.W., and Veech, R.L., *Anal. Biochem.*, 132(2), 418-423 (1983).
- 7. Dupret, B. et al., MethodsX, 5, 244-256 (2018).
- Anderson, N.L. *et al.*, "Two-Dimensional Electrophoresis as Applied to Problems of in Vitro and in Vivo Toxicology", in *Methods in Toxicology, Volume 1: Part B – In Vitro Toxicity Indicators* (C.A. Tyson and J.M. Frazier, eds.). Academic Press (San Diego, CA), pp. 463-474 (1994).
- Ramos-Vara, J.A., "Immunohistochemical Methods", in *Making And Using Antibodies: A Practical Handbook, Second Edition* (G.C. Howard and M.R. Kaser, eds.). CRC Press (Boca Raton, FL), Chapter 11, pp. 303-342 (2014).
- Ko, E., and Martin, R., "Fluorescence in situ Hybridization (FISH) Analysis in Human Sperm Cells", in *Diagnostic Cytogenetics* (R.-D. Wegner, ed.). Springer-Verlag (Berlin, Heidelberg), pp. 335-355 (1999).

Trizma is a registered trademark of Sigma-Aldrich Co., LLC.

GCY,RXR,MAM 12/18-1