

3050 Spruce Street Saint Louis, Missouri 63103 USA Telephone 800-325-5832 • (314) 771-5765 Fax (314) 286-7828 email: techserv@sial.com sigma-aldrich.com

# ProductInformation

### **Glucose Oxidase Aspergillus niger**

Product Number **G 6891** Storage Temperature 2-8 °C

## **Product Description**

Enzyme Commission (EC) Number: 1.1.3.4 CAS Number: 9001-37-0 Molecular Weight: 160 kDa (gel filtration)<sup>1</sup> Isoelectric Point:  $4.2^2$ Extinction coefficient:  $E^{1\%} = 16.7 (280 \text{ nm})^3$ Synonyms: GOD, Gox,  $\beta$ -D-Glucose: oxygen 1-oxidoreductase

Glucose oxidase from *Aspergillus niger* is a dimer consisting of 2 equal subunits with a molecular weight of 80 kDa each. Each subunit contains one mole of flavin adenine dinucleotide and one mole of iron. The enzyme is a glycoprotein containing approximately 16% neutral sugar and 2% amino sugars.<sup>1</sup> The enzyme also contains 3 cysteine residues and 8 potential sites for N-linked glycosylation.<sup>4</sup>

Glucose oxidase is capable of oxidizing Daldohexoses, monodeoxy-D-glucoses, and methyl-Dglucoses at varying rates. D-glucose, 2-deoxy-Dglucose, 4-O-methyl-D-glucoses, 6-deoxy-D-glucose, 4-deoxy-D-glucose, 3-deoxy-D-glucose and 3-Omethyl-D-glucose are oxidized at decreasing rates and in the order listed. The pH optimum for glucose oxidase is 5.5, while it has a broad activity range of pH 4-7.<sup>2</sup> Glucose oxidase is specific for  $\beta$ -D-glucose with a K<sub>M</sub> of 33-110 mM.<sup>5,6</sup> Glucose oxidase does not require any activators, but it is inhibited by Ag<sup>+</sup>, Hg<sup>2+</sup>, Cu<sup>2+</sup>, phenylmercuric acetate and p-chloromercuribenzoate. It is not inhibited by the nonmetallic SH reagents: N-ethylmaleimide, iodoacetate, and iodoacetamide.<sup>7</sup>

Glucose oxidase can be utilized in the enzymatic determination of D-glucose in solution. As glucose oxidase oxidizes  $\beta$ -D-glucose to D-gluconolactate and hydrogen peroxide, horseradish peroxidase is often used as the coupling enzyme in glucose determinations. Although glucose oxidase is specific for  $\beta$ -D-glucose, solutions of D-glucose containing  $\alpha$ -D-glucose will mutorotate to  $\beta$ -D-glucose as the  $\beta$ -D-glucose is consumed by the enzymatic reaction.<sup>8</sup>

### **Precautions and Disclaimer**

For Laboratory Use Only. Not for drug, household or other uses.

### Storage/Stability

This product is offered as a solution in 100 mM sodium acetate buffer, pH approximately 4, containing 0.002% thimerosal. It is not recommended that dilute solutions be used as stock solutions.

#### References

- Tsuge, H. J., et al., Purification, properties, and molecular features of glucose oxidase from *Aspergillus niger*. J. Biochem., **78**, 835-843 (1975).
- Pazur, J. H., and Kleppe, K., The oxidation of glucose and related compounds by glucose oxidase from *Aspergillus niger*. Biochemistry, **3**, 578-583 (1964).
- CRC Handbook of Biochemistry and Molecular Biology, Fasman, G. D., ed., CRC Press (Boca Raton, FL:1990), pp. 244.
- 4. Frederick, K. R., et al., Glucose oxidase from *Aspergillus niger*. Cloning, gene sequence, secretion from *Saccharomyces cerevisiae* and

kinetic analysis of a yeast-derived enzyme. J. Biol. Chem., **265**, 3793-3802 (1990).

- 5. Swoboda, B. E. P., and Massey, V., Purification and properties of glucose oxidase from *Aspergillus niger*. J. Biol. Chem., **240**, 2209-2215 (1965).
- Gibson, Q. H., et al., Kinetics and mechanism of action of glucose oxidase. J. Biol. Chem., 239, 3927-3934 (1964).
- Nakamura, S., and Ogura, Y., Mode of inhibition of glucose oxidase by metal ions. J. Biochem., 64, 439-447 (1968).
- 8. Methods of Enzymatic Analysis, 3rd ed., vol. 2, Bergmeyer, H. U., ed., Weinheim (Deerfield Beach, Fl: 1983), pp. 201-202.

TMG/AJH 1/03

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.