



Product Information

RESTRICTION ENDONUCLEASE Hpa II

Product No. **R 0629**

Store at 0 °C to -20 °C

Product Summary

Recognition Sequence: 5'C/CGG'3

Activity: 10,000 units/ml

Cutting: 100%

Ligation: >95%

Recutting: >95%

No degradation detected with >40 units for 16 hrs.

Fold over digestion: 640 (40 units x 16 hrs.)

Package Size: 500 units

Unit Definition

One unit is the enzyme activity that completely cleaves 1 µg λDNA in 1 hr. at 37 °C in a total volume of 25 µl of Buffer SL for restriction endonucleases.

Specificity

Hpa II recognizes the sequence C/CGG and generates fragments with 5'-cohesive termini.¹ Hpa II is sensitive to methylation at the internal cytosine. Hpa II does not cleave ^mCCGG and C^mCGG.

Comments

Digestion Buffer SL is supplied as a 10x concentrate. Information for heat inactivation Hpa II is not available.

Hpa II Storage and Dilution Buffer Composition

20 mM Tris-HCl

50 mM KCl

0.5 mM EDTA

5 mM 2-mercaptoethanol

50% (v/v) glycerol

pH 7.5

1x Digestion Buffer SL (B 3782) Composition for Hpa II: 100% Digestion at 37 °C.

10 mM Tris-HCl

10 mM MgCl₂

1 mM dithioerythritol (DTE)

pH 7.5

Quality Control Testing

Absence of unspecific endonuclease activities:

1 µg λDNA is incubated for 16 hrs. in 50 µl buffer SL with excess of Hpa II.

Ligation and Recutting Assay

Hpa II fragments, obtained by complete digestion of 1 µg λ DNA, are adjusted to pH 7.5 at 20 °C. The Hpa II fragments are then ligated with 0.5 unit T4-DNA ligase at pH 7.5 at 20 °C. A 10 µl reaction mixture, incubated for 16 hrs. at 20 °C, contained 0.5 unit T4-DNA ligase, 66 mM Tris-HCl, 5 mM MgCl₂, 1 mM dithiothreitol and 1 mM ATP.

The degree of ligation and subsequent recutting with Hpa II to yield the typical pattern of λ x Hpa II fragments is determined.

References

1. Garfin, D., et al., Biochem. Biophys. Res. Commun., **59**, 108 (1974).