

ProductInformation

POLYMYXIN B - AGAROSE Sigma Prod. No. P1411

CAS NUMBER: N/A

PHYSICAL DESCRIPTION:

Matrix: Cross-linked 4% beaded agarose (activated by cyanogen bromide)

Attachment: Polymyxin B (Sigma product no. P1004) is attached through an amino group, with a 1-atom

spacer

Ligand immobilized: Approximately 1 mg per mL

Binding capacity: 200-500 μg lipopolysaccharide from Escherichia coli serotype 0128:B12 per mL

Form: suspension in 50% glycerol containing 0.02% sodium azide

STORAGE / STABILITY AS SUPPLIED:

The resin should be stable at least one year stored at 2-8°C. Freezing the suspension will damage the agarose beads.

GENERAL REMARKS:

Polymyxin B is a strongly cationic cyclic antibiotic isolated from Bacillus polymyxa. It inhibits growth of microorganisms, including E. coli. It is a variable mixture of B1 and B2, with approximate molecular weight 1450. Polymyxins combine with cell membranes and disrupt normal permeability to small molecules. Polymyxin B and the other polymyxin antibiotics act primarily by binding membrane phospholipids and disrupting the cytoplasmic membrane...Polymyxin B has a bactericidal action on most Gram-negative bacilli, except Proteus spp....not active against Neisseria species, ... most fungi and Gram-positive bacteria. Polymyxin B has a bacterial properties and applications have been reported Including its use in a bacterial permeabilization reagent.

Usage of Polymyxin B Sulfate is often reported in terms of unit activity: "one unit is contained in 0.000119 mg of the second International Standard Preparation (1969) of Polymyxin B Sulfate which contains 8403 units per mg."²

Because of its affinity for bacterial cell membranes (which are lipopolysaccharides, also called endotoxins), Polymyxin B has been immobilized on agarose and used to remove endotoxins from solutions.^{7,8} Histamine agarose has been used for the same purpose.^{9,10}

P1411 is tested for its ability to bind LPS from E. coli serotype 0128:B12. A general protocol for use of polymyxin B-resins to remove endotoxins is given below.

POLYMYXIN B - AGAROSE Sigma Prod. No. P1411

SUGGESTED PROTOCOL:11

- 1. The resin should be washed with 3-5 column volumes of buffer to remove the glycerol storage solution. (Endotoxin-free 0.1 M ammonium bicarbonate buffer pH 8.0 may be used.) Pack immobilized polymyxin B (10 mL gel) in a disposable polypropylene column.
- 2. Equilibrate the column with 100 mL buffer.
- 3. Dissolve sample, for example, bovine serum albumin (BSA, 400 mg, 80 E.U./mg) in 2.0 mL buffer, and apply to the column.
- 4. After the protein solution has completely entered the gel, elute with endotoxin-free buffer at pH 7.8, at a flow rate of approximately 0.2 ml/minute, while collecting 1.0-mL fractions.
- 5. Collect the fractions containing protein (BSA) and immediately lyophilize. Endotoxin testing of purified BSA should indicate a level of less than 1 E.U./mg protein. If the actual level is higher, recycle the protein through this resin to remove further quantities of endotoxin.

The polymyxin B agarose can be regenerated at least 10 times. For regeneration, wash the column with 50 mL 1% sodium deoxycholate solution, followed by 50 mL endotoxin-free water. The resin can be stored in 25% ethanol at 4°C.

REFERENCES:

- 1. Data for Biochemical Research, 3rd Ed., Dawson et al. eds. (Oxford Press, 1986), P. 301.
- 2. *Martindale: The Extra Pharmacopoeia,* 30th Ed., J. Reynolds et al., eds. (Pharmaceutical Press, 1993), p. 194.
- 3. J. Gen. Microbiol., 10, 491 (1954).
- 4. Cornu, J., *Ann. Microbiol.*, 131B, 121 (1980). (In French).
- 5. Lerner, H.R. et al., *Physiol. Plant*, 57, 90 (1983).
- 6. Schupp, J.M. et al., *Biotechniques*, 19, 18 (1995).
- 7. Molvig, H. and Baek, L., Scand. J. Immunol., 26, 611-619 (1987).
- 8. Issekutz, A.C., J. Immunol. Methods, 61, 275-281 (1983).
- 9. Hou, K.C. and Zaniewski, R., *Biochim. Biophys. Acta*, 1073, 149-154 (1991).
- 10. Sofer, G., *Bio/Technology*, Dec. 1984, p. 1035.
- Hermanson, G.T. et al., Immobilized Affinity Ligand Techniques (Academic Press, 1992), pp. 354-355.

Sigma warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.