

Product Information

D-Luciferin

Synthetic

L9504

Product Description

Synonyms: (S)-2-(6-Hydroxy-2-benzothiazolyl)-2-thiazoline-4-carboxylic acid, 4,5-Dihydro-2-(6-hydroxy-2-benzothiazolyl)-4-thiazolecarboxylic acid, Firefly Luciferin

Molecular Formula: C₁₁H₈N₂O₃S₂

Molecular Weight: 280.3

CAS Number: 2591-17-5

 λ_{max} : 268 nm, 330 nm¹Extinction Coefficient: E^{mm}= 7.04 (268 nm), 18.2 (330 nm)¹Solvent: N₂-sparged ethanol

D-Luciferin is a naturally occurring compound in organisms capable of bioluminescence, such as fireflies. The enzyme luciferase oxidizes D-luciferin to produce a light-emitting excited state molecule.^{2,3}

The excitation and emission spectra for D-Luciferin have been published. The excitation is pH-dependent, with a maximum of 327 nm at pH 4 and 385 nm at pH 11. The emission profile is identical at both pH's, with a maximum at 537 nm.⁴ The dependence of the bioluminescence of the luciferase-luciferin system on Zn²⁺ concentration has been published.⁵

ATP can be measured with a reagent made up of luciferin and luciferase from firefly. Discussions of extraction buffers for releasing ATP from bacteria and tissues have been published.^{6,7} If D-luciferin is used for assaying the concentration of ATP in cell lysates, it is important to know if ATPases are present. These enzymes must be inactivated in the extraction process so that the ATP is not destroyed. Heat or low pH are usually used and do not affect the integrity of the ATP.

Several theses⁸ and dissertations⁹⁻¹⁵ cite use of product L9504 in their research protocols.

Precautions and Disclaimer

This product is for R&D use only. Not for drug, household, or other uses. Please consult the Safety Data Sheet for information regarding hazards and safe handling practices.

Solubility

D-luciferin is tested for solubility in DMSO at 10 mg/mL.

To solubilize luciferin free base in water:

- Bubble nitrogen through distilled water.
- Stir in the luciferin. It will not dissolve.
- Add as close to 1 equivalent as possible of NaHCO₃. Allow the NaHCO₃ to stir slowly at room temperature. It will take about 30 minutes to solubilize.
- The solution should be a faint yellow solution at a pH of approximately 6.5. If too much sodium bicarbonate is added, the solution will be too alkaline, and the luciferin will oxidize and form a green solution.

Preparation Instructions

Do **not** use Tris-HCl to prepare solutions of D-luciferin. If the application requires a Tris buffer, Tris-acetate buffer is recommended instead.

Solutions of D-luciferin in ice-cold 0.1 M Tris-acetate buffer, pH 7.5-7.75, prepared at 4 °C and protected from light (amber bottle in ice bath) are stable for 8-24 hours at 4 °C. In general, it is not recommended to store D-luciferin solutions frozen, although several publications report freezing of luciferin stock solutions under conditions that we have not ourselves tested:

- 10 mM stock solution in water:methanol (10 mg of D-luciferin in 2.57 mL water + 1 mL of MeOH), stored at -20 °C¹⁶
- 100 mM stock solution at -80 °C¹⁷

References

1. White, E. H. et al., *J. Am. Chem. Soc.*, 85(3), 337-343 (1963).
2. Adams, S. T., Jr., and Miller, S. C., *Curr. Opin. Chem. Biol.*, 21, 112-120 (2014).
3. Syed, A. J., and Anderson, J. C., *Chem. Soc. Rev.*, 50(9), 5668-5705 (2021).
4. Bowie, L. J., *Methods Enzymol.*, 57, 15-28 (1978).
5. DeLuca, M., and McElroy, W. D., *Methods Enzymol.*, 57, 3-15 (1978).
6. Chappelle, E. W. et al., *Methods Enzymol.*, 57, 65-72 (1978).
7. Karl, D. M., *Methods Enzymol.*, 57, 85-94 (1978).
8. Afrin, Farzana, "Investigation of the RNA-interference pathway of *Toxoplasma gondii*". University of Windsor, M.Sc. thesis, 29 (2018).
9. Causse, Sébastien, "Etude de la dynamique d'activation de la transcription des gènes par l'ARN Polymerase 2" ("Study of the activation dynamics of gene transcription by RNA polymerase 2"). Université Pierre et Marie Curie - Paris VI, Ph.D. dissertation, p. 158 (2011).
10. Han, Dongjun, "Role of Oct4 in pXEN cell differentiation and MET process". Humboldt-Universität zu Berlin, Dr. rer. nat. dissertation, p. 29 (2019).
11. Ahmed, Mohammed, "Modulation du microenvironnement tumoral pour contourner la résistance du glioblastome aux traitements". ("Modulation of the tumor microenvironment to overcome glioblastoma treatment resistance"). Université Paris-Saclay, Ph.D. dissertation, p. 126 (2020).
12. Capri, Joseph Randolph, "Integrative and Quantitative Mass Spectrometry (IQ-MS): Developing a unified sample workup for multi-omics analysis". University of California Los Angeles, Ph.D. dissertation, pp. 27 (2020).
13. Sattler, Lisa-Marie, "Das Polyzystische Ovarialsyndrom als Autoimmunerkrankung - sind natürliche Gonadotropin-Releasing-Hormon-Rezeptor-Autoantikörper von diagnostischer Relevanz?" ("Polycystic ovary syndrome as an autoimmune disease - are natural gonadotropin-releasing hormone receptor autoantibodies of diagnostic relevance?"). Charité-Universitätsmedizin Berlin, Dr. med. dissertation, p. 19 (2022).
14. Damm, K., *Methods in Neurosciences*, 22, 265-276 (1994).
15. Lee, B.-H. et al., Chapter 43, "High-Throughput Screening of Arabidopsis Mutants with Deregulated Stress-Responsive Luciferase Gene Expression", in *Luminescence Biotechnology: Instruments and Applications* (Van Dyke, K. et al., eds.). CRC Press (Boca Raton, FL), p. 562 (2002).

Notice

We provide information and advice to our customers on application technologies and regulatory matters to the best of our knowledge and ability, but without obligation or liability. Existing laws and regulations are to be observed in all cases by our customers. This also applies in respect to any rights of third parties. Our information and advice do not relieve our customers of their own responsibility for checking the suitability of our products for the envisaged purpose.

The information in this document is subject to change without notice and should not be construed as a commitment by the manufacturing or selling entity, or an affiliate. We assume no responsibility for any errors that may appear in this document.

Technical Assistance

Visit the tech service page at [SigmaAldrich.com/techservice](https://www.sigmaaldrich.com/techservice).

Terms and Conditions of Use

Warranty, use restrictions, and other conditions of sale may be found at [SigmaAldrich.com/terms](https://www.sigmaaldrich.com/terms).

Contact Information

For the location of the office nearest you, go to [SigmaAldrich.com/offices](https://www.sigmaaldrich.com/offices).

The life science business of Merck KGaA, Darmstadt, Germany operates as MilliporeSigma in the U.S. and Canada.

MilliporeSigma, and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

© 2024 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.