

- Removes >65% of serum albumin
- Lower non-specific binding than traditional Cibacron Blue
- Easy, 15-minute protocol
- Can be used for multiple species without protocol change



Improves resolution of lower abundance proteins by electrophoresis, chromatography, or MS

The presence of albumin in serum or plasma samples can decrease the resolution of one-dimensional electrophoresis (1DE), two-dimensional electrophoresis (2DE), chromatography, or mass spectrometry systems.

The Montage Albumin Deplete Kit provides a fast, convenient, and reproducible method for improving resolution by reducing albumin levels in serum or plasma samples.

Improve Resolution of Lower Abundance Proteins

The kit's centrifugal columns are prepacked with an affinity resin that removes >65% of albumin from typical human serum samples. The affinity resin has been specially formulated to bind albumin without binding significant amounts of other serum or plasma proteins. This reduction in albumin levels can improve the detection of lower abundance proteins by subsequent analytical techniques.

Process Samples in Just 15 Minutes


Simply equilibrate the affinity column; add the sample; spin to bind the albumin; wash the column; and collect the eluate. If desired, the bound fraction (albumin and other retained proteins) can be recovered using a stripping solution supplied with the kit.

The kit is optimized for the depletion of human albumin. However, depletion has been demonstrated with other mammalian species including bovine, canine, goat, mouse, rabbit, and rat samples with no protocol changes.

Ensure the Reproducibility of Your Results

Preparing reagents and affinity columns in the lab increases the potential for variability in your results. The Montage Albumin Deplete Kit comes ready-to-use, with prepared reagents and columns. Each lot of kits is QC tested for reproducible albumin depletion and low non-specific binding.





Equilibrate column




Add sample



Deplete albumin



Wash to reduce non-specific binding



Elute bound fraction (optional)


= Albumin

△ = Other serum/plasma proteins

Depletion of Albumin from Human Plasma

Typical Human Albumin Depletion Measured by Radial Immunodiffusion Assay

Human Plasma Volume (µL)	Albumin Depletion (%)
20	85
50	75
75	65
100	58



Lane 1: Molecular Weight Markers


Lane 3: Unprocessed Plasma

Lane 5: Flow-through (Albumin-depleted)
Plasma

Lane 7: Bound Fraction

Figure 1. Human plasma sample (20 µL) was diluted with 180 µL of equilibration buffer and processed with the Montage Albumin Deplete Kit. Samples of the flow-through and bound fractions were then separated by 1D SDS-PAGE and stained with colloidal Coomassie™ blue. To accurately assess albumin removal, all sample loads were volume normalized to the undepleted starting sample. Radial immunodifussion assay revealed 85% albumin reduction. Gel densitometry revealed 10% non-specific protein binding to the resin.

Reproducibility



Lane 1: Molecular Weight Markers

Lanes 2-11: Replicates


Figure 2. Ten replicates of a human plasma sample (75 μL) were diluted with 125 μL of equilibration buffer and individually processed with the Montage Albumin Deplete Kit. Flow-through (albumin-depleted) fractions were collected; acetone-precipitated to further concentrate the sample; and then resuspended in SDS-PAGE sample buffer. Samples were separated by 1D SDS-PAGE and stained with colloidal Coomassie blue.

Radial immunodetection assay revealed an average of 66% albumin reduction across the 10 separate depletions with a standard deviation of 3.2 and a CV of 2.1%.

Depletion of Albumin from Animal Sera

Typical Animal Albumin Depletion Measured by Gel Densitometry

Species	Albumin Depletion (%)	Non-Specific Binding (%)
Mouse	42.0	18.7
Rat	65.0	17.8
Bovine	39.0	9.7
Rabbit	46.0	12 1



Rat Serum

Lane 1: Molecular Weight Markers

Lane 3: Unprocessed Rat Serum


Lane 5: Flow-through (Albumin-depleted)

Serum

Lane 7: Bound Fraction

Figure 3. Plasma samples (75 μ L) from a range of animal species were diluted with 125 μ L equilibration buffer and processed using the Montage Albumin Deplete Kit. Bound material was eluted from the columns using the stripping buffer supplied with the kit and separated by 1D SDS-PAGE.

Two-Dimensional Electrophoresis of Human Plasma






Figure 4. A 20 μL plasma sample was dissolved in IPG loading buffer and analyzed by 2DE (left). The same plasma sample was then processed using the Montage Albumin Deplete Kit and analyzed (right).

Magnification of selected gel regions demonstrates the increased visualization of lower abundance proteins achieved with the Montage Albumin Deplete Kit.



Specifications	
Kit Contents	Affinity columns, micro-centrifuge tubes, buffers, and reagents for processing either 4 or 24 samples
Storage	2 to 8 °C
Resin Volume	200 µL packed bed resin/column

Binding Capacity

Resin binding capacity: ~ 40 mg of human serum albumin per gram of moist

resin under static equilibrium binding conditions

Dynamic binding

Typically 2 mg human serum albumin capacity/column: Non-specific binding: Typically < 14% from human serum

Ordering Information		
Description	Qty/Pk	Catalogue No.
Montage Albumin Deplete Kit	Process 4 samples Process 24 samples	LSKA D00 04 LSKA D00 24

To Place an Order or Receive **Technical Assistance**

For additional information call your nearest Millipore office: In the U.S. and Canada, call toll-free 1-800-MILLIPORE (1-800-645-5476)

In the U.S., Canada and Puerto Rico, fax orders to 1-800-MILLIFX (1-800-645-5439)

Internet: www.millipore.com

Tech Service:

www.millipore.com/techservice



Now you can buy Millipore products online @

www.millipore.com/purecommerce

Millipore Worldwide

AUSTRALIA

Tel. 1 800 222 111 or (02) 9888 8999 Fax (02) 9878 0788

Tel. (01) 877 89 26 Fax (01) 877 16 54

BALTIC COUNTRIES Tel. +358 9 804 5110 Fax +358 9 256 5660

BELGIUM AND LUXEMBOURG

Tel. +32 2 726 88 40 Fax +32 2 726 98 84

Tel. (011) 5548-7011 Fax (011) 5548-7923

CANADA

Tel. 1-800-645-5476 Fax 1-800-645-5439

CHINA, PEOPLE'S REPUBLIC OF

Beijing:

(8610) 8519 1250 (8610) 6518 1058 Fax (8610) 8519 1255

Guangzhou: Tel. (8620) 8755 4049 Fax (8620) 8752 0172

Hong Kong: Tel. (852) 2803 9111 Fax (852) 2513 0313 Shanghai:

Tel. (8621) 5306 9100 Fax (8621) 5306 0838

CZECH REPUBLIC Tel. 02-2051 3841 Fax 02-2051 4298

DENMARK Tel. 70 10 00 23

Fax 70 10 13 14 EASTERN EUROPE,

C.I.S., AFRICA, MIDDLE EAST AND GULF Life Sciences Division:

Tel. +33 3 88 38 9536 Fax +33 3 88 38 9539 BioPharmaceutical Division: Tel. +43 1 877-8926 Fax +43 1 877-1654

Laboratory Water Division: Tel. +33¹ 30 12 7000 Fax +33 1 30 12 7180

FINLAND

Tel. (09) 804 5110 Fax (09) 256 5660

FRANCE

Tel. 01 30 12 7000 Fax 01 30 12 7180

GERMANY

Tel. (06196) 494-0 Fax (06196) 43901

HUNGARY

Tel. 01 381 0433 01 381 0434 01 209 3232 Fax 01 209 0295

Tel. (91) 80-839 46 57 Fax (91) 80-839 63 45 **IRELAND**

Life Sciences and Laboratory Water Divisions: Tel. +44[']1923 816375 Fax +44 1923 818297 BioPharmaceutical Division: Tel. (021) 883 666 Fax (021) 883 048

ITALY

Vimodrone (Milano): Tel. (02) 25.07.81 Fax (02) 26.50.324 Roma:

Tel. (06) 52.03.600 Fax (06) 52.95.735

JAPAN

Tel. (03) 5442-9711 Fax (03) 5442-

9736 Life Sciences 9737 BioPharm. 9734 Jah Water

KORFA

Tel. (822) 551-0990 Fax (822) 551-0228

MALAYSIA

Tel. 603-7957-1322 Fax 603-7957-1711

MEXICO

Tel. (525) 576-9688 Fax (525) 576-8706

THE NETHERLANDS Tel. 076-5022000 Fax 076-5022436

NORWAY Tel. 22 67 82 53 Fax 22 66 04 60 **POLAND**

Tel. 22-669 12 25 22-663 70 31 Fax 22-663 70 33

PUERTO RICO

Tel. (787) 273-8495 Fax (787) 747-6553

SINGAPORE

Tel. (65) 842 1822 Fax (65) 842 4988

SPAIN AND PORTUGAL Tel. +34 917 283 960 Fax +34 917 292 909

SWEDEN

Tel. 08-628 6960 Fax 08-628 6457

SWITZERLAND

Tel. (01) 908-30-60 Fax (01) 908-30-80

TAIWAN

Tel. 886-2-2792-9333 Fax 886-2-2792-6555

Tel. 01923 816375 Fax 01923 818297

U.S.A

Tel. (781) 533-6000 Fax (781) 533-3110

IN ALL OTHER COUNTRIES

Millipore Intertech (U.S.A.) Tel. +1 (781) 533-8622 Fax +1 (781) 533-8630

Millipore and Montage are trademarks of Millipore Corporation.

Coomassie is a trademark of BASF Aktiengesellschaft. Lit. No. PF1770EN00 8/02 02-188 Printed in U.S.A.

© 2002 Millipore Corporation, Bedford, MA. All rights reserved.

