

Automation of Small Volume PCR Purification

Chris Barbagallo¹, Libby Kellard¹, Peter J. Rapiejko¹, John McCool², and Jennifer Goode²

¹Millipore Corporation, Life Sciences Division, Danvers, MA 01923; and ²Tecan USA, Research Triangle Park, NC 27709

Abstract

The expansion of small-scale laboratory processes to industrial-scale environments has necessitated the development of new instrument and product technologies. Aspects of importance include small volume sample recovery, compatibility with laboratory automation and the means to deliver a high throughput process. The polymerase chain reaction (PCR) is a key element of many current and emerging genomics applications, including sequencing, microarrays and genotyping. Each of these technologies is dependent upon rapid and robust post-PCR clean-up methods that effectively remove salts, primers and dNTPs. Building upon Millipore's size-exclusion based purification technology, the MontageTM PCR_{μ96} plate has been developed specifically to meet the demands of modern PCR-based applications. The tapered well design of the Montage PCR_{μ96} plate enables 20 μl recovery volumes for a wide range of input reaction sizes. The Montage PCR_{μ96} plate rapidly delivers the high quality PCR products required for the most demanding of downstream applications while maintaining high DNA recovery. Full automation of the PCR clean-up is demonstrated using a Tecan Genesis[®] Te-MOTM Multi-Channel Pipetting Option, a stand-alone 96 channel pipettor.


Materials and Method

- 20 μL aliquots of the PCR samples were diluted with 80 μL TE buffer
- 100 μL from each well were then transferred to the Montage PCR_{μ96} plate on the vacuum manifold
- Vacuum was applied at 20" Hg until the wells were dry (7 minutes)
- Samples were then resuspended in 20 μL TE buffer via incubation for 5 minutes and subsequent mixing of 15 μL for 75 cycles
- Purified samples were then transferred to a V-bottom collection plate

*This protocol was carried out using two vacuum manifolds on the Te-MO deck.

** In all cases, some volume is lost to the membrane during resuspension and transfer to the collection plate (approximately 1 - 2 μL).

PCR Product Recovery by SYBR[®] Quantitation and Primer Removal




Starting Concentrations			
	301 bp	657 bp	1159 bp
	8.6 ng/μL	25.5 ng/μL	21.7 ng/μL
Percent Recovery			
	301 bp	657 bp	1159 bp
Avg	84.9%	94.8%	99.2%
CV	8.7%	8.5%	6.7%
n	96	96	96

Primer Removal	
>	>99%
n	183

Recovery data was calculated using a fluorometric assay with SYBR Green I nucleic acid gel stain (Molecular Probes, Inc). Primer removal was performed by measuring the retention of a fluoresceinated 20-mer in the purified PCR samples.

Agarose Gel Analysis of Purified PCR Samples




1.5% agarose gel analysis of PCR products after purification using Montage-PCR_{μ96} on the Te-MO. 5 μL of samples loaded on gel.


The Montage PCR_{μ96} Plate Designed for Automation



Montage PCR_{μ96} Method in GeminiTM Version 4.0



Deck setup for PCR Cleanup on the Tecan Genesis Te-MO



Conclusions

- Fast processing time utilizing the Te-MO's 96-channel pipetting capabilities (under 20 minutes to purify two plates of 96 samples)
- Small volume recovery for sample concentration
- Efficient removal of PCR primers and salts
- Micro-well format compatible with sequencing, genotyping, and microarray production.