Supelco_®

Analytical Products





Melting Point Standards

Melting point (MP) is defined as the temperature at which a chemical substance changes from solid to the liquid state. Melting point determination is widely used for the characterization of pure chemicals and pharmaceutical drugs. It is also used in R & D labs for chemical quality control (QC), and quality assurance (QA) to identify solid crystalline substances and to check their purity. Analytical laboratories involved in QC/QA need to calibrate their melting point instrumentation regularly to ensure that their instruments are in accordance with the specific requirements defined by their local, national, and international standards laboratories. Standards used in the calibration of melting point instrumentation are called melting point standards.

We offer melting point standards ranging from Phenyl salicylate, melting point 41.2-43.2 °C to Sodium methanesulfonate, melting point 353.1-355.1 °C.

Technical Benefits:

- Traceable to primary standards (LGC, London)
- Grade: Analytical Standard
- Shelf Life: 2 years
- Provided with Certificates of Analysis and Safety Data Sheet
- Two modes for melting point evaluation:
 - The Pharmacopeia mode: neglects that the furnace temperature is different (higher) during the heating process than the sample temperature, meaning that the furnace temperature is measured rather than the sample temperature. Therefore, the pharmacopeia melting point depends strongly on the heating rate.
 - The Thermodynamic mode: The thermodynamic melting point is the physically correct melting point. This value does not depend on heating rate or other parameters.



Melting Point Determination as per European Pharmacopeia's (2.2.14.) Requirements:

The melting point is determined by Capillary method as described in European Pharmacopeia 10.5 (2.2.14.). The pharmacopeia's requirements for melting point determination are listed:

- Use capillaries with outer diameters ranging from 1.3-1.5 mm and wall thicknesses from 0.1-0.3 mm.
- Apply a constant heating rate of 1 °C/min.
- The recorded temperature represents the temperature of the heating stand, which can be an oil bath or a metal block, in which the thermocouple is positioned.

METTLER TOLEDO Melting Point Standards

METTLER TOLEDO standards are provided with melting points specific to METTLER TOLEDO instrument measured in Thermodynamic mode and Pharmacopeia modes. Also, these standards are useful for the calibration of METTLER TOLEDO melting point instruments. The value for the melting point is an average of 6 to 12 measurements with a METTLER TOLEDO MP 90 and MP 70 excellence melting point instrument (the measuring range is from ambient temperature up to 400 °C) that is calibrated with primary standards.

Table 1: METTLER TOLEDO Calibration substances

| Cat. No. | Description | Composition | Melting Point | Pack size | |
|----------|---------------------------------------------------------------------|----------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------|-----------|--|
| 44770 | METTLER TOLEDO Calibration substance ME 30034252, Phenyl salicylate | Phenyl salicylate | 41.2 °C (Standard Deviation = 0.2 °C) (For measurements in thermodynamic mode) | | |
| | | | 42.3 °C (Standard Deviation = 0.2 °C) (For measurements with a METTLER TOLEDO thermosystem FP in pharmacopeia mode at 1.0 °C /min heating rate) | 5g | |
| | | | 43.2 °C (Standard Deviation = 0.2 °C) (For measurements with a METTLER TOLEDO MP melting point system in pharmacopeia mode at 1.0 °C /min heating rate) | | |
| 73664 | METTLER TOLEDO Calibration substance ME 18870, Benzophenone | Benzophenone | 47.4 °C (Standard Deviation = 0.2 °C) (For measurements in thermodynamic mode at 0.2 °C /min heating rate) | | |
| | | | 48.5 °C (Standard Deviation = 0.2 °C) (For measurements with a METTLER TOLEDO thermosystem FP900 in pharmacopeia mode at 1.0 °C /min heating rate) | 5g | |
| | | | 49.4 °C (Standard Deviation = 0.2 °C) (For measurements with a METTLER TOLEDO MP melting point system in pharmacopeia mode at 1.0 °C /min heating rate) | | |
| 77634 | METTLER TOLEDO Calibration substance ME 51143093, Vanillin | Vanillin | 80.9 °C (Standard Deviation = 0.2 °C) (For measurements in thermodynamic mode) | | |
| | | | 82.0 °C (Standard Deviation = 0.2 °C) (For measurements with a METTLER TOLEDO thermosystem FP in pharmacopeia mode at 1.0 °C /min heating rate) | 5g | |
| | | | 82.9 °C (Standard Deviation = 0.2 °C) (For measurements with a METTLER TOLEDO MP melting point system in pharmacopeia mode at 1.0 °C /min heating rate) | | |
| 73983 | METTLER TOLEDO Calibration substance ME 18555, Benzoic acid | Benzoic acid | 121.9 °C (Standard Deviation = 0.2 °C) (For measurements in thermodynamic mode) | | |
| | | | 123.0 °C (Standard Deviation = 0.2 °C) (For measurements with a METTLER TOLEDO thermosystem FP in pharmacopeia mode at 1.0 °C /min heating rate) | 5g | |
| | | | 123.9 °C (Standard Deviation = 0.2 °C) (For measurements with a METTLER TOLEDO MP melting point system in pharmacopeia mode at 1.0 °C /min heating rate) | | |

Table 1: METTLER TOLEDO Calibration substances (continued).

| Cat. No. | Description | Composition | Melting Point | Pack size |
|----------|-------------------------------------------------------------------------------------------|-----------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------|------------|
| 15809 | METTLER TOLEDO Calibration substance ME 30130597, p-Toluamide, | p-Toluamide | 159.2 °C (Standard Deviation = 0.3 °C) (For measurements in thermodynamic mode) | |
| | | | 160.2 °C (Standard Deviation = 0.3 °C (For measurements with a METTLER TOLEDO thermosystem FP in pharmacopeia mode at 1.0 °C /min heating rate) | 5g |
| | | | 161.2 °C (Standard Deviation = 0.3 °C) (For Measurements with a METTLER TOLEDO MP melting point system in pharmacopeia mode at 1.0 °C /min heating rate) | |
| 41131 | METTLER TOLEDO Calibration substance ME 51143091, Saccharin | Saccharin | 228.3 °C (Standard Deviation = 0.3 °C) (For measurements in thermodynamic mode) | 5g |
| | | | 229.3 °C (Standard Deviation = 0.3 °C) (For measurements with a METTLER TOLEDO thermosystem FP in pharmacopeia mode at 1.0 °C /min heating rate) | |
| | | | 230.3 °C (Standard Deviation = 0.3 °C) (For measurements with a METTLER TOLEDO MP melting point system in pharmacopeia mode at 1.0 °C /min heating rate) | |
| 94993 | METTLER TOLEDO Calibration substance ME 30130598, Methyltriphenylphosphoniumbromide | Methyltriphenylphosphoniumbromide | 232.6 °C (Standard Deviation = 0.3 °C) (For measurements in thermodynamic mode) | |
| | | | 233.7 °C (Standard Deviation = 0.3 °C) (For measurements with a METTLER TOLEDO thermosystem FP900 in pharmacopeia mode at 1.0 °C /min heating rate) | 5g |
| | | | 234.6 °C (Standard Deviation = 0.3 °C) (For measurements with a METTLER TOLEDO MP melting point system in pharmacopeia mode at 1.0 °C /min heating rate) | |
| 75035 | METTLER TOLEDO Calibration substance ME 18872, Caffeine | Caffeine | 236.2 °C (Standard Deviation = 0.2 °C) (For measurements in thermodynamic mode) | |
| | | | 237.3 °C (Standard Deviation = 0.2 °C) (For measurements with a METTLER TOLEDO thermosystem FP in pharmacopeia mode at 1.0 °C /min heating rate) | 5g |
| | | | 238.2 °C (Standard Deviation = 0.2 °C (For measurements with a METTLER TOLEDO MP melting point system in pharmacopeia mode at 1.0 °C /min heating rate) | |
| 42123 | METTLER TOLEDO Calibration substance ME 30130599, Sodium acetate anhydrous | Sodium acetate anhydrous | 329.1 °C (Standard Deviation = 0.4 °C) (For measurements in thermodynamic mode) | |
| | | | 330.1 °C (standard deviation = 0.4 °C) (For measurements with a METTLER TOLEDO thermosystem FP in pharmacopeia mode at 1.0 °C /min heating rate) | 5 g |
| | | | 331.1 °C (Standard Deviation = 0.4 °C) (For measurements with a METTLER TOLEDO MP melting point system in pharmacopeia mode at 1.0 °C /min heating rate) | |

Table 1: METTLER-TOLEDO Calibration substances (continued).

| Cat. No. | Description | Composition | Melting Point | Pack size |
|----------|---------------------------------------------------------------------------|-------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------|
| 49143 | METTLER TOLEDO Calibration substance ME 51143095, Potassium nitrate | Potassium nitrate | 333.2 °C (Standard Deviation = 0.5 °C) (For measurements in thermodynamic mode) | |
| | | | 334.3 °C (Standard Deviation = 0.5 °C) (For measurements with a METTLER TOLEDO thermosystem FP in pharmacopeia mode at 1.0 °C /min heating rate) | - 5g |
| | | | 335.2 °C (Standard Deviation = 0.5 °C) (For measurements with a METTLER TOLEDO MP melting point system in pharmacopeia mode at 1.0 °C /min heating rate) | |
| 04229 | METTLER TOLEDO Calibration substance ME 30130610, Sodium methanesulfonate | Sodium methanesulfonate | 353.1 °C (Standard Deviation = 0.4 °C) (For measurements in thermodynamic mode) | |
| | | | 354.2 °C (Standard Deviation = 0.4 °C) (For measurements with a METTLER TOLEDO thermosystem FP900 in pharmacopeia mode at 1.0 °C /min heating rate) | - 5g |
| | | | 355.1 °C (Standard Deviation = 0.4 °C) (For measurements with a METTLER TOLEDO mp melting point system in pharmacopeia mode at 1.0 °C /min heating rate) | |

Table 2. Melting point standards: The melting point is measured with an average of 6 - 12 measurements with a Buechi B-545 instrument that is calibrated with primary standard

| Cat. No. | Description | Composition | Melting Point | Pack size |
|----------|-----------------------------------|------------------------------|--------------------------------------------|-----------|
| 50296 | Melting point standard 47-49 °C | Benzophenone | 47.87 °C (±0.3 °C) (Thermodynamic Mode) | 1 g |
| 30230 | | Benzophenone | | 5 g |
| 77402 | Melting Point Standard 69-71 °C | 1-Heptadecanecarboxylic acid | 69.2 °C (±0.3 °C) (Thermodynamic Mode) | 1 g |
| 77402 | | | | 5 g |
| 01422 | Melting point standard 79-81 °C | Naphthalene | 80.0 °C (±0.3 °C) (Thermodynamic Mode) | 1 g |
| 76170 | Melting point standard 121-123 °C | Benzoic acid | 122.0 °C (±0.3 °C) (Thermodynamic Mode) | 5 g |
| 42183 | Melting point standard 182-184 °C | 4-Methoxybenzoic acid | 183.2 °C (±0.3 °C) (Thermodynamic Mode) | 250 mg |
| 42103 | Melting point standard 102-104 C | 4-Methoxyberizoic acid | | 1 g |
| 41019 | Melting point standard 235-237 °C | 1,3,7-Trimethylxanthine | 236.0 °C (±0.3 °C) (Thermodynamic Mode) | 1 g |
| | | 1,5,7 mmetryradiume | | 5 g |
| 67372 | Melting point standard 283-286 °C | Anthraguinone | 284.4 °C (±0.3 °C) (Thermodynamic Mode) | 250 mg |
| | | Anumaquinone | | 1 g |

To place an order or receive technical assistance in the U.S. and Canada, call toll-free 1-800-645-5476 For other countries across Europe and the world, please visit: **EMDMillipore.com/offices**For Technical Service, please visit: **EMDMillipore.com/techservice**

MilliporeSigma 400 Summit Drive Burlington, MA 01803

EMDMillipore.com



