

Product Information

Acetyl coenzyme A lithium salt

≥93% (HPLC)

A2181


Product Description

CAS Registry Number: 32140-51-5 ('x' lithium salt);
75520-41-1 (trilithium salt)Molecular Formula: C₂₃H₃₈N₇O₁₇P₃S (free acid);
C₂₃H₃₅N₇O₁₇P₃SLi₃ (trilithium salt)Formula Weight: 809.57 (free acid basis); 827.37
(trilithium salt basis)

Synonyms: Acetyl CoA Lithium; Acetyl-S-CoA Lithium

E^M (260 nm) = 16,400 (water)¹E^M (259 nm) = 15,400 (0.1 M phosphate buffer, pH 7)E^M (232 nm) = 8,700 (water)¹

Structure:



Coenzymes comprise a class of molecules, generally derived from vitamins, which function catalytically in enzyme systems.² The acetic acid moiety, which is bound by a high-energy bond (free energy 34.3 kJ/mol) to the -SH group of Coenzyme A, is a precursor to fatty acids, steroids and other naturally occurring compounds, such as terpenes and acetogenins present in plants.^{3,4} In the transfer reaction by Acetyl CoA of the C2 acetyl fragment, either the carboxyl group or the methyl group may react (electrophilic vs. nucleophilic reaction, respectively). The biosynthetic pathways for Acetyl CoA have been illustrated.³

Several methods of preparation and methods for the determination of Acetyl CoA and other CoA derivatives have been described in the literature.⁵⁻⁷ Several theses⁸ and dissertations⁹⁻¹⁷ have cited use of A2181 in their research protocols.

Precautions and Disclaimer

For R&D use only. Not for drug, household, or other uses. Please consult the Safety Data Sheet for information regarding hazards and safe handling practices.

Storage/Stability

It is suggested to store this product at -20 °C, under desiccator conditions, as this product is moisture-sensitive.

Solubility

This product is soluble in water at 100 mg/mL.

Acetyl CoA is generally stable in neutral and moderately acidic solutions.¹ Acetyl CoA hydrolyzes in strong acid, and hydrolyzes more rapidly in alkaline solutions. Although we have not tested long-term storage of acetyl CoA solutions, aqueous solutions of sodium acetyl CofA have been reported to be stable for at least 6 months at -20 °C, with a suggestion to store at -70 °C "if possible".¹⁸

References

1. Dawson, R.M.C. *et al.*, *Data for Biochemical Research*, 3rd edition. Clarendon Press (Oxford, UK), pp. 116-117 (1993).
2. Meyers, R.A., *Molecular Biology and Biotechnology: A Comprehensive Desk Reference*. VCH Publishers, Inc. (New York, NY), p. 960 (1995).
3. Lehninger, A.L., *Biochemistry*. Worth Publishers, Inc., p. 534 (1972).

4. Scott, T., and Eagleson, M., *Concise Encyclopedia: Biochemistry*, 2nd edition. Walter de Gruyter Press (New York, NY), p. 4 (1988).
5. Simon, E.J., and Shemin, D., *J. Am. Chem. Soc.*, **75(10)**, 2520 (1953).
6. Stadtman, E.R., *Meth. Enzymol.*, **3**, 931-941 (1957).
7. Bergmeyer, H.-U. et al., *Methods of Enzymatic Analysis*, 2nd ed., Vol. 4 (Bergmeyer, H.-U., ed.). Verlag Chemie/Academic Press, pp. 1998-2000 (1974).
8. Manoharan, Naveenan, "Characterization of mitochondrial respiration in Atlantic salmon (*Salmo salar*) fed with glutamate and succinate enhanced diet". Norwegian University of Science and Technology, M.Sc. thesis, p. 19 (2018).
9. Doepel, Lorraine Stephania, "Peripartum Performance and Amino Acid Metabolism of the Dairy Cow". University of Alberta, Ph.D. dissertation, p. 228 (2002).
10. Workman, Aspen M., "Regulation of Bovine Herpesvirus 1 (BHV-1) Productive Infection by Cellular Transcription Factors". University of Nebraska Lincoln, Ph.D. dissertation, pp. 37, 69, 140 (2011).
11. Moreno, Alison E., "Structural and Biochemical Studies of Histone Acetyltransferases and Sirtuin Enzymes". Johns Hopkins University, Ph.D. dissertation, p. 93 (2015).
12. Gonzalez, Rhiannon Wold, "Development of Immunological Methods for the Detection of Micropollutants in Fresh Water Samples". Texas Woman's University, Ph.D. dissertation, p. 33 (2016).
13. Cliff, Timothy Scott, "Defining Roles of Myc as a Developmental Regulator of Cell Fate Through Metabolic Control". University of Georgia, Ph.D. dissertation, p. 150 (2017).
14. Kappler, Lisa Claudia Charlotte, "Method development for valid high-resolution profiling of mitochondria and Omics investigation of mitochondrial adaptations to excess energy intake and physical exercise". Eberhard Karls Universität Tübingen, Dr. rer. nat. dissertation, p. 15 (2018).
15. Beebe, Jennifer Diane, "Analysis of a Novel Hdac8-H1.3 Complex in Several Human Carcinoma Cell Lines". Indiana University, Ph.D. dissertation, p. 38 (2019).
16. Lacerda de Albuquerque, Ralph, "Predictors of Maximal Oxygen Consumption During Exercise, and Ecological and Behavioral Correlates in Lizards and Mammals". University of California Riverside, Ph.D. dissertation, p. 68 (2019).
17. Weidling, Ian Walter, "Respiratory chain deficiency alters cellular proteostasis and triggers Alzheimer's disease-like tau alterations". University of Kansas, Ph.D. dissertation, p. 177 (2019).
18. Lake, R.A., and Owen, M.J., *Methods Mol. Biol.*, **7**, 23-33 (1991).

Notice

We provide information and advice to our customers on application technologies and regulatory matters to the best of our knowledge and ability, but without obligation or liability. Existing laws and regulations are to be observed in all cases by our customers. This also applies in respect to any rights of third parties. Our information and advice do not relieve our customers of their own responsibility for checking the suitability of our products for the envisaged purpose.

The information in this document is subject to change without notice and should not be construed as a commitment by the manufacturing or selling entity, or an affiliate. We assume no responsibility for any errors that may appear in this document.

Technical Assistance

Visit the tech service page at [SigmaAldrich.com/techservice](https://www.sigmaaldrich.com/techservice).

Standard Warranty

The applicable warranty for the products listed in this publication may be found at [SigmaAldrich.com/terms](https://www.sigmaaldrich.com/terms).

Contact Information

For the location of the office nearest you, go to [SigmaAldrich.com/offices](https://www.sigmaaldrich.com/offices).

The life science business of Merck KGaA, Darmstadt, Germany operates as MilliporeSigma in the U.S. and Canada.

MilliporeSigma, and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.
© 2022 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.

A2181pis Rev 03/22 MAC,GCY

**MILLIPORE
SIGMA**