

3050 Spruce Street, St. Louis, MO 63103 USA Tel: (800) 521-8956 (314) 771-5765 Fax: (800) 325-5052 (314) 771-5757 email: techservice@sial.com sigma-aldrich.com

# **Product Information**

# TrueGel3D Polymer

Modified for fast gelling (FAST-DEXTRAN), allows cell recovery

Catalog Number **TRUEDEXF** Storage Temperature –70 °C

### **Product Description**

FAST-DEXTRAN polymer is a natural degradable polymer functionalized with maleimide thiol-reactive groups, which react rapidly with crosslinkers (PEG or CD cell-degradable crosslinker) to form biomimetic hydrogels in few seconds. The TrueGel3D buffer, pH 5.5, provided in the kit ensures homogenous mixing of polymer and crosslinker by slowing down gel formation and simultaneously controlling osmotic conditions. Most cells are not affected by a short exposure to buffer at pH 5.5. TrueGel3D hydrogel formed from dextran polymers are dissolved by treatment of TrueGel3D enzymatic cell recovery solution for post culture analysis.

The chemically defined hydrogel formed from FAST-DEXTRAN polymers allows complete control over gel stiffness. The polymers are transparent and mimic the natural extracellular matrix environment. FAST-DEXTRAN polymers are used when cells need to obtain appropriate medium within a few minutes, or when the application requires fast gelation, as in the case of bioprinting.

## Components

FAST-DEXTRAN
3 × 170 μL
Lyophilized
Each tube contain 30 mmol/L of reactive groups
Catalog Number TRU-FDE

 TrueGel3D buffer, pH 5.5, 10× 600 μL Catalog Number TRUBUF-55PH

#### **Precautions and Disclaimer**

This product is for R&D use only, not for drug, household, or other uses. Please consult the Safety Data Sheet for information regarding hazards and safe handling practices.

#### **Preparation Instructions**

- Centrifuge the vial to make sure entire material is at the bottom of the tube.
- Add 175 µL of water to make a 30 mmol/L concentration of maleimide groups.
- Vortex until all material is dissolved.
- Incubate the tube on ice for 5 minutes.
- Briefly vortex and centrifuge the tube. Note: Keep on ice while in use.
- FAST-DEXTRAN is now ready to use.

### Storage/Stability

- The lyophilized powders may be stored unopened in the original bottles at -70 °C for up to one year.
- FAST-DEXTRAN is stored at -70 °C after reconstitution.
- Buffers are stored at 4 °C for short term (<2 months) and between -20 °C to -70 °C for long term.
- Water can be stored between –70 °C and room temperature.

BG,MAM 01/18-1