Instructions for Use

Wright-Giemsa Stain, Modified

Procedure No. WG

Intended Use

Wright-Giemsa solution is intended for use in staining blood films or bone marrow films. Solutions are for "In Vitro Diagnostic Use." For professional use only. The data obtained from this manual, qualitative procedure demonstrates cellular elements in blood and bone marrow films of human specimens. This data when reviewed in conjunction with other diagnostic tests and information may be used as an aid to diagnosis of viral and bacterial infections, allergenic reactions, and other pathological states.

Wright-Giemsa Stain is a modified Romanowsky stain used for differentially staining the cellular elements of blood. When blood films are stained as described in the procedure, the white blood cell nucleus and cytoplasm take on the characteristic blue or pink coloration. The combination of purified eosin and thiazine dyes in the product eliminates inconsistent staining and yields lot-to-lot reproducible chromogenic response.

Reagents

Wright-Giemsa Stain, Modified (Cat. No. WG: WG16-500ML; WG32-1L; WG80-2.5L; WG128-4L)
Wright-Giemsa stain, modified, 0.4% w/v, C.I. 45380 and 52010, buffered at pH 6.8, in methanol

Special Materials Required but Not Provided

- Phosphate Buffer (Cat. No. P3288-1VL or P3288-12VL)
A mixture of sodium phosphate and potassium phosphate, 0.0083 M/L, pH 7.2
- Methanol, Acetone Free (Cat. No. M1775-1GA)
- Microscope / Slides / Coverslips

Storage and Stability

Store Wright-Giemsa solutions at room temperature (18-26°C). Reagent label bears expiration date.

Store Phosphate Buffer and Methanol at room temperature (18-26°C). Store working phosphate buffer solution at 2-8°C.

Deterioration

Discard Wright-Giemsa stain solution if a precipitate develops or water artifacts appear in red cells. Discard the working phosphate solution if turbidity or visible bacterial growth is present.

Preparation

Wright-Giemsa solution is supplied ready to use.

The phosphate buffer should be prepared by diluting 1 vial of buffer to 1 gallon or 3.8 liters of deionized water. Mix well to dissolve.

Methanol is ready to use.

Precautions

These IVDs are intended for in vitro diagnostic use in a clinical laboratory environment. These IVDs are for professional use by qualified personnel only. Sigma-Aldrich IVDs may be operated by laboratory personnel who are trained to handle human specimens that can be infectious, use microscopes and other laboratory equipment and have color perception and visual acuity to distinguish colors and other objects under a microscope.

Normal precautions exercised in handling laboratory reagents should be followed. Dispose of waste observing all local, state, provincial or national regulations.

Procedure

Specimen Collection

No known test method can offer complete assurance that blood samples or tissue will not transmit infection. Therefore, all blood derivatives or tissue specimens should be considered potentially infectious.

Fresh whole blood films or fresh films from blood anticoagulated with EDTA should be used. Thicker films (e.g., bone marrow) will usually require longer staining times.

Notes

- For greater cellular detail, staining time may be increased. Color (tones of blue or red) may be varied by increasing or decreasing time in deionized water.
- Rapid (30 seconds) staining is not recommended for bone marrow. For those preparations 1-3 minutes in stain and 2-10 minutes in deionized water gives satisfactory results.
- For batch staining, slide racks and dishes such as supplied by Sakura Finetek for Tissue-Tek® are recommended since this system allows for vertical placement of slides.
- Color can be varied by increasing or decreasing time in deionized water. Bone marrows should be stained for at least 3 minutes and buffered for 3-10 minutes.
- Staining times outlined in the above procedures have given satisfactory results in our laboratories. Individual preferences may dictate adjustment of these times.
- The deionized water recommended for this procedure should be between pH 6.8 to 7.2. If the deionized water is more acidic than pH 6.8, a phosphate buffer is recommended.
- In case the slides cannot be immediately stained, the slides can be fixed in methanol for 1 minute and dried. Stain as soon as possible.
- Positive control slides should be included in each run.

Procedure

Dip Method (Rapid)

- Place approximately 50 mL Wright-Giemsa Stain in a Coplin jar.
- Fill another Coplin jar with water or phosphate buffer.
- Place thoroughly dried blood film, feather edge DOWN, in Wright-Giemsa Stain for approximately 30 seconds.
NOTE: Rapid dipping for 5-10 seconds may reduce water artifacts on films that are not thoroughly dried.
- Remove slide from stain and place in deionized water or phosphate buffer, pH 6.8-7.2, feather edge DOWN, for approximately 1-10 minutes. DO NOT AGITATE SLIDE WHILE IT IS IN DEIONIZED WATER.
- Rinse briefly in running deionized water and air dry thoroughly before evaluation.

Horizontal Staining Method

- Place thoroughly dried blood film on an appropriate staining rack.
- Flood slide with 1-2 mL Wright-Giemsa Stain.
- After 1 minute, add an equal volume of deionized water or phosphate buffer, pH 6.8-7.2, and mix thoroughly by gently blowing on slide.
- After 1-3 minutes, thoroughly rinse with deionized water and air dry.

Performance Characteristics

Nuclei will be varying shades of purple. Cytoplasmic staining will be varying shades of blue to light pink. Fine reddish to lilac granules may be present in cytoplasm of some cell types. Basophils will demonstrate dark blue black granules in the cytoplasm. Eosinophils will demonstrate orange granules in the cytoplasm. Red blood cells should be pink to orange.¹

If observed results vary from expected results, please contact Sigma-Aldrich Technical Service for assistance.

Analytical Performance Characteristics

The analytical performance results for the given tests conducted on all target structures, confirm 100% sensitivity, specificity, and repeatability.

Cat. No	Product Description	Target	Intra-assay Specificity	Intra-assay Sensitivity	Inter-assay Specificity	Inter-assay Sensitivity
WG	Wright-Giemsa Stain, Modified	Nuclei	3 of 3	3 of 3	3 of 3	3 of 3
		Cytoplasm	3 of 3	3 of 3	3 of 3	3 of 3
		Granules	3 of 3	3 of 3	3 of 3	3 of 3

Warnings and Hazards

Refer to Safety Data Sheet and product labeling for any updated risk, hazard or safety information.

WG16, WG32, WG80, WG128:

H225 Highly flammable liquid and vapor.

H301 + H311 + H331 Toxic if swallowed, in contact with skin or if inhaled.

H319 Causes serious eye irritation.

H370 Causes damage to organs (Eyes, Central nervous system).

P210 Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.

P280 Wear protective gloves/ protective clothing/ eye protection/ face protection/ hearing protection.

P301 + P310 IF SWALLOWED: Immediately call a POISON CENTER/ doctor.

P303 + P361 + P353 IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water.

P304 + P340 + P311 IF INHALED: Remove person to fresh air and keep comfortable for breathing. Call a POISON CENTER/ doctor.

P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

If during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.

Symbol Definitions

Symbols as defined in EN ISO 15223-1:2021

	Manufacturer		Catalogue Number
	Consult Instructions for Use		Batch Code
	Authorized Representative in the European Community/ European Union		European Union Declaration of Conformity (defined in IVDR 2017/746)
	Use-by Date		In vitro diagnostic medical device
	Temperature Limit		Caution
	Date of Manufacture		Importer
	Indicates the Authorised Representative in Switzerland		

References

- Hematology: Principles and Procedures, Sixth Edition, Brown AB, Lea & Febiger, Philadelphia 1993 p100-102

Contact Information

To place an order, please visit our web site at SigmaAldrich.com. For Technical Service, please visit the tech service page on our web site at SigmaAldrich.com/techservice.

Revision History

Rev. 4.0	2014
Rev. 5.0	2016
Rev. 6.0	2023

Transferred to new template with current branding. Specified for professional use in intended use and precautions. Revised intended use to align with IVDR guidelines. Updated Material Safety Data Sheet to Safety Data Sheet. Updated contact information. Removed instruction to follow CLSI for specimen collection. Removed EN 980 and changed to EN ISO 15223-1:2021 for symbols. Added adverse event contact information. Updated literature references. Added a revision history table. Added Warnings and Hazards. Added CH-REP information.

Sigma-Aldrich, Inc.,
3050 Spruce Street,
St. Louis, MO 63103 USA
an affiliate of Merck KGaA,
Darmstadt, Germany
+1(314) 771-5765

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Merck KGaA,
64271
Darmstadt,
Germany

MDSS CH GmbH
Laurenzenvorstadt 61
5000 Aarau
Switzerland

The Initial M and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

© 2023 Merck KGaA, Darmstadt, Germany and/or its affiliates. All rights reserved.

The life science business of Merck operates as MilliporeSigma in the U.S. and Canada.

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MDSS GmbH