

#### sigma-aldrich.com

 3050 Spruce Street, St. Louis, MO 63103 USA

 Tel: (800) 521-8956
 (314) 771-5755

 Fax: (800) 325-5052
 (314) 771-5757

 email: techservice@sial.com sigma-aldrich.com

# **Product Information**

Glucose Oxidase, Type VII from Aspergillus niger

Catalog Number **G2133** Storage Temperature –20 °C

CAS RN 9001-37-0 EC 1.1.3.4 Synonyms: Gox; β-D-Glucose:oxygen 1-oxidoreductase

# **Product Description**

Glucose oxidase catalyzes the oxidation of  $\beta$ -D-glucose to form D-glucono-1,5-lactone and hydrogen peroxide.

 $\beta$ -D-glucose +  $O_2 \rightarrow$  D-glucono-1,5-lactone +  $H_2O_2$ 

Glucose oxidase can be utilized for the enzymatic determination of D-glucose in solution. As glucose oxidase oxidizes  $\beta$ -D-glucose to D-gluconolactate and hydrogen peroxide, horseradish peroxidase is often used as the coupling enzyme in glucose determinations. Although glucose oxidase is specific for  $\beta$ -D-glucose, solutions of D-glucose can be quantified, as  $\alpha$ -D-glucose will mutorotate to  $\beta$ -D-glucose as the  $\beta$ -D-glucose is consumed by the enzymatic reaction.<sup>1</sup>

Molecular mass:<sup>2</sup> 160 kDa (gel filtration)

Glucose oxidase from *Aspergillus niger* is a dimer consisting of 2 equal subunits each with a molecular mass of 80 kDa. Each subunit contains one molecule of flavin adenine dinucleotide (FAD) and one molecule of iron. The enzyme is a glycoprotein containing approximately 16% neutral sugar and 2% amino sugar.<sup>2</sup> The enzyme also contains 3 cysteine residues and 8 potential sites for *N*-linked glycosylation.<sup>3</sup>

Extinction coefficient:<sup>4</sup>  $E^{1\%} = 16.7$  (280 nm)

Isoelectric point:<sup>5</sup> 4.2

Optimal pH:<sup>2</sup> 5.5 (broad activity range of pH 4-7)

Inhibitors:  $Ag^+$ ,  $Hg^{2+}$ , and Cu<sup>2+</sup> ions, phenylmercuric acetate and *p*-chloromercuribenzoate inhibit glucose oxidase. Nonmetallic sulfhydryl reagents, such as *N*-ethylmaleimide, iodoacetate, and iodoacetamide, are not inhibitors.<sup>6</sup>

Substrates: Glucose oxidase is relatively specific for  $\beta$ -D-glucose (K<sub>M</sub> of 33–110 mM)<sup>7,8</sup> It also oxidizes D-aldohexoses, monodeoxy-D-glucoses, and methyl-D-glucoses at varying rates. The following substrates are listed in decreasing order of oxidation rate: D-glucose, 2-deoxy-D-glucose, 4-O-methyl-D-glucose, 6-deoxy-D-glucose, 4-deoxy-D-glucose, 3-deoxy-D-glucose, 3-O-methyl-D-glucose

One publication has examined kinetics data and structural data to postulate on the chemical mechanism of action of glucose oxidase from *A. niger.*<sup>9</sup> The role of FAD in the oxidation of glucose, as catalyzed by glucose oxidase from *A. niger*, has been investigated using electrospray ionization mass spectrometry.<sup>10</sup> The crystal structure of a partially deglycosylated form of glucose oxidase from *A. niger* has been reported.<sup>11</sup>

This product is supplied as a lyophilized powder containing phosphate buffer salts and sodium chloride.

Protein content: ≥60% protein

Specific activity: ≥100,000 units/g solid (without added oxygen). If the reaction mixture is saturated with oxygen, the activity may increase by up to 100%.

Unit definition: One unit will oxidize 1.0  $\mu$ mole of  $\beta$ -D-glucose to D-gluconolactone and H<sub>2</sub>O<sub>2</sub> per minute at pH 5.1 at 35 °C.

Other enzyme activities:

Catalase: ≤10 Sigma units/mg protein Amylase, maltase, glycogenase, invertase, and galactose oxidase: lot-specific results are reported on the Certificate of Analysis.

# **Precautions and Disclaimer**

This product is for R&D use only, not for drug, household, or other uses. Please consult the Safety Data Sheet for information regarding hazards and safe handling practices.

## **Preparation Instructions**

This enzyme is soluble (0.2 mg/mL) in 50 mM sodium acetate buffer, pH 5.1, yielding a clear solution. One publication reports preparation of stock solutions of this product in 50 mM sodium acetate, pH 5.0, at 3.7 mg/mL.<sup>12</sup> Another publication reports preparation of 2 mg/mL stock solutions of this product in a buffer of 80 mM PIPES, pH 6.9, 1 mM EGTA, and 1 mM MgCl<sub>2</sub>, and storage of stock solutions in aliquots at -20 °C.<sup>13,14</sup> We have not tested this latter condition.

### Storage/Stability

This product is stored at -20 °C and is stable for at least 5 years, when properly stored and unopened.

### References

- Bergmeyer, H.U., in *Methods of Enzymatic* Analysis, Volume 2 (3<sup>rd</sup> ed.). Academic Press (Deerfield Beach, FL), pp. 201-202 (1983).
- Tsuge, H. et al., J. Biochem., 78(4), 835-843 (1975).
- Frederick, K.R. et al., J. Biol. Chem., 265(7), 3793-3802 (1990).
- Fasman, G.D., in CRC Handbook of Biochemistry and Molecular Biology, CRC Press (Boca Raton, FL), p.244 (1990).

- Pazar, J.H., and Kleppe, K., *Biochemistry*, 3(4), 578-583 (1964).
- Nakamura, S., and Ogura, Y., *J. Biochem.*, 64(4), 439-447 (1968).
- Swoboda, B.E.P., and Massey, V., *J. Biol. Chem.*, 240(5), 2209-2215 (1965).
- Gibson, Q.H., et al. J. Biol. Chem., 239(11), 3927-3934 (1964).
- 9. Wohlfahrt, G., *et al. Mol. Cell. Biochem.*, **260(1)**, 69-83 (2004).
- 10. Wang, Y., Chem. Sci., 9(3), 594-599 (2018).
- 11. Hecht, H.J. *et al.*, *J. Mol. Biol.*, **229(1)**, 153-172 (1993).
- 12. Lyubimova, A. *et al.*, *Nat. Protoc.*, **8(9)**, 1743-1758 (2013).
- Bechstedt, S., and Brouhard, G.J., "Fluorescence-Based Assays for Microtubule Architecture", Chapter 21, in *Microtubules, In Vitro* (2<sup>nd</sup> ed). *Methods in Cell Biology*, Vol. 115 (J.J. Correia and L. Wilson, eds.), Academic Press, pp. 343-354 (2013).
- Uppalapati, M., et al., "Microtubule Motors in Microfluidics", Chapter 13, in *Methods in Bioengineering: Biomicrofabrication and Biomicrofluidics* (J.D. Zahn, ed.). Artech House (Boston & London), p. 316 (2010).

VLR,RBG,GCY,MAM 10/18-1

©2018 Sigma-Aldrich Co. LLC. All rights reserved. SIGMA-ALDRICH is a trademark of Sigma-Aldrich Co. LLC, registered in the US and other countries. Sigma brand products are sold through Sigma-Aldrich, Inc. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see product information on the Sigma-Aldrich website at www.sigmaaldrich.com and/or on the reverse side of the invoice or packing slip.