

Product Information

COMPETENT CELLS, LOW EFFICIENCY JM109

Product No. **C 2582**

Store at -70°C

Product Description

JM109 cells are made competent for transformation by a modification of the method of Hanahan.

Provides $>10^7$ colony forming units per μg of control DNA (pGEM-3Z diluted to $0.1\text{ ng}/\mu\text{l}$ in TE buffer).

Storage

Always avoid multiple freeze-thaw cycles or exposure to frequent temperature changes. These fluctuations can greatly alter product stability. Do not re-freeze unused thawed aliquot.

Genotype

endA1, *recA1*, *gyrA96*, *thi*, *hsdR17* (r_K^- , m_K^+), *relA1*, *supE44*, *D(lac-proAB)*, [*F=*, *traD36*, *proAB*, *lacI^qZΔM15*].

Storage Buffer

100	mM KCl
45	mM MnCl_2
10	mM CaCl_2
3	mM HACoCl_3
10	mM KOAc
10	% glycerol (w/v)
5	% sucrose (w/v)
pH is adjusted to 6.4 using dilute HCl	

Transformation Protocol

1. Remove a tube of frozen cells from -70°C and thaw on ice.
2. To a chilled 17 x 100 mm polypropylene tube, add $50\ \mu\text{l}$ of competent cells. Note: Use of standard microcentrifuge tubes will reduce the transformation efficiency by 2-fold.
3. Add 1-50 ng of plasmid DNA ligation mixture (in a volume that should not exceed $10\ \mu\text{l}$) per $100\ \mu\text{l}$ of competent cells. Mix by gently tapping the tube.
4. Immediately place the tube on ice for 10 minutes.
5. Heat shock the cells for 45-50 seconds in a 42°C water bath (critical step). **DO NOT SHAKE.**
6. Immediately place the tubes on ice for 2 minutes.
7. Add $900\ \mu\text{l}$ of SOC medium (Product No. S1797, available as part of #COMP-T, competent cell preparation kit) and incubate for 60 minutes at 37°C with shaking (approximately 225 rpm).
8. Dilute the experimental reaction as necessary and spread $100\text{-}200\ \mu\text{l}$ on plates containing the appropriate antibiotic.

Reference

Hanahan, D., J. Mol. Biol. **166**, 557 (1983)