

3050 Spruce Street Saint Louis, Missouri 63103 USA Telephone 800-325-5832 • (314) 771-5765 Fax (314) 286-7828 email: techserv@sial.com sigma-aldrich.com

ProductInformation

Pyruvic acid

Product Number **1 07360** Storage Temperature 2-8 °C

Replacement for Product Number P1656

Product Description

Molecular Formula: $C_3H_4O_3$ Molecular Weight: 88.06 CAS Number: 127-17-3 Density: 1.27 g/ml¹ Melting Point: 11.8 C¹ pK_A : 2.49² Based on the density 1.27 g/ml and a purity of 100%, the molarity would be 14.4. A purity of 98% has a

molarity of 14.1. Synonyms: 2-oxopropanoic acid; α-ketopropionic acid; acetylformic acid; pyroracemic acid

Pyruvic acid is the end product of the glycolysis pathway, whereby glucose is converted to pyruvate with the production of ATP. In the mitochondria of aerobic organisms, pyruvate is converted to acetyl coenzyme A, which in turn is oxidized completely to CO₂. When oxygen is not present in sufficient quantities, pyruvate is metabolized to lactate. In anaerobic organisms such as yeast, pyruvate is converted to ethanol. In gluconeogenesis, pyruvate is converted to glucose.⁴ Other metabolic fates of pyruvic acid include conversion to alanine by transamination and to oxaloacetate by carboxylation.⁵

Pyruvic acid is utilized as a component in culture broths and media.^{6,7} The use of sodium pyruvate in Wallen fermentation medium to enhance the conversion of oleic acid to 10-ketostearic acid by *Bacillus sphaericus* has been described.⁸ A protocol that uses sodium pyruvate to establish stable transfection of human B cell lines has been published.⁹ Pyruvic acid polymerizes and decomposes on standing unless pure and kept in air-tight containers. In aqueous solution, 54-71% of pyruvic acid is in hydrated (gem-diol) form, but the anion is only 3-5% hydrated. There is no enol form present in solution.³

Precautions and Disclaimer

For Laboratory Use Only. Not for drug, household or other uses.

References

- 1. The Merck Index, 12th ed., Entry #8205.
- Lange's Handbook of Chemistry, 13th ed., Dean, J.A., Ed., McGraw-Hill (New York, NY: 1985), pp. 5-52.
- Data for Biochemical Research, 3rd ed., Dawson, R. M. C., et al., Oxford University Press (New York, NY: 1986), pp. 50-51.
- 4. Biochemistry, 3rd ed., Stryer, L., W. H. Freeman (New York, NY: 1988), pp. 349-394.
- 5. Textbook of Biochemistry with Clinical Correlations, Devlin, T. M., ed., Wiley-Liss (New York, NY: 1992), p. 248.
- Geshi, M., et al., Effects of sodium pyruvate in nonserum maturation medium on maturation, fertilization, and subsequent development of bovine oocytes with or without cumulus cells. Biol. Reprod., 63(6), 1730-1734 (2000).
- Jiang, X., and Doyle, M. P., Growth supplements for *Helicobacter pylori*. J. Clin. Microbiol., **38(5)**, 1984-1987 (2000).
- Kuo, T. M., et al., Conversion of fatty acids by Bacillus sphaericus -like organisms. Curr. Microbiol., 45(4), 265-271 (2002).
- Brielmeier, M., et al., Improving stable transfection efficiency: antioxidants dramatically improve the outgrowth of clones under dominant marker selection. Nucleic Acids Res., 26(9), 2082-2085 (1998).

MWM/NSB 11/03

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.