Sigma-Aldrich®

Instructions for Use

Periodic Acid-Schiff (PAS) Staining System

Procedure No. 395

MERCK

Intended Use

The PAS staining system offers standard and microwave procedures for the histochemical demonstration of lymphocytes and mucopolysaccharides. The PAS reaction in tissue sections is useful for the demonstration of mucopolysaccharides. The Diastase (α -amylase) digestion procedure, followed by the PAS stain is useful as an aid in the diagnosis of glycogen storage diseases. Periodic Acid Schiff (PAS) reagents are for "In Vitro Diagnostic Use". For professional use only. The data obtained from this manual, qualitative procedure identifies lymphocytes and mucopolysaccharides in blood or bone marrow films as well as Formalin-Fixed Paraffin Embedded tissue from human specimens. This data when reviewed in conjunction with other diagnostic tests and information may be used as an aid to diagnosis of fungal infections and glycogen storage diseases. Histological visualization of lymphocytes and mucopolysaccharides is a unique technique currently commonly used in medicine.

When treated with periodic acid, glycols are oxidized to aldehydes. After reaction with Schiff's reagent, the released fuchsin stains the cellular components containing the oxidizable compounds.³ This reaction can be performed on blood or bone marrow films, tissue touch preparations or tissue sections.^{1,2}

Diastase (α -amylase) digestion may be employed as an aid in the diagnosis of glycogen storage disease. Diastase hydrolyzes starch, glycogen and degradation products originating in these polysaccharides present in tissue. The resultant by-products of the digestion process are rinsed away prior to PAS staining.⁴

A technique for rapid PAS staining using microwave ovens is included.⁵⁻⁷

Reagents

Periodic Acid Solution (Cat. No. 3951-100ML and 395132-1L)

Periodic acid, 1 g/dL

Schiff's Reagent (Cat. No. 3952-50ML and 3952016-500ML)

Pararosaniline HCl, 1%, and sodium metabisulfite, 4%, in hydrochloric acid, 0.25 mol/L. Danger. Causes severe skin burns and eye damage. May cause cancer.

Hematoxylin Solution, Gill No. 3 (Cat. No. GHS3-100ML)

Certified hematoxylin, 6 g/L, C.I. 75290 sodium iodate, 0.6 g/L, aluminum sulfate, 52.8 g/L and stabilizer

Special Materials Required but Not Provided

- Formaldehyde solution, 37%
- Reagent Alcohol
- Positive PAS control slides, such as PAS TISSUE-TROL™ (Cat. No. TTR009-25EA) should be included in each run
- Xylene or xylene substitute

Microwave Procedures Only

- Microwave Oven
- Coplin jar with vented lids
- Scott's Tap Water Substitute Concentrate (Cat. No. S5134-6x100ML)
- α -Amylase (for Diastase Extraction Procedure only) (Cat. No. A3176)

Storage and Stability

Store Periodic Acid Solution, and Schiff's Reagent in refrigerator (2–8°C). Store Hematoxylin Solution, Gill No. 3 at room temperature (18–26°C). Reagents are stable until expiration date shown on labels. Formation of a precipitate in Schiff's Reagent does not affect performance.

Deterioration

Discard Hematoxylin Solution, Gill No. 3, if solution turns brown (over oxidized from air) or purple (loss of acidity) or when time required for suitable staining exceeds the time recommended in the procedure by more than 5 minutes.

Preparation

Periodic Acid Solution, Schiff's Reagent and Hematoxylin Solution, Gill No. 3 are supplied ready to use. Warm reagents to room temperature (18–26°C) before use.

Filter Schiff's Reagent if precipitate is present. NOTE: A slight pink coloration due to leaching of color from slides does not interfere with test.

Scott's Tap Water Substitute working solution is prepared by diluting 1 part of Scott's Tap Water Substitute Concentrate with 9 parts deionized water.

Formalin-Ethanol Fixative Solution is prepared by mixing 5 mL of formaldehyde with 45 mL of 95% ethanol (Reagent Alcohol). Prepare fresh daily and keep tightly capped.

Precautions

The IVDs included in this kit are intended for in vitro diagnostic use in a clinical laboratory environment. These IVDs are for professional use by qualified personnel only. Sigma-Aldrich IVDs may be operated by laboratory personnel who are trained to handle human specimens that can be infectious, use microscopes and other laboratory equipment and have color perception and visual acuity to distinguish colors and other objects under a microscope.

Normal precautions exercised in handling laboratory reagents should be followed. Dispose of waste observing all local, state, provincial or national regulations.

PAS TISSUE-TROL™ control slides are paraffin embedded human tissue containing PAS and should be considered potentially infectious.

Procedure

Specimen Collection

No known test method can offer complete assurance that blood samples or tissue will not transmit infection. Therefore, all blood derivatives or tissue specimens should be considered potentially infectious.

Freshly prepared whole, EDTA or heparinized blood or bone marrow films are used. Fix as soon as possible.³

For polysaccharides, tissue fixed in 10% neutral buffered formalin, Zenker's or Bouin's may be used.¹ It should be noted that some carbohydrates are water soluble.² For the demonstration of glycogen, Carnoy's fluid, Gendre's fluid or acid alcoholic formalin are recommended.¹ Time required for diastase extraction may be prolonged when tissue is fixed in a picric acid containing fixative.¹ Cut tissue sections at 5 microns.

Notes

If a microwave oven is used, please see the Owner's Manual for instructions.

Blood films prepared from clinically healthy individuals may be included for control purposes. Polymorphonuclear leukocytes will show an intense red cytoplasmic stain. Tissue sections known to be PAS positive and/or contain glycogen should be included each time a stain sequence is performed.

Procedure

Blood, Bone Marrow, or Tissue Touch Preparations

Standard Procedure

- Fix air dried films for **1 minute** at **room temperature** in Formalin-Ethanol Fixative Solution.
- Rinse slides **1 minute** in slowly running tap water.
- Immerse slides in Periodic Acid Solution for **5 minutes** at **room temperature**.
- Rinse slides in several changes of distilled water.
- Immerse slides in Schiff's Reagent for **15 minutes** at room temperature.
NOTE: Immediately after use, cap Schiff's Reagent and return to refrigerator (2–8°C).
- Wash slides in running tap water for **5 minutes**.
- Counterstain slides in Hematoxylin Solution, Gill No. 3, for **90 seconds**.
- Rinse slides in running tap water for **15–30 seconds**, air dry and examine microscopically under oil immersion (900x) lens. Slides may be mounted in toluene or xylene based mounting media.

Microwave Procedure

- Fix air dried films at **room temperature** for **1 minute** in Formalin-Ethanol Fixative Solution.
- Rinse slides for **1 minute** in slowly running tap water.
- Place slides in **40 mL** of Periodic Acid Solution contained in a plastic Coplin jar.
- Microwave on **800 watts** for **10 seconds**.
- Rinse well in several changes of deionized water.
- Place slides in **40 mL** Schiff's Reagent contained in a plastic Coplin jar.
- Microwave at **800 watts** for **15 seconds**. Mix solution with a beral pipet or applicator stick and let incubate for **1 minute**.
- Rinse in warm, gently running tap water for **5 minutes**.
- Place slides in **40 mL** of Hematoxylin Solution Gill No. 3 contained in a plastic Coplin jar.
- Microwave on **800 watts** for **10 seconds**.
- Rinse in running tap water for **1–2 minutes**, then blue in Scott's Tap Water Substitute working solution at **room temperature**.
- Rinse in running tap water. Air dry.
- Slides may be mounted in toluene or xylene based mounting media.

Tissue Sections

Standard Procedure

- Deparaffinize and hydrate sections to deionized water.
- Immerse slides in Periodic Acid Solution for **5 minutes** at **room temperature** (18–26°C).
- Rinse slide in several changes of distilled water.
- Immerse slides in Schiff's Reagent for **15 minutes** at **room temperature** (18–26°C).
- Wash slides in running tap water for **5 minutes**.
- Counterstain slides in Hematoxylin Solution, Gill No. 3, for **90 seconds**.
- Rinse slides in running tap water.
- Dehydrate, clear and mount sections in toluene or xylene based mounting media.

Microwave Procedure

- Deparaffinize and hydrate to deionized water.
- Place slides in **40 mL** of Periodic Acid Solution contained in a plastic Coplin jar. Loosely cover jar with lid, or use lids with holes drilled into them.
- Microwave at **800 watts** for **10 seconds**.
- Rinse well in several changes of deionized water.
- Place slides in **40 mL** Schiff's Reagent contained in a plastic Coplin jar.
- Microwave at **800 watts** for **15 seconds**. Mix solution with a beral pipet or applicator stick and let incubate for **1 minute**.
- Rinse in warm, gently running tap water for **5 minutes**.
- Place slides in Hematoxylin Solution, Gill No. 3, contained in a plastic Coplin jar.
- Microwave on **800 watts** for **10 seconds**.
- Rinse in running tap water for **1–2 minutes**, then blue in Scott's Tap Water Substitute working solution at **room temperature**. Rinse in running tap water. Dehydrate, clear and mount.

Microwave Procedure for Diastase (α -Amylase) Digestion

- Use duplicate test slides. Label one for digestion with diastase and one for PAS staining only. NOTE: Slides coated with a tissue adhesive are recommended. Do not collodionize sections when doing diastase digestion.¹
- Deparaffinize and hydrate slides to deionized water.
- Prepare Diastase (α -amylase) Working Solution by dissolving 0.2 g α -Amylase, (Cat. No. A3176), to 40 mL of deionized water. Mix well and place in plastic Coplin jar. Prepare just prior to use.
- Microwave at **600 watts** for **25 seconds**.
- Remove slides from Coplin jar and rinse digested slide in running tap water for **5 minutes**.
- Using both the digested and the undigested slides, proceed with the Microwave Procedure for Tissue, step 2.

Performance Characteristics

PAS positive substances stain pink to red and nuclei are blue. A Diastase (α-Amylase) Extraction slide will have no visible PAS staining of glycogen when compared to the undigested glycogen positive control slide.

If observed results vary from expected results, please contact Sigma-Aldrich Technical Service for assistance.

Analytical Performance Characteristics

The analytical performance results for the given tests conducted on all target structures, confirm 100% sensitivity, specificity, and repeatability.

Cat. No	Product Description	Target	Intra-assay Specificity	Intra-assay Sensitivity	Inter-assay Specificity	Inter-assay Sensitivity
3951	Periodic Acid Solution	Monocytes	3 of 3	3 of 3	3 of 3	3 of 3
		Granulocytes	3 of 3	3 of 3	3 of 3	3 of 3
		Glycogen	3 of 3	3 of 3	3 of 3	3 of 3
		Basement Membrane	3 of 3	3 of 3	3 of 3	3 of 3
3952	Schiff's Reagent	Monocytes	3 of 3	3 of 3	3 of 3	3 of 3
		Granulocytes	3 of 3	3 of 3	3 of 3	3 of 3
		Glycogen	3 of 3	3 of 3	3 of 3	3 of 3
		Basement Membrane	3 of 3	3 of 3	3 of 3	3 of 3
GHS3-50ML	Hematoxylin Solution, Gill No. 3	Nuclei	3 of 3	3 of 3	3 of 3	3 of 3

Warnings and Hazards

Refer to Safety Data Sheet and product labeling for any updated risk, hazard or safety information.

395B:

- H290: May be corrosive to metals.
- H302: Harmful if swallowed.
- H314: Causes severe skin burns and eye damage.
- H350: May cause cancer.
- H373: May cause damage to organs (Kidney) through prolonged or repeated exposure if swallowed.
- H412: Harmful to aquatic life with long lasting effects.
- P273: Avoid release to the environment.
- P280: Wear protective gloves/ protective clothing/ eye protection/ face protection.
- P301 + P312: IF SWALLOWED: Call a POISON CENTER/ doctor if you feel unwell.
- P303 + P361 + P353: IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water.
- P304 + P340 + P310: IF INHALED: Remove person to fresh air and keep comfortable for breathing. Immediately call a POISON CENTER/ doctor.
- P305 + P351 + P338: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

If during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.

Symbol Definitions

Symbols as defined in EN ISO 15223-1:2021

	Manufacturer		Catalogue Number
	Consult Instructions for Use		Batch Code
	Authorized Representative in the European Community/ European Union		European Union Declaration of Conformity (defined in IVDR 2017/746)
	Use-by Date		In vitro diagnostic medical device
	Temperature Limit		Caution
	Date of Manufacture		Importer
	Indicates the Authorised Representative in Switzerland		

References

1. Sheehan DC, Hrapchak BB: Theory and Practice Histotechnology, 2nd ed. CV Mosby, St. Louis, (MO), pp 52, 164–167, 1980

2. Culling CFA, Allison RT, Barr WT: Cellular Pathology Technique, 4th ed. Butterworths, pp 216–220, 1985

3. Davey FR, Nelson DA: Periodic Acid Schiff (PAS) Stain. IN Hematology, 2nd ed. WJ Williams, E Buetler, AJ Erslev, RW Rundles, McGraw-Hill, New York, pp 1630–1632, 1977

4. Thompson SW: Selected His tochemi cal and Histopathological Methods, CC Thomas, Springfield, (IL), pp 520–539, 1966

5. Leong AS-Y, Milios J: Rapid immunoperoxidase staining of lymphocyte antigens using microwave irradiation. J Pathol 148:183, 1986

6. Brinn NT: Rapid metallic histologic staining using the microwave oven. J Histotechnol 6:125, 1983

7. Valle S: Special stains in the microwave oven. J. Histotechnol 9:237, 1986

Contact Information

To place an order, please visit our web site at [SigmaAldrich.com](https://www.sigmaaldrich.com). For Technical Service, please visit the tech service page on our web site at [SigmaAldrich.com/techservice](https://www.sigmaaldrich.com/techservice).

Revision History

Rev. 4.0	2016
Rev. 5.0	2022
Rev. 6.0	2022
	Transferred to new template with current branding. Specified for professional use in intended use and precautions. Moved aid to diagnosis statement to intended use. Revised intended use to align with IVDR guidelines. Updated Material Safety Data Sheet to Safety Data Sheet. Updated contact information. Removed instruction to follow CLSI for specimen collection. Removed EN 980 and changed to EN ISO 15223-1:2021 for symbols. Added adverse event contact information. Added Warnings and Hazards.
Rev. 7.0	2024
	Updated naming conventions. Revised Inteded Use statement. Added CH-REP and importer information

Sigma-Aldrich, Inc.,
3050 Spruce Street,
St. Louis, MO 63103 USA
an affiliate of Merck KGaA,
Darmstadt, Germany
+1(314) 771-5765

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The Initial M, TISSUE-TROL, and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

© 2024 Merck KGaA, Darmstadt, Germany and/or its affiliates. All rights reserved.

The life science business of Merck operates as MilliporeSigma in the U.S. and Canada.

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Bibliografia

1. Sheehan DC, Hrapchak BB: Theory and Practice Histotechnology, 2nd ed. CV Mosby, St. Louis, (MO), pp 52, 164–167, 1980
2. Culling CFA, Allison RT, Barr WT: Cellular Pathology Technique, 4th ed. Butterworths, pp 216–220, 1985
3. Davey FR, Nelson DA: Periodic Acid Schiff (PAS) Stain. IN Hematology, 2nd ed. WJ Williams, E Buetler, AJ Erslev, RW Rundles, McGraw-Hill, New York, pp 1630–1632, 1977
4. Thompson SW: Selected Histochemical and Histopathological Methods, CC Thomas, Springfield, (IL), pp 520–539, 1966
5. Leong AS-Y, Milios J: Rapid immunoperoxidase staining of lymphocyte antigens using microwave irradiation. J Pathol 148:183, 1986
6. Brinn NT: Rapid metallic histologic staining using the microwave oven. J Histotechnol 6:125, 1983
7. Valle S: Special stains in the microwave oven. J. Histotechnol 9:237, 1986

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Referencias

- 1. Sheehan DC, Hrapchak BB: Theory and Practice Histotechnology, 2nd ed. CV Mosby, St. Louis, (MO), pp 52, 164–167, 1980
- 2. Culling CFA, Allison RT, Barr WT: Cellular Pathology Technique, 4th ed. Butterworths, pp 216–220, 1985
- 3. Davey FR, Nelson DA: Periodic Acid Schiff (PAS) Stain. IN Hematology, 2nd ed. WJ Williams, E Buetler, AJ Erslev, RW Rundles, McGraw-Hill, New York, pp 1630–1632, 1977
- 4. Thompson SW: Selected His tochemi cal and Histopathological Methods, CC Thomas, Springfield, (IL), pp 520–539, 1966
- 5. Leong AS-Y, Milios J: Rapid immunoperoxidase staining of lymphocyte antigens using microwave irradiation. J Pathol 148:183, 1986
- 6. Brinn NT: Rapid metallic histologic staining using the microwave oven. J Histotechnol 6:125, 1983
- 7. Valle S: Special stains in the microwave oven. J. Histotechnol 9:237, 1986

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Referencer

- Sheehan DC, Hrapchak BB: Theory and Practice of Histotechnology, 2nd ed. CV Mosby, St. Louis, (MO), pp 52, 164–167, 1980
- Culling CFA, Allison RT, Barr WT: Cellular Pathology Technique, 4th ed. Butterworths, pp 216–220, 1985
- Davey FR, Nelson DA: Periodic Acid Schiff (PAS) Stain. IN Hematology, 2nd ed. WJ Williams, E Buetler, AJ Erslev, RW Rundles, McGraw-Hill, New York, pp 1630–1632, 1977
- Thompson SW: Selected Histochemical and Histopathological Methods, CC Thomas, Springfield, (IL), pp 520–539, 1966

- Leong AS-Y, Milios J: Rapid immunoperoxidase staining of lymphocyte antigens using microwave irradiation. J Pathol 148:183, 1986
- Brinn NT: Rapid metallic histologic staining using the microwave oven. J Histotechnol 6:125, 1983
- Valle S: Special stains in the microwave oven. J. Histotechnol 9:237, 1986

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Referências

- 1. Sheehan DC, Hrapchak BB: Theory and Practice Histotechnology, 2nd ed. CV Mosby, St. Louis, (MO), pp 52, 164–167, 1980
- 2. Culling CFA, Allison RT, Barr WT: Cellular Pathology Technique, 4th ed. Butterworths, pp 216–220, 1985
- 3. Davey FR, Nelson DA: Periodic Acid Schiff (PAS) Stain. IN Hematology, 2nd ed. WJ Williams, E Buetler, AJ Erslev, RW Rundles, McGraw-Hill, New York, pp 1630–1632, 1977
- 4. Thompson SW: Selected Histological and Histopathological Methods, CC Thomas, Springfield, (IL), pp 520–539, 1966
- 5. Leong AS-Y, Milios J: Rapid immunoperoxidase staining of lymphocyte antigens using microwave irradiation. J Pathol 148:183, 1986
- 6. Brinn NT: Rapid metallic histologic staining using the microwave oven. J Histotechnol 6:125, 1983
- 7. Valle S: Special stains in the microwave oven. J. Histotechnol 9:237, 1986

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1. Sheehan DC, Hrapchak BB: Theory and Practice Histotechnology, 2nd ed. CV Mosby, St. Louis, (MO), pp 52, 164–167, 1980
2. Culling CFA, Allison RT, Barr WT: Cellular Pathology Technique, 4th ed. Butterworths, pp 216–220, 1985
3. Davey FR, Nelson DA: Periodic Acid Schiff (PAS) Stain. IN Hematology, 2nd ed. WJ Williams, E Buetler, AJ Erslev, RW Rundles, McGraw-Hill, New York, pp 1630–1632, 1977

4. Thompson SW: Selected Histological and Histopathological Methods, CC Thomas, Springfield, (IL), pp 520–539, 1966
5. Leong AS-Y, Milios J: Rapid immunoperoxidase staining of lymphocyte antigens using microwave irradiation. J Pathol 148:183, 1986
6. Brinn NT: Rapid metallic histologic staining using the microwave oven. J Histotechnol 6:125, 1983
7. Valle S: Special stains in the microwave oven. J. Histotechnol 9:237, 1986

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Referanser

- Sheehan DC, Hrapchak BB: Theory and Practice Histotechnology, 2nd ed. CV Mosby, St. Louis, (MO), pp 52, 164–167, 1980
- Culling CFA, Allison RT, Barr WT: Cellular Pathology Technique, 4th ed. Butterworths, pp 216–220, 1985
- Davey FR, Nelson DA: Periodic Acid Schiff (PAS) Stain. IN Hematology, 2nd ed. WJ Williams, E Buetler, AJ Erslev, RW Rundles, McGraw-Hill, New York, pp 1630–1632, 1977
- Thompson SW: Selected Histochemical and Histopathological Methods, CC Thomas, Springfield, (IL), pp 520–539, 1966
- Leong AS-Y, Milios J: Rapid immunoperoxidase staining of lymphocyte antigens using microwave irradiation. J Pathol 148:183, 1986
- Brinn NT: Rapid metallic histologic staining using the microwave oven. J Histotechnol 6:125, 1983
- Valle S: Special stains in the microwave oven. J. Histotechnol 9:237, 1986

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Lähteet

- Sheehan DC, Hrapchak BB: Theory and Practice Histotechnology, 2nd ed. CV Mosby, St. Louis, (MO), pp 52, 164–167, 1980
- Culling CFA, Allison RT, Barr WT: Cellular Pathology Technique, 4th ed. Butterworths, pp 216–220, 1985
- Davey FR, Nelson DA: Periodic Acid Schiff (PAS) Stain. IN Hematology, 2nd ed. WJ Williams, E Buetler, AJ Erslev, RW Rundles, McGraw-Hill, New York, pp 1630–1632, 1977
- Thompson SW: Selected His tochemi cal and Histopathological Methods, CC Thomas, Springfield, (IL), pp 520–539, 1966
- Leong AS-Y, Milios J: Rapid immunoperoxidase staining of lymphocyte antigens using microwave irradiation. J Pathol 148:183, 1986
- Brinn NT: Rapid metallic histologic staining using the microwave oven. J Histotechnol 6:125, 1983
- Valle S: Special stains in the microwave oven. J. Histotechnol 9:237, 1986

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Referenties

1. Sheehan DC, Hrapchak BB: Theory and Practice Histotechnology, 2nd ed. CV Mosby, St. Louis, (MO), pp 52, 164–167, 1980
2. Culling CFA, Allison RT, Barr WT: Cellular Pathology Technique, 4th ed. Butterworths, pp 216–220, 1985
3. Davey FR, Nelson DA: Periodic Acid Schiff (PAS) Stain. IN Hematology, 2nd ed. WJ Williams, E Buetler, AJ Erslev, RW Rundles, McGraw-Hill, New York, pp 1630–1632, 1977
4. Thompson SW: Selected Histological and Histopathological Methods, CC Thomas, Springfield, (IL), pp 520–539, 1966
5. Leong AS-Y, Milios J: Rapid immunoperoxidase staining of lymphocyte antigens using microwave irradiation. J Pathol 148:183, 1986
6. Brinn NT: Rapid metallic histologic staining using the microwave oven. J Histotechnol 6:125, 1983
7. Valle S: Special stains in the microwave oven. J. Histotechnol 9:237, 1986

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Šaltiniai

1. Sheehan DC, Hrapchak BB: Theory and Practice Histotechnology, 2nd ed. CV Mosby, St. Louis, (MO), pp 52, 164–167, 1980
2. Culling CFA, Allison RT, Barr WT: Cellular Pathology Technique, 4th ed. Butterworths, pp 216–220, 1985
3. Davey FR, Nelson DA: Periodic Acid Schiff (PAS) Stain. IN Hematology, 2nd ed. WJ Williams, E Buetler, AJ Erslev, RW Rundles, McGraw-Hill, New York, pp 1630–1632, 1977
4. Thompson SW: Selected His tochemi cal and Histopathological Methods, CC Thomas, Springfield, (IL), pp 520–539, 1966
5. Leong AS-Y, Milios J: Rapid immunoperoxidase staining of lymphocyte antigens using microwave irradiation. J Pathol 148:183, 1986
6. Brinn NT: Rapid metallic histologic staining using the microwave oven. J Histotechnol 6:125, 1983
7. Valle S: Special stains in the microwave oven. J. Histotechnol 9:237, 1986

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Odkazy

- Sheehan DC, Hrapchak BB: Theory and Practice Histotechnology, 2nd ed. CV Mosby, St. Louis, (MO), pp 52, 164–167, 1980
- Culling CFA, Allison RT, Barr WT: Cellular Pathology Technique, 4th ed. Butterworths, pp 216–220, 1985
- Davey FR, Nelson DA: Periodic Acid Schiff (PAS) Stain. IN Hematology, 2nd ed. WJ Williams, E Buetler, AJ Erslev, RW Rundles, McGraw-Hill, New York, pp 1630–1632, 1977
- Thompson SW: Selected Histological and Histopathological Methods, CC Thomas, Springfield, (IL), pp 520–539, 1966

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MDSS GmbH
MDSS CH GmbH
Chemie GmbH