MultiScreen₃₈₄-PCR Filter Plate

Automation-friendly 384-well plate for high-throughput purification of PCR* products



Use MultiScreen₃₈₄-PCR for one-step cleanup prior to sequencing or microarray preparation.

High Throughput PCR Purification

 $MultiScreen_{384}$ -PCR filter plates use size selective membranes to separate PCR products from contaminants. Designed with the needs of the high throughput laboratory in mind, MultiScreen₃₈₄-PCR meets SBS standards and is ideally suited to "walk away" automation.

Quick and Easy Protocol MultiScreen₃₈₄-PCR minimizes the procedural complexity involved in traditional precipitation or bind and elute methods. Use of the ${\sf MultiScreen}_{{\tt 384}}\text{-PCR filter plate elimi-}$ nates centrifugation, time consuming incubations, ethanol-drying steps and manifold disassembly routines. Resuspension in your choice of buffer or water adds flexibility and facilitates integration with downstream steps.

The protocol is simple. Filter until wells are empty and then resuspend PCR products free of dNTP's, salt, and 90% free of primers in your choice of buffer.

Simply Filter, Reconstitute, and Recover



1. Introduce samples into wells (up to 100 µL).



2. Filter using MultiScreen₃₈₄ vacuum manifold for 10 to 20 minutes or until wells are empty. DNA is retained on the surface of the membrane while smaller contami nants are filtered to waste



3. Add resuspension buffer to each well (typically 10 to 100 µl) and recover purified DNA from the membrane surface by mixing or shaking. Retrieve purified samples by aspiration.



Design and Protocol Saves Time and Money

The 20 minute protocol does not require special reagents that can contribute to increased cost, process time and automation complexity.

Moreover, there is no risk of exonuclease degradation of the PCR products and microarray preparation is facilitated by the removal of nucleotides and salt.

Features:

- Automation friendly
- No practical capacity limits
- Resuspension at a desired concentration and in an aqueous buffer of choice allows for greater flexibility when preparing samples for downstream applications
- $-10-100~\mu L$ working volume better for processing smaller samples

Performance




Figure 1. Representative electropherogram from a 1159 bp PCR product purified using the MultiScreen $_{\it 384}$ -PCR filter plate.

Specifications

Maximum Well Volume: 130 μL Working Well Volume: 100 μL Non-recoverable volume: < 2 μL Distance from bottom of plate to membrane surface: 3 mm

Materials of construction

Plastic Housing: Clear Polystyrene Membrane: Proprietary

Automation Platform Recommendations

MultiScreen₃₈₄-PCR filter plates perform best with 96 or 384 tip liquid handlers (e.g., Beckman's Multimek) capable of pipetting ≥ 25 µL/second or with plate shakers with 0.12 inch orbit @ 1100 rpm. Other liquid handlers or shakers operating at slower speeds are equally effective, but resuspension times are longer.

Contact Millipore Technical Services or consult with one of the Millipore Application Scientists for help with your automation needs (1-800-MILLIPORE).

PCR Product Recovery




Figure 2. MultiScreen₃₈₄-PCR filter plate's average recovery determined using a fluorescence-based assay (n=96). Good recoveries are observed for PCR products >300 bp with salts, dNTPs and >90% of primers being removed in a single step.

Multiple Resuspension Options





Figure 3. 20 μL PCR products were purified on the MultiScreen₃₈₄-PCR plate, resuspended in 20 μL TE buffer as indicated and percent recovery was determined by gel electrophoresis, ethidium bromide staining and densitometry. Good recoveries are observed with liquid handlers using disposable or fixed tips. Shaking too was shown to be effective and is noteworthy for laboratories who are not equipped with 96 pipette liquid handling systems.

Chemical Compatibility

Chemical	,	Chemical	
EDTA	Е	Ethanol 100%	Q
Formamide 70%	Χ	Ethanol 70%	Е
Formamide 20%	Е	Ethanol 20%	Е
DMSO 100%	Χ	Sodium Citrate	Е
DMSO 50%	Χ	2-Pyrrolidinone 20%	Q
DMSO 20%	Q		



Ordering Information

Description	Pack Size	Catalogue No.
MultiScreen ₃₈₄ -PCR filter plate	10	S384 PCR 10
MultiScreen ₃₈₄ -PCR filter plate	50	S384 PCR 50
MultiScreen ₃₈₄ Manifold	1	SAVM 384 01

To Place an Order or Receive Technical Assistance

For additional information call your nearest Millipore office: In the U.S. and Canada, call toll-free 1-800-MILLIPORE

(1-800-645-5476) In the U.S., Canada and Puerto Rico, fax orders to 1-800-MILLIFX (1-800-645-5439)

E-mail: tech service@millipore.com



Now you can buy Millipore products online @

www.millipore.com/purecommerce

MILLIPORE

Millipore and MultiScreen are trademarks of Millipore Corporation or an affiliated company. MultiScreen₃₈₄-PCR plate is patent pending.

Multimek is a trademark of Beckman Instruments, Inc. PCR is covered by patents issued to Hoffman-LaRoche, Inc.

Millipore Lit. No. PF1546EN00 Printed in U.S.A. 6/00 00-160

Copyright 2000 Millipore Corporation or an affiliated company. All rights reserved.

Millipore Worldwide

AUSTRALIA

Tel. 1 800 222 111 or (02) 9888 8999 Fax (02) 9878 0788

AUSTRIA, CENTRAL EUROPE, C.I.S., AFRICA, MIDDLE EAST AND GUIF

Tel. (43) 1 877-8926 Fax (43) 1 877-1654

BALTIC COUNTRIES Tel. +358 9 804 5110

Fax +358 9 256 5660

BELGIUM AND LUXEMBOURG

Tel. (02) 726 88 40 Fax (02) 726 98 84

Tel. (011) 548-7011 Fax (011) 548-7923

CANADA

Tel. 1-800-645-5476 Fax 1-800-645-5439

CHINA, PEOPLE'S REPUBLIC OF

Beiiina Tel. (8610) 8519-1250

(8610) 8519-1058 Fax (8610) 8519-1255 Guangzhou:

Tel. (8620) 8755-4021 Fax (8620) 8755-4350

Hong Kong: Tel. (852) 2803-9111 Fax (852) 2513-0313 Shanghai

Tel. (8621) 5306-9100

Fax (8621) 5306-0838

CZECH REPUBLIC Tel. 02-2051 3841

02-2051 3842 Fax 02-2051 4298

DENMARK

Tel. 70 10 00 23 Fax 70 10 13 14

FINLAND

Tel. (09) 804 5110 Fax (09) 256 5660

FRANCE

Tel. (01) 30.12.70.00 Fax (01) 30.12.71.80

GERMANY

Tel. (06196) 494-0 Fax (06196) 43901

HUNGARY

Tel. 1-205 9784 Fax 1-205 9792

INDIA

Tel. (91) 80-839 46 57 Fax (91) 80-839 63 45

ITALY

Vimodrone (Milano): Tel. (02) 25.07.81 Fax (02) 26.50.324 Roma:

Tel. (06) 52.03.600 Fax (06) 52.95.735

JAPAN

Tel. (03) 5442-9711 Fax (03) 5442-. 9736 Analytical 9737 BioProcess 9734 Lab Water

KOREA

Tel. (822) 551-0990 Fax (822) 551-0228

MALAYSIA

Tel. 603-757 1322 Fax 603-757 1711

MEXICO

Tel. (525) 576-9688 Fax (525) 576-8706

THE NETHERLANDS Tel. 076-5022000

Fax 076-5022436 **NORWAY**

Tel. 22 67 82 53

Fax 22 66 04 60

POLAND

Tel. 22-669 12 25 22-663 70 31 Fax 22-663 70 33

PUERTO RICO

Tel. (787) 273-8495 Fax (787) 747-6553

SINGAPORE

Tel. (65) 842 1822 Fax (65) 842 4988

SPAIN AND PORTUGAL

Tel. 917 283 960 Fax 917 292 909 Barcelona: Tel. 934 525 530

Fax 934 516 048

Madrid:

SWEDEN

Tel. 08-628 6960 Fax 08-628 6457

SWITZERLAND

Volketswil-Tel. (01) 908-30-60 Fax (01) 908-30-80 lausanne

Tel. (021) 641 2550 Fax (021) 641 2551

TAIWAN

Taipei: Tel. (886-2) 700-1742 Fax (886-2) 755-3267 Hsin Chu City: Tel. (886-3) 571-0178

Fax (886-3) 572-9520

U.K. AND IRELAND

Tel. (01923) 816375 Fax (01923) 818297

U.S.A. Tel. (781) 533-6000 Fax (781) 533-3110

IN ALL OTHER COUNTRIES

Millipore Intertech (U.S.A.) Tel. +1 (781) 533-8622 Fax +1 (781) 533-8630