

POLY-L-LYSINE COATED SLIDES FOR MICROARRAYS

ProductInformation

Product Code S 1313 Store at Room Temperature

TECHNICAL BULLETIN

Product Description

Poly-L-Lysine Coated Slides for Microarrays are an optimal substrate for printing DNA microarrays. The slides are prepared by coating stringently clean glass slides with poly-L-lysine resulting in a slide surface densely populated with positive charges. The availability of the positively charged amine groups at the surface translates into high DNA binding and strong hybridization signals. The negatively charged phosphate backbone of DNA is electrostatically attached to the positively charged amines. Treatment with ultraviolet irradiation results in a more stable coupling.

Each lot of slides undergoes rigorous quality testing with high standards for surface energy, coating uniformity, autofluorescence, particulation, and array printing/hybridization to ensure performance and lot-tolot consistency.

Slide Dimensions: 75 mm X 25 mm X 1 mm

Precautions and Disclaimer

Poly-L-Lysine Coated Slides for Microarrays are for laboratory use only, not for drug, household or other uses.

Storage/Stability

Store Poly-L-Lysine Coated Slides for Microarrays at room temperature sealed with desiccant in the storage bag provided. Use the slides in a clean environment. Particles can interfere with the printing (arraying) process and can lead to background signal during detection. In addition, slide handling should be minimized. Touching of the printing surface, except by the printer, should be avoided.

Product Code

A 3095

Reagents and Equipment Required but Not Provided

(Sigma product codes have been given where appropriate)

Product Name ArrayHyb[™] LowTemp Microarray Hybridization Buffer

Standard Microarray Spotting M 1435 Solution Humid Chamber H 6644 Wheaton 20 slide staining dish S 6141 22 x 22 mm Coverslips C 9802 E 7148 95% Ethanol Belly-Dancer/Hybridization Water Z36,762-1 Bath (115 V or 220/240 V) Z36,764-8 BioLink BLX Crosslinker Z37.537-3 Z37.538-1

Wash Solutions

Heat Block

(115 V or 230 V)

Use ready-to-use wash solutions from the Microarray Hybridization Wash Pack (Product Code M 2185) or prepare solutions according to the following table:

		Wash Buffers		
		1	2	3
	0.5% SDS	1X SSC,	0.2X	0.05X
	Solution	0.03% SDS	SSC	SSC
10% SDS	50 ml	3 ml		
(L 4522)				
20X SSC		50 ml	10 ml	2.5 ml
S 6639)				
H ₂ O	to 1L	to 1L	to 1L	to 1L

Procedure

(

(ŀ

A. DNA Printing

Dilute double-stranded DNA (0.1-1 μ g/ μ l in water) or oligonucleotide (100-150 µM in water) samples with an equal volume of the Standard Microarray Spotting Solution. Print slides according to arrayer manufacturer's or standard protocol. Store printed slides desiccated at room temperature until ready for pre-hybridization treatment and hybridization.

B. Post Printing/Pre-Hybridization Treatment

Slides must be pre-treated prior to hybridization after printing DNA or oligonucleotide samples. The following procedure has been optimized for hybridizations using ArrayHyb LowTemp Microarray Hybridization Buffer. Overall signal and consistency may be compromised if slides are pre-treated using other procedures. This procedure denatures spotted, double-stranded DNA, making it available for hybridization, while eluting any non-crosslinked nucleic acid. Multiple slides may be pre-treated simultaneously. Do not write on the slide at this stage; many inks are soluble in ethanol and may contribute to fluorescent backgrounds.

- 1. Heat enough distilled, deionized water to 95-100 °C to cover the slides in a slide staining rack.
- Fill a humid chamber (H 6644) with approximately 50 ml of water and warm to 65 °C in a water bath. Other chambers may require more or less volume. The chamber bottom must be covered with water to facilitate equal hydration of the slides.
- 3. Place the arrayed slides, DNA side down, over the humid chamber. Allow spots to hydrate for approximately 5 seconds.
- 4. Dry the arrayed slide (DNA side up) on a 95-100 °C surface/heat block for approximately 5 seconds.
- 5. UV crosslink the DNA with the slide surface by placing slides, DNA side, up on a rack (an empty humid chamber) and supply 65 mJ of 254 nm irradiation.
- 6. Place slides in a slide rack.
- Incubate the slides for approximately 2 minutes in a slide staining dish filled with 0.5% SDS solution on an orbital shaker. This step removes the unbound nucleic acids from the arrays and helps block nonspecific binding of nucleic acids.
- 8. Rinse the slides by dipping the slide rack several times in a fresh dish of distilled, deionized water.
- Gently plunge the slide rack into 95-100 °C distilled, deionized water for 2 minutes. <u>Note</u>: This step is not necessary for single-stranded oligonucleotides.
- 10. Rinse the slides by plunging the slide rack several times in a staining dish filled with 95+% ethanol.

- 11. Quickly transfer the slides to a centrifuge equipped with a swinging bucket rotor for multiwell plates (place paper towels below the rack to absorb liquid) and spin at 50-100 x g for 5 to 10 minutes.
- 12. Use the treated slides immediately or store in a slide box desiccated at 2-8 °C.
- C. Hybridization
- 1. Equilibrate pre-treated slides to room temperature.
- Prepare a labeled cDNA mixture by combining labeled nucleic acid and blocking agent(s), if desired (see below), in a microcentrifuge tube. It may be necessary to first concentrate the mixture by ethanol precipitation or vacuum centrifugation. Suspend the probe mixture/pellet in an appropriate volume of ArrayHyb LowTemp Microarray Hybridization Buffer for the size cover slips being utilized (2.5 to 3 µl/cm² or 12.5 to 15 µl for 22 x 22 mm cover slips).

	Final
Common Blocking Agents	Concentration
	(µg/ml)
Single-stranded DNA (D 8899)	100
Poly (dA) (P 0887)	400
Yeast tRNA (R 9001)	200
CoT-1 DNA [®]	400

- 3. Heat the ArrayHyb LowTemp/probe mixture at 60 °C for 5 minutes.
- 4. Centrifuge the contents to the bottom of the microcentrifuge tube. Pipette the hybridization solution onto a cover slip.
- 5. Slowly lower the slide, array side down, until surface tension allows the cover slip to be raised with the slide, taking care not to introduce bubbles.
- 6. Incubate the slides from 6 hours to overnight at 50 °C in a humidity-controlled environment. This can be achieved by placing slides in an empty humid chamber placed in a shallow hybridization water bath. Alternatively, hybridization may be carried out in one of several commercially available hybridization chambers immersed in a temperature controlled water bath.

D. Post-Hybridization Wash

<u>Note</u>: Do not allow the slides to dry out at any point during this procedure. This can cause non-specifically bound probe to become irreversibly attached to the slide, resulting in high backgrounds.

- 1. Fill separate wash dishes with Wash Buffers 1 (two dishes), 2, and 3.
- Remove the slides one at a time from the humidity controlled incubator. Immerse the slides into the Wash Buffer 1 dish and gently remove the cover slips. Place the slides into a slide rack in the second Wash Buffer 1 dish. Incubate with gentle mixing on an orbital shaker for 5 minutes at room temperature.
- Rinse the slides by plunging the slide rack into fresh distilled, deionized water dish several times. Remove excess liquid from rack and transfer to Wash Buffer 2. Incubate for 5 minutes at room temperature with gentle mixing on an orbital shaker.
- Rinse the slides by plunging the slide rack into fresh distilled, deionized water dish several times. Remove excess liquid from rack and transfer to Wash Buffer 3. Incubate for 5 minutes at room temperature with gentle mixing on an orbital shaker.
- After the third washing place the slide rack in a centrifuge equipped with a swinging bucket rotor for multiwell plates. Centrifuge at 50-100 x g for 5 to 10 minutes. Immediately remove the slides from the centrifuge and store in a light-proof slide box at room temperature.
- 6. Slides are now ready for scanning.

Related Products

Product Name	Product Code
ArrayHyb™ Hybridization Buffer	A 7718
ArrayHyb™ LowTemp Microarray	A 3095
Hybridization Buffer	
GenElute [™] Mammalian Total RNA Kits	RTN-10
	RTN-70
	RTN-350
GenElute [™] mRNA from Total RNA Kits	MRN-10
	MRN-70
GenElute™ PCR Purification Kit	GEN-PCR
SigmaSpin™ Post-Reaction Clean-up Columns	S 5059
Microarray Hybridization Wash Pack	M 2185
Hybridization Water Bath	Z36,765-6
(115V or 220V)	Z36,766-4
BioLink BLX UV Crosslinker	Z37,537-3
(115V or 230V)	Z37,538-1
Belly-Dancer Orbital Shaker	Z36,760-5
(115V or 220/240V)	Z36,761-3
Poly (dA) _n	P 0887
Yeast tRNA	R 5636

References

- Schena, M., *et al.* Parallel human genomic analysis: microarray-based expression monitoring of 1000 genes. Proc. Natl. Acad. Sci. USA, **93**, 10614-10619 (1996).
- 2. Schena, M., *et al.* Quantitative monitoring of gene expression patterns with a complementary DNA microarray. Science, **270**, 467-470 (1995).
- Schena, M. (Ed.), <u>Microarray Biochip Technology</u>, Eaton Publishing, Natick, MA, 2000, Product No. M 4309.
- Schena, M. (Ed), <u>DNA Microarrays, A Practical</u> <u>Approach</u>, Oxford University Press, Oxford, England, 1999, Product No. D 6187.



Problem	Cause	Solution
Low Signal	Inefficient binding of nucleic acids to the slide during printing	DNA must be cleaned properly prior to spotting. Silica matrix columns, such as the GenElute PCR DNA Purification Kit, are recommended for purification of the DNA.
	Hybridization temperature is too high.	Decrease hybridization temperature (typically by 5-10 °C) to achieve optimal hybridization signal.
	Probe was not labeled efficiently.	Check the probe for labeling efficiency. If poor label incorporation is observed, remake the probe. Protect labeled probes from exposure to light.
Background Fluorescence	Drying of hybridization solution at the edges of the cover slip.	Be sure to hybridize the slides in a humidified chamber. The use of humid chambers will prevent this problem.
	Probe was allowed to dry to the slide during transfer to Wash Buffer 1.	Use extra care to prevent any drying of the probe solution on the slides.
	Dust has accumulated on the slide.	Dust particles will show up as isolated spots of very high fluorescence. Protect slides from general lab air as much as possible. Some dust may be removed by the use of compressed air to "dust" the slide prior to scanning.
Non-specific Hybridization Signal	Hybridization temperature is too low.	Increase hybridization temperature (typically by 5-10 °C) to eliminate non-specific hybridization signal.
Irregular spot morphology	Poor printing of DNA	Dissolve the DNA in Standard Microarray Spotting Solution or another validated printing solution such as 3X SSC.
	Deformed spotting pins	Replace spotting pins.
Spots appear as streaks or comets.	Recommended pre-treatment procedure was not followed.	When concentrated DNA (0.5 mg/ml) is spotted on the slide, only a small percentage becomes bound to the surface. The remaining unbound DNA must be washed away. This is accomplished in the recommended pre-treatment procedure by washing the slides in 0.5% SDS solution. This step removes unbound DNA from the spots and prevents the DNA
		by washing the slides in 0.5% SDS solution. This step removes unbound DNA from the spots and prevents the from binding to the surrounding slide surfaces.

BW 8/01

Sigma brand products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.