

3050 Spruce Street, St. Louis, MO 63103 USA
Tel: (800) 521-8956 (314) 771-5765 Fax: (800) 325-5052 (314) 771-5757
email: techservice@sial.com sigma-aldrich.com

Product Information

Western Blot Procedure

Preparation Instructions

 $1\times$ PBS - 10 mM NaH $_2PO_4\cdot H_2O$ and 130 mM NaCl, pH 7

 $1 \times PBST - 1 \times Phosphate$ buffer saline with 0.2% (v/v) TWEEN® 20

Blocking Buffer – Add 5 g of non-fat milk to 100 mL of $1 \times PBST$

Procedure

Perform appropriate SDS-PAGE and transfer protein to a PVDF membrane. Keep the membrane at –20 °C no longer than a week, if the membrane will not be used immediately.

- Add an adequate volume of Blocking Buffer and leave the membrane at room temperature for 1 hour or overnight at 2–8 °C.
- Remove the membrane from the Blocking Buffer, add the primary antibody diluted in Blocking Buffer to the membrane, and incubate at 2–8 °C overnight.

- 3. Wash the membrane with $1 \times PBST$ for 10 minutes. Repeat the wash 3 times.
- 4. Add the secondary antibody-HRP conjugate and the membrane at room temperature for 1 hour.
- 5. Wash the membrane $1 \times PBST$ for 10 minutes. Repeat the wash 4 times.
- After washing the membrane, place it in a sealable bag and add sufficient freshly prepared chemiluminescent substrate reagents to coat the entire membrane.
- 7. Take photographs immediately with CCD camera at 5 seconds, 20 seconds, 1.5 minutes, and 5 minutes to obtain proper image.

TWEEN is a registered trademark of Croda International PLC.

JX.LPG.MAM 12/08-1