Chrysamine G and BTA-1: New probes for β-amyloid aggregates

Prod. Nos. C 1115 and B 9934

BTA-1

The pathological signature of Alzheimer's disease is the deposition of β -amyloid protein (A β). Its cleavage products, such as A β 40 and A β 42, form amyloid fibrils and plaques in the brains of affected individuals. Compounds that have affinity for A β have the ability to prevent neurotoxicity by inhibiting aggregation of amyloid fibrils. In addition, these molecules can also serve to quantify amyloid deposits *in vivo* and thus may serve as diagnostic tools [1].

Chrysamine G is a lipophilic analog of **Congo red** (Prod. Nos. **C 6767** and **C 6277**). The compound displays a K_i value of 25.3 nM [2] and labels both high- ($K_d = 0.2 \ \mu\text{M}$; $B_{max} = 1.13 \ \text{moles}$ per mole of A β 40) and low- ($K_d = 39 \ \mu\text{M}$) affinity binding sites for synthetic A β [1]. At concentrations of 0.1 to 1 μ M, Chrysamine G increases the number of surviving neurons following exposure to both A β 25-35 and A β 40 [3].

BTA-1 is a fluorescent **Thioflavin-T** (Prod. No. **T 3516**) derivative. The compound exhibits high affinity for synthetic aggregated A β 40 fibrils, displaying a K $_{\rm i}$ value of 11 nM vs. 580 nM for Thioflavin-T [4]. In addition, BTA-1 crosses the blood-brain barrier and displays 50-fold higher affinity than Thioflavin-T. It selectively stains cerebral plaques and cerebrovascular amyloid deposits in the brains of PS1/APP transgenic mice, as well as A β fibrils in postmortem brain tissue obtained from Alzheimer's disease patients [4].

Both Chrysamine G and BTA-1 are high affinity probes for $A\beta$ that cross the blood-brain barrier. They will thus serve as useful tools in the study of $A\beta$ -induced neurotoxicity, as well as in the development of *in vivo* imaging tools to detect formation of amyloid fibrils and plaques.

References

- Kang, J. and Han, K., The amide derivatives of chrysamine G bind to the β-amyloid fibril., Bull. Korean Chem. Soc., 22, 1065-1066 (2001).
- Klunk, W.E., et al., Imaging Aβ plaques in living transgenic mice with multiphoton microscopy and methoxy-X04, a systemically administered Congo red derivative., J. Neuropathol. Exp. Neurol., 61, 797-805 (2002).
- 3. Ishii, K., et al., Chrysamine G and its derivative reduce amyloid β -induced neurotoxicity in mice., Neurosci. Lett., 333, 5-8 (2002).
- Mathis, C.A., et al., A lipophilic thioflavin-T derivative for positron emission tomography (PET) imaging of amyloid in brain., *Bioorg. Med. Chem. Lett.*, 12, 295-298 (2002).

VDM 11: Potent and selective anandamide membrane transporter inhibitor

Prod. No. V 3264

Anandamide (Prod. No. **A 0580**) is a putative endogenous agonist at CB_1 cannabinoid and TPRV1 vanilloid receptors [1]. Its physiological effects are terminated by a two-step process that involves reuptake via the anandamide membrane transporter (AMT), followed by hydrolysis catalyzed by the enzyme fatty acid amide hydrolase [2]. A derivative of anandamide, VD M11 has recently been described as a potent and selective AMT inhibitor [3].

VDM 11 inhibits AMT activity, displaying IC_{50} values of ~10 μ M in C6 glioma and RBL-2H3 cells [3]. However, although equipotent to another AMT inhibitor, **AM404** (Prod. No. **A-262**) VDM 11 exhibits negligible agonist activity at the human TRPV1 vanil-loid receptor, as well as CB₁ and CB₂ cannabinoid receptors [3].

VDM 11 is therefore a new and potentially superior pharmacological tool for studying the role of AMT in the physiological termination of the effects of anandamide.

References

- De Petrocellis, L., et al., The activity of anandamide at vanilloid VR1 receptors requires facilitated transport across the cell membrane and is limited by intracellular metabolism., J. Biol. Chem., 276, 12856-12863 (2001).
- Ueda, N., et al., The fatty acid amide hydrolase (FAAH)., Chem. Phys. Lipids, 108, 107-121 (2000).
- De Petrocellis, L., et al., Overlap between the ligand recognition properties
 of the anandamide transporter and the VR1 vanilloid receptor: inhibitors
 of anandamide uptake with negligible capsaicin-like activity., FEBS Lett.,
 483, 52-56 (2000).

URB597: Potent and selective fatty acid amide hydrolase (FAAH) inhibitor

Prod. No. **U 4133**

The enzyme fatty acid amide hydrolase (FAAH) has been shown to hydrolyze endogenous cannabinoids, such as anandamide (Prod. No. A 0580), 2-arachidonylglycerol (Prod. No. A 8973) and oleamide (Prod. No. O 2136). In mice lacking FAAH, it has been

shown that a decrease in FAAH activity leads to reduced sensation of pain and enhanced signaling by endogenous cannabinoids [1]. Blocking FAAH may thus serve to enhance the analgesic effects of endogenous cannabinoids [2].

URB597 has recently been described as a potent and selective FAAH inhibitor, displaying IC $_{50}$ values of 4.6 nM and 0.5 nM in rat cortical neuron membranes and intact neurons, respectively [3]. This compound does not interfere with the binding of anandamide to CB $_1$ and CB $_2$ cannabinoid receptors (IC $_{50}$ > 100 mM), but has been shown to exhibit anxiolytic effects in rats (ID $_{50}$ 0.15 mg/Kg) [3].

URB597 is thus a promising tool with which to explore the development of novel anxiolytic and analgesic drugs.

References

- Cravatt, B.F., et al., Supersensitivity to anandamide and enhanced endogenous cannabinoid signaling in mice lacking fatty acid amide hydrolase., Proc. Natl. Acad. Sci. USA, 98, 9371-9376 (2001).
- Cravatt, B.F., et al., Molecular characterization of an enzyme that degrades neuromodulatory fatty-acid amides., Nature, 384, 83-87 (1996).
- Kathuria, S., et al., Modulation of anxiety through blockage of anandamide hydrolysis., Nat. Med., 1, 76-81 (2003).



Celltransmissions Vol 20, No. 2, 2004