



PRAK Substrate Peptide
Product Number **P 0240**

Synonyms: P38 Regulated/Activated
Protein Kinase Substrate Peptide

Product Description

The 1171 dalton synthetic decapeptide, KKLRTLSVA, is provided as a substrate for p38-Regulated/Activated Protein Kinase or PRAK, active (P 0365). PRAK is a 471 amino acid serine/threonine kinase that is activated in response to cellular stress and proinflammatory cytokines.^{1,2,3} Activity of PRAK is regulated by p38 α and p38 β (members of the MAP kinase superfamily). Activated PRAK will phosphorylate PRAK Substrate Peptide and small heat shock protein 27 (HSP27).¹

Supplied as a frozen solution in water.

Precautions and Disclaimer

For laboratory use only. Not for drug, household or other uses. Please consult Material Safety Data Sheet for handling recommendations.

Storage/Stability

Undiluted solution is stable for one year upon delivery if properly stored at -20 °C. Aliquot to avoid repeated freeze-thaw cycles.

Protein Kinase Assay Procedure

Stock Solutions

1. Dilute PRAK Substrate Peptide to 0.17 $\mu\text{g}/\mu\text{l}$ in Assay Dilution Buffer (ADB) (20 mM MOPS, pH 7.2, 25 mM β -glycerol phosphate, 5 mM EGTA, 1 mM sodium orthovanadate, 1 mM dithiothreitol). A final concentration of 30 μM is used per assay (10 $\mu\text{l}/\text{assay}$).
2. Dilute PRAK, active (P0365) to 20 ng/ μl in ADB. A volume of 10 μl will be used per assay.

Product Information

3. Prepare 10 μl aliquots of a 1mCi/100 μl stock of [γ -³²]ATP (100 $\mu\text{Ci}/\text{vial}$). Before beginning the assay, dilute an aliquot to 1 $\mu\text{Ci}/\mu\text{l}$ with 90 μl of 75 mM MgCl_2 and 500 μM cold ATP.

Assay Procedure:

1. Add 20 μl ADB to a microcentrifuge tube.
2. Add 10 μl of 0.17 $\mu\text{g}/\mu\text{l}$ PRAK Substrate Peptide.
3. Add 10 μl of diluted PRAK, active (200 ng purified enzyme per assay).
4. Add 10 μl of the diluted [γ -³²P]ATP solution.
5. Incubate with agitation for 10 minutes at 30 °C.
6. Transfer a 35 μl aliquot onto the center of a 2 cm by 2 cm P81 paper square.
7. Wash the paper squares three times with 0.75% phosphoric acid.
8. Wash the paper squares once for 3 minutes with acetone.
9. Transfer the paper squares to a scintillation vial and add 5 ml scintillation cocktail.
10. Read in a scintillation counter.

References

1. New, L. *et al.*, PRAK, a novel protein kinase regulated by the p38 MAP Kinase. *EMBO J.*, **17**, 3372-3384 (1998).
2. Raingeaud, J., *et al.*, Pro-inflammatory cytokines and environmental stress cause p38 MAP kinase activation by dual phosphorylation on tyrosine and threonine. *J. Biol. Chem.*, **270**, 7420-7426 (1995).
3. Foltz, I., *et al.*, Hemopoietic growth factors with the exception of interleukin-4 activate the p38 mitogen activated protein kinase pathway. *J. Biol. Chem.*, **272**, 3296-3301 (1997).

LCM/PHC 12/01

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.