

For life science research only.
Not for use in diagnostic procedures.



Digoxigenin-3-O- methylcarbonyl- ϵ -aminocaproic acid-N-hydroxysuccinimide ester

 **Version: 19**
Content Version: June 2021

Cat. No. 11 333 054 001 5 mg

Store product at -15 to -25°C .

1.	General Information	3
1.1.	Contents	3
1.2.	Storage and Stability	3
	Storage Conditions (Product)	3
	Reconstitution	3
1.3.	Additional Equipment and Reagent required	3
1.4.	Application	3
2.	How to Use this Product	4
2.1.	Before you Begin	4
	Working Solution.....	4
	Preparation of working solution for the labeling of proteins	4
2.2.	Protocols	4
	Labeling of oligonucleotides.....	4
	Isolation of labeled oligonucleotides.....	5
	Labeling of proteins.....	5
	Purification of labeled proteins.....	5
2.3.	Parameters	6
	Chemical Formula.....	6
	Chemical Name.....	6
	Structural formula.....	6
	Molecular Weight	6
3.	Additional Information on this Product	6
3.1.	Test Principle	6
4.	Supplementary Information	7
4.1.	Conventions.....	7
4.2.	Changes to previous version.....	7
4.3.	Trademarks.....	8
4.4.	License Disclaimer.....	8
4.5.	Regulatory Disclaimer.....	8
4.6.	Safety Data Sheet.....	8
4.7.	Contact and Support.....	8

1. General Information

1.1. Contents

Vial / Bottle	Label	Function / Description	Content
1	Digoxigenin-3-O-methylcarbonyl- ϵ -aminocaproic acid-N-hydroxysuccinimide ester	Lyophilized, white powder	1 vial, 5 mg

1.2. Storage and Stability

Storage Conditions (Product)

When stored at -15 to -25°C , the product is stable through the expiry date printed on the label.

Vial / Bottle	Label	Storage
1	Digoxigenin-3-O-methylcarbonyl- ϵ -aminocaproic-acid-N-hydroxysuccinimide ester	Store at -15 to -25°C .

Reconstitution

Proceed as described in the Protocols section, **Labeling of oligonucleotides**, **Labeling of proteins**.

1.3. Additional Equipment and Reagent required

For labeling of oligonucleotides

- Double-distilled water
 - Dimethylformamide (DMF)
 - Sodium borate, 100 mM, pH 8.5
 - Sodium acetate, 3 M, pH 8.5
 - Ethanol, A.R.
- i** Use fresh DMF to increase labeling efficiency.

For isolation of labeled oligonucleotide

- HPLC column, RP 18, 5 μm , 25 \times 0.4 cm
- Buffer A: Triethylammoniumacetate, 100 mM, pH 7
- Buffer B: Triethylammoniumacetate, 100 mM, acetonitrile (1:1), pH 7

For labeling of proteins

- Ethanol, A.R.
- Dimethylsulfoxide (DMSO), A.R.
- Phosphate-buffered saline (PBS), pH 8.5

For purification of labeled proteins

- Sephadex G25 column

1.4. Application

Digoxigenin-3-O-methylcarbonyl- ϵ -aminocaproic acid-N-hydroxysuccinimide ester can be used for labeling of:

- Oligonucleotides
- Proteins

2. How to Use this Product

2.1. Before you Begin

Working Solution

Preparation of working solution for the labeling of proteins

With stirring, dissolve 0.327 mg Digoxigenin-3-O-methylcarbonyl- ϵ -aminocaproic acid-N-hydroxysuccinimide ester in 8.18 μ l DMSO or ethanol.

2.2. Protocols

Labeling of oligonucleotides

i Typically, the oligonucleotide is prepared by chemical solid phase synthesis. It should be reacted in a final synthesis cycle with, for example, "Aminolink II", generating a 5'-terminal amino function. After cleavage of the protection group by concentrated ammonia, the oligomer solution is concentrated in a vacuum.

- 1 Dissolve the solid residue of the crude oligonucleotide in a mixture consisting of 300 μ l double-distilled water and 30 μ l 3 M sodium acetate.

- 2 Transfer the solution to a reaction vial, add 0.9 ml ethanol, mix, and chill 2 to 3 hours at -15 to -25°C .

- 3 Spin 15 minutes at $20,000 \times g$ and decant the supernatant.

- 4 Wash the pellet with 100 μ l ice-cold ethanol.

- 5 Spin 5 minutes at $20,000 \times g$ and decant supernatant.

- 6 Dissolve ethanol precipitated oligonucleotide (approximately 100 nM) in 50 μ l 100 mM sodium borate buffer, pH 8.5.

- 7 Dissolve 2 μ mol (1.3 mg) Digoxigenin-3-O-methylcarbonyl- ϵ -aminocaproic acid-N-hydroxysuccinimide ester in 50 μ l fresh DMF.
– Add this mixture to the oligonucleotide solution.
i A clear solution should be obtained. If necessary, add several drops of buffer or DMF to homogenize.

- 8 Incubate overnight at $+15$ to $+25^{\circ}\text{C}$.

Isolation of labeled oligonucleotides

i Since the labeling reaction is not quantitative, separation of the labeled oligonucleotide from the unlabeled compound must be performed. This can be done by reverse phase HPLC.

1 Concentrate the mixture from the labeling reaction under vacuum.

i Alternatively, perform an ethanol precipitation.

2 Dissolve the residue in 1 ml double-distilled water and apply onto an HPLC column, RP 18, 5 μ m, 25 \times 0.4 cm.

3 Perform the elution as follows:

Buffer	Description
A	Triethylammoniumacetate, 100 mM, pH 7
B	Triethylammoniumacetate, 100 mM, acetonitrile (1:1), pH 7

– Use a gradient in 30 minutes to 80% B.

i The digoxigenin-labeled oligonucleotide is eluted in the last fraction.

4 Concentrate the appropriate fraction under vacuum and desalt as usual, for example, gelfiltration on Sephadex or dialysis in a SPECTAPOR 1000).

Labeling of proteins

The following steps describe the labeling of 1 mg polyclonal antibody.

i The molar reaction mixture of antibody: Digoxigenin-3-O-methylcarbonyl- ϵ -aminocaproic acid-N-hydroxysuccinimide ester is 1:70.

1 Dissolve 1 mg polyclonal antibody in 1 ml PBS, pH 8.5.

2 Add Digoxigenin-3-O-methylcarbonyl- ϵ -aminocaproic acid-N-hydroxysuccinimide ester working solution, see section **Working Solution**.

3 Incubate 2 hours at +15 to +25°C.

Purification of labeled proteins

To remove unbound Digoxigenin-3-O-methylcarbonyl- ϵ -aminocaproic acid-N-hydroxysuccinimide ester and the solvent, dialyze the protein solution against PBS or use a Sephadex G25 column.

3. Additional Information on this Product


2.3. Parameters

Chemical Formula



Chemical Name

Structural formula



Molecular Weight

658.8 Da

3. Additional Information on this Product

3.1. Test Principle

The activated ester reacts with amino groups under mild conditions, therefore the digoxigenin moiety can be introduced, for example, into proteins or 5'-amino substituted oligonucleotides.

4. Supplementary Information

4.1. Conventions

To make information consistent and easier to read, the following text conventions and symbols are used in this document to highlight important information:

Text convention and symbols

 *Information Note: Additional information about the current topic or procedure.*

 **Important Note: Information critical to the success of the current procedure or use of the product.**

① ② ③ etc. Stages in a process that usually occur in the order listed.

① ② ③ etc. Steps in a procedure that must be performed in the order listed.

* (Asterisk) The Asterisk denotes a product available from Roche Diagnostics.

4.2. Changes to previous version

Layout changes.

Editorial changes.

4. Supplementary Information

4.3. Trademarks

All product names and trademarks are the property of their respective owners.

4.4. License Disclaimer

For patent license limitations for individual products please refer to:

List of biochemical reagent products.

4.5. Regulatory Disclaimer

For life science research only. Not for use in diagnostic procedures.

4.6. Safety Data Sheet

Please follow the instructions in the Safety Data Sheet (SDS).

4.7. Contact and Support

To ask questions, solve problems, suggest enhancements or report new applications, please visit our **Online Technical Support Site.**

To call, write, fax, or email us, visit **sigma-aldrich.com**, and select your home country. Country-specific contact information will be displayed.

