



**PROTEIN A-10 NM COLLOIDAL GOLD LABELED**  
Product Number **P 1039**

## Product Information

### Preparation Instructions

Product should be diluted for most applications. It is recommended that diluent buffer contain 0.15 M saline buffered at pH 6 to 8, plus 0.5% albumin (A 7638) and 0.05% Tween 20 to minimize background (additional buffer supplement may be required for certain applications e.g., see "dot blot" diluent). It is also recommended that, prior to application, the diluted conjugate be allowed to equilibrate at least 20 minutes in lower glycerol content. Optimum concentration of the conjugate must be determined empirically dependent on specific usage and generally may range from final  $A_{520} = 1.0$  to 0.05 (1:5 – 1:100 dilution) with incubation times ranging from 30 minutes to 12 hours.

### Storage/Stability

This undiluted product may be stored for extended periods at  $-20^{\circ}\text{C}$ . Diluted samples should not be stored below  $0^{\circ}\text{C}$  as freezing may cause aggregation of the colloid.

### Results

Binding is evaluated by a "dot blot" assay modified from the method of Brada and Roth.<sup>2</sup> Serial dilutions are

prepared from a 1 mg/ml positive control protein solution. One micro-liter ( $1\ \mu\text{l}$ ) of each solution is adsorbed onto a nitrocellulose membrane and allowed to dry. The gold conjugate is diluted to  $A_{520} = 0.25$  (approx. 1:20) with 0.15 M NaCl, 0.01 M sodium phosphate, pH 7.0, 5 mg/ml albumin and 0.05% Tween 20. The spotted membranes are incubated with the gold for 1 hour at  $25^{\circ}\text{C}$ . The detection limit is the minimum amount of protein that can be detected as a pink-red spot on the membrane.

### References

1. Ackerman, G.A., *et al.*, J. Histochem. Cytochem., **31**, 433 (1983).
2. Brada, D. and Roth, J. Anal. Biochem., **142**, 70 (1984).
3. Bendayan, M., J. Electron Microscopy Technique, **1**, 243 (1984).
4. Bendayan, M., J. Histochem. Cytochem., **30**, 81 (1981).

KMR/JWM 01/02

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.