3050 Spruce Street, St. Louis, MO 63103 USA
Tel: (800) 521-8956 (314) 771-5765 Fax: (800) 325-5052 (314) 771-5757
email: techservice@sial.com sigma-aldrich.com

# **Product Information**

# SILu™MAb Mouse - Stable Isotope Labeled Monoclonal Antibody Standard

mouse, recombinant, expressed in CHO cells

Catalog Number **MSQC10** Storage Temperature –20 °C

## **Product Description**

SILu™MAb Mouse is a recombinant, stable isotope-labeled, mouse monoclonal antibody, which incorporates [¹³C<sub>6</sub>, ¹⁵N<sub>4</sub>]-Arginine and [¹³C<sub>6</sub>, ¹⁵N<sub>2</sub>]-Lysine. Expressed in CHO cells, SILuMAb Mouse is designed to be used as a universal internal standard for bioanalysis of monoclonal murine antibodies in animal pharmacokinetics (PK) studies.

SILuMAb Mouse is an IgG1 antibody with a kappa light chain. Because of overlap with common sequences in the Fc region with candidate murine antibodies, SILuMAb Mouse provides universal utility, and thus, eliminates the need to produce candidate-specific internal standards.

Each vial of SILuMAb Mouse contains the labeled antibody lyophilized from a solution of phosphate buffered saline. Vial content was determined by measuring  $A_{280}$  and using an extinction coefficient ( $E^{0.1\%}$ ) of 1.4.

# Sequence Information

# SILuMAb Mouse Heavy Chain:

QVQLQQSAAELARPGASVKMSCKASGYSFTTYTIHWVKQRPGQGL EWIGYINPSSGYAAYNQNFKDETTLTADPSSSTAYMELNSLTSED SAVYYCAREKFYGYDYWGQGATLTVSSAKTTPPSVYPLAPGSAAQ TNSMVTLGCLVKGYFPEPVTVTWNSGSLSSGVHTFPAVLQSDLYT LSSSVTVPSSTWPSETVTCNVAHPASSTKVDKKIVPRDCGCKPCI CTVPEVSSVFIFPPKPKDVLTITLTPKVTCVVVDISKDDPEVQFS WFVDDVEVHTAQTQPREEQFNSTFRSVSELPIMHQDWLNGKEFKC RVNSAAFPAPIEKTISKTKGRPKAPQVYTIPPPKEQMAKDKVSLT CMITDFFPEDITVEWQWNGQPAENYKNTQPIMDTDGSYFVYSKLN VOKSNWEAGNTFTCSVLHEGLHNHHTEKSLSHSPG

## SILuMAb Mouse Light Chain:

DVLMTQIPLSLPVSLGDQASISCRSSQSIVHRNGNTYLEWYLLKP GQSPKLLIYKVSNRFSGVPDRFSGSGSGTDFTLKISRVEAEDLGV YYCFQGSHVPYTFGGGTKLEIRRADAAPTVSIFPPSSEQLTSGGA SVVCFLNNFYPKDINVKWKIDGSERQNGVLNSWTDQDSKDSTYSM SSTLTLTKDEYERHNSYTCEATHKTSTSPIVKSFNRNEC

#### **Precautions and Disclaimer**

This product is for R&D use only, not for drug, household, or other uses. Please consult the Safety Data Sheet for information regarding hazards and safe handling practices.

# **Preparation Instructions**

SILuMAb Mouse recovery is maximized when 0.1% formic acid is used for reconstitution of the lyophilized product.

Reconstitution with other solvents may reduce recovery. Do not freeze after reconstitution.

- 1. Briefly centrifuge the vial at  $\sim$ 10,000  $\times$  g to collect the product at the bottom of the vial.
- Add 500 μL of ultrapure water containing 0.1% formic acid to the vial.
- 3. Mix the contents by gently inverting the vial a minimum of 5 times.
- 4. Allow the vial to stand at room temperature for at least 15 minutes and repeat mixing by inversion.

# Storage/Stability

Store the lyophilized product at -20 °C.

SILu is a trademark of Sigma-Aldrich Co. LLC.

JK,MAM 11/16-1