

Product Information

N-26 Sequencing Primer

Catalog Number **P7832**

Store Temperature $-20\text{ }^{\circ}\text{C}$

Product Description

Nucleotide sequence:

5'-OH-CAT-CAT-AAC-GGT-TCT-GGC-AAA-TAT-TC-
-3-OH'

The N-26 sequencing primer is a chemically synthesized, 26 base, single strand oligodeoxy-ribonucleotide designed for double-stranded or single-stranded DNA sequencing of FLAG® fusion junctions corresponding to the N-terminus of FLAG fusion proteins expressed by *E. coli* pFLAG expression vectors. The N-26 primer is especially useful for sequencing N-terminal FLAG fusion junctions encoded by inserts cloned into the pFLAG-MAC™, pFLAG-ATS™, pFLAG-1, pFLAG-2, pFLAG-Shift™₁₂, and pFLAG-Shift™_{12c} N-terminal expression vectors. For all pFLAG expression vectors, the N-26 primer is complementary to the plus strand map positions 1 to 26.

The N-26 sequencing primer is supplied at an initial concentration of 5 μM in 0.1 \times TE buffer (1 mM Tris, pH 8.0, 0.1 mM EDTA).

Reagents Required but Not Provided

- 1 \times TE buffer (10 mM Tris, pH 8.0, 1.0 mM EDTA)
- 5 M NaOH
- 3 M potassium/5 M acetate: To 60 ml of 5 M potassium acetate add 11.5 ml of glacial acetic acid and 28.5 ml of water. The resulting solution is 3 M with respect to potassium and 5 M with respect to acetate.
- Isopropanol

Procedure

- A. Preparation of Denatured pFLAG-1 DNA Template
(For sequencing single-stranded DNA, proceed to Section B.)

The following protocol is designed to allow sufficient template for two sets of 4 sequencing reactions (2G, 2A, 2T and 2C reactions).

1. Denature 7 micrograms of highly purified pFLAG-1 supercoiled DNA in 70 μl of 1 \times TE by adding 3 μl (1/25 volume) of 5 M NaOH and incubating at 37 $^{\circ}\text{C}$ for 5 minutes.
2. Precipitate the denatured pFLAG-1 supercoiled DNA for 30 minutes at room temperature by adding 150 μl (2 volumes) of 3 M potassium/5 M acetate: isopropanol (1:3 ratio mixture).
3. Collect the precipitated and denatured pFLAG-1 DNA by centrifugation at 10,000 $\times g$ for 5 minutes.
4. Wash the pFLAG-1 DNA pellet with 1 ml ethanol and dry.
5. Resuspend the pFLAG-1 DNA in 20 μl of 1 \times TE. The final concentration should be 0.35 $\mu\text{g}/\mu\text{l}$.

The irreversibly denatured pFLAG-1 DNA template can be stored at $-20\text{ }^{\circ}\text{C}$ at this point if desired.

B. Priming pFLAG-1 DNA Template with N-26 Sequencing Primer

The following protocol uses 10 μ l or half of the preceding preparation of irreversibly denatured pFLAG-1 DNA template. This is sufficient for one set of 4 sequencing reactions (G, A, T and C).

1. Dilute a 3 μ l aliquot of the N-26 sequencing primer with 6 μ l of 1 \times TE to make a final concentration of 1.67 pmol/ μ l.
2. Add 2-3 μ l (3.5-5 pmoles) of N-26 sequencing primer to 10 μ l (1 pmole) of denatured pFLAG-1 DNA template.
3. To 12-13 μ l of N-26 primer/pFLAG-1 DNA template add an appropriate volume of sequencing buffer to make the buffer concentration 1 \times . Heat at 70 $^{\circ}$ C in a wet temperature block for 2 minutes.
4. Slowly cool the mixture to 45 $^{\circ}$ C by placing the temperature block at room temperature for about 20 minutes.
5. Distribute the N-26 primer/pFLAG-1 DNA template to four tubes to be used in the G, A, T and C DNA sequencing reaction.

The N-26 primer/pFLAG-1 DNA template is now ready for supercoil sequencing of the DNA sequence corresponding to the N-terminal FLAG fusion junction.

The DNA sequence corresponding to the fusion junction will be approximately 150 bases away from the 3' end of the N-26 sequencing primer reading from the bottom of an autoradiogram.

FLAG is a registered trademark of Sigma-Aldrich[®] Biotechnology LP and Sigma-Aldrich Co. pFLAG-MAC, pFLAG-ATS, and pFLAG-Shift are trademarks of Sigma-Aldrich Biotechnology LP and Sigma-Aldrich Co.

MAM 09/10-1

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.