

Detection Platform Compatibility with SmartFlare™ Probes

Technology	Imaging						Non-Imaging				
	High Content	Microscopy				Imaging Flow Cytometry	Flow Cytometers	Fluorescence-Activated Cell Sorters	Plate Readers*		
		Wide field/ Epi-fluorescence	Confocal	Two-photon Confocal	Imaging plate readers*						
Examples	Biotek Cytation™ 3 PerkinElmer Opera™	Evos® Station	Nikon C2	Zeiss LSM 510-NLO	PerkinElmer EnSight™	Merck Millipore Amnis® brand FlowSight® and ImageStream®X Systems	Merck Millipore Muse® Cell Analyzer	Merck Millipore guava easyCyte™ Flow Cytometer	FACSAria™	Tecan Infinite® M1000	
SmartFlare™ Fit	HIGH	MED	HIGH	HIGH	MED	HIGH	HIGH	HIGH	HIGH	LOW*	
Dye Compatibility	Cy 3 Ex 550 Em 570	Yes (mono or filters)	Yes (530 nm excitation via lightcube)	Yes (561 nm excitation)	Yes (543 nm excitation)	Yes (mono)	Yes (561 nm excitation)	Yes (532 nm excitation)	Not optimal (488 nm excitation only)	Yes (561 nm excitation)	Yes (530 nm excitation)
	Cy 5 Ex 625 Em 670	Yes (mono or filters)	Yes (via optional Cy5 lightcube)	Yes (633 nm excitation)	Yes (633 nm excitation)	Yes (mono)	Yes (642 nm excitation)	No	Yes (640 nm excitation)	Yes (633 nm excitation)	Yes (635 nm excitation)
	Cy5 / Cy3 cross-compatible	Yes (mono)	No (Bleed over from channels can occur if no PMT)	Yes	Yes	Yes (+25% detection time for two-color)	Yes (with 561 nm and 642 nm laser options)	No	Not optimal	Yes	Yes
Applications	<ul style="list-style-type: none"> HTS screens Time-lapse experiments Live cell assays 	<ul style="list-style-type: none"> Live cell applications (using heating tray) 	<ul style="list-style-type: none"> Morphology - Timecourse/ Live cell imaging with incubator 	<ul style="list-style-type: none"> Morphology in whole cells Timecourse/ Live cell imaging with incubator 	<ul style="list-style-type: none"> Cell population imaging Identification of positive cell events within population 	<ul style="list-style-type: none"> Single-cell imaging, morphology Identification of positive cell events within population Detection of rare events in populations 	<ul style="list-style-type: none"> Single cell analysis 	<ul style="list-style-type: none"> Single cell analysis 	<ul style="list-style-type: none"> Separation of cell populations from heterogeneous samples Enrichment tool 	<ul style="list-style-type: none"> Top and bottom well readings 	
Typical sample size (# cells)	Up to 384-well plate, T-25 flask, up to 100 mm dishes; cell number varies	Live cell culture vessels (dishes up to 100 mm, flasks up to T-75); cell number varies	8-well or 4-well chamber slides; 35 mm dishes and multiwell plates (~10 ⁵ cells per well)	8-well or 4-well chamber slides; 35 mm dishes and multiwell plates (~10 ⁵ cells per well)	Up to 384-well plates	Both systems will image ~95% of cells in a 200 µL sample, regardless of the number of cells present.	Single-sample reads in tubes (200 µL) 30,000 cells per tube	96-well plate reads (~200 µL per well) 30,000 cells per tube	5 mL polypropylene tube (1 - 3 mL volume) up to ~10 ⁷ cells per sample	Up to 1536 well plate; volume and cell numbers vary	

*Usage on non-imaging plate readers as not been validated in our labs but is likely provided the platform has the correct Ex and Em wavelength capability and can perform bottom. Most Non-imaging readers only allow for whole-well fluorescence reads, resulting in population averages, not single cell resolution which SmartFlare™ probes provide.