

Product Information

Endoglycoceramidase II recombinant, from *Rhodococcus sp.* expressed in *E. coli*

Product Number **E 9030**
 Storage Temperature $-20\text{ }^{\circ}\text{C}$

E.C. 3.2.1.123
 Synonyms: ECGase II; Oligoglycosylglucosylceramide
 glycohydrolase

Product Description

Endoglycoceramidase II (ECGase II) is an endo-glycohydrolase. The enzyme is useful for structural studies of glycosphingolipids. In the presence of detergents, it hydrolyzes the linkage between the oligosaccharide and ceramide of acidic and neutral glycosphingolipids, resulting in intact oligosaccharides and ceramides. The enzyme does not act on galactosyl- and glucosyl-ceramides. It does not hydrolyze cerebrosides, sulfatides, and sphingomyelin. The galactosylceramide linkage of gala-type glycosphingolipids is not hydrolyzed. Globo-type glycosphingolipids are strongly resistant to hydrolysis by ECGase II.

Molecular weight: 58.9 kDa

Optimal pH: 5.5

Inhibition: 1 mM of Hg^{2+} , Zn^{2+} , and Cu^{2+}

This enzyme preparation is isolated and purified from *E. coli* carrying the gene for the *Rhodococcus sp.* enzyme. The product is supplied as 0.1 unit of ECGase II in 20 mM sodium acetate buffer, pH 6.0, containing 0.2% BSA and 0.1% Lubrol PX. This preparation contains detergent and does not contain Activator II.

Purity: This preparation is essentially free from the following exoglycosidase and other enzyme activities:

α -galactosidase, β -galactosidase, α -mannosidase, α -N-acetylgalactosaminidase, α -fucosidase, β -N-acetylgalactosaminidase, sialidase, β -N-acetylglucosaminidase, glycopeptidase, endo- β -N-acetylglucosaminidase, proteinase, and sphingomyelinase.

Substrate Specificity:
 (1.5 milliunits of the enzyme was incubated with 10 nmole of the substrate at $37\text{ }^{\circ}\text{C}$ for 16 hours.)

Substrate	Hydrolysis (%)
Ganglio series	
G _{T1b}	100
G _{D1a}	100
G _{M1a}	100
G _{M3}	100
Asialo G _{M1}	100
Globo Series	
Gb5Cer	19
Gb4Cer	10
Neolacto series	
IV ³ NeuAc α -nLC ₄	100
IV ⁶ NeuAc α -nLC ₄	100
IV ³ NeuAc α III ³ Fuc α -nLC ₄ (SialylLewis X)	100
III ³ Fuc α -nLC ₄ (SialylLewis X)	100
Lacto series	
IV ³ NeuAc α -LC ₄	100
IV ⁶ NeuAc α -LC ₄	100
Lactosylceramide	100
Cerebrosides	
Glucosyl ceramide	0
Galactosyl ceramide	0
Sulfatide	0

Unit Definition: One unit will catalyze the hydrolysis of 1 μ mole of asialo-G_{M1} per minute at $37\text{ }^{\circ}\text{C}$ at pH 5.0.

Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

Storage/Stability

The product ships on dry ice and storage at or below $-20\text{ }^{\circ}\text{C}$ is recommended until use. Store the solution in aliquots at $2-8\text{ }^{\circ}\text{C}$ once thawed and use within 2-3 days.

References

1. Ito, M., and Yamagata, T., J. Biol. Chem., **261**, 14278-14282 (1986).
2. Ito, M., and Yamagata, T., J. Biol. Chem., **264**, 9510-9519 (1989).
3. Ito, M. et al, J. Biol. Chem., **266**, 7919-7926 (1991).
4. Ito, M. et al, J. Biochem., **110**, 328-332 (1991).
5. Ito, M. et al, Eur. J. Biochem., **218**, 637-643 (1993).
6. Ito, M. et al, Eur. J. Biochem., **218**, 645-649 (1993).
7. Muramoto, K. et al, Biochem. Biophys. Res. Comm., **202**, 398-402 (1994).

RBG/MAM 08/04

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.