Instructions for Use

Leukocyte Peroxidase
(Myeloperoxidase) Kit

Procedure No. 391


Intended Use

Sigma-Aldrich Peroxidase (Myeloperoxidase) kit is intended for use as a general laboratory histochemical stain for the demonstration of leukocyte peroxidase. The Peroxidase (Myeloperoxidase) kit is for professional "In Vitro Diagnostic Use" only. This manual, qualitative procedure identifies polymorphonucleocytes in blood or bone marrow films of human specimens. Histological visualization and identification of polymorphonucleocytes is a unique technique currently commonly used in medicine.

Diaminobenzidine (DAB), as a benzidine substitute for peroxidase (myeloperoxidase) cytochemistry, has less tinctorial power than benzidine.^{1,2} This made DAB much less attractive to histochemists. However, Hanker and associates improved DAB methodology, making it more suitable for differentiating granulocytes, their precursors and monocytes from cells of lymphoid origin.³⁻⁶ According to their modification, the brown reaction product is first intensified with copper salts followed by application of Gill's modified Papanicolaou stain, resulting in intense grey-black granules at sites of neutrophil and monocyte myeloperoxidase.⁷

The Sigma-Aldrich procedure is similar to that developed by Hanker, involving the following reactions:

When treated with Gill's modified Papanicolaou stain, the reaction product is further intensified and characteristic color imparted to neutrophils, eosinophils and basophils. It should be noted that evidence indicating the carcinogenicity of DAB is equivocal and good laboratory practice should preclude any potential hazard.⁸

Reagents

Diaminobenzidine (Cat. No. 3911-1VL)

Diaminobenzidine tetrahydrochloride, 25 mg/vial (10 vials). Danger. Suspected of causing genetic defects. May cause cancer. Obtain special instructions before use. Use personal protective equipment as required. If exposed or concerned: Get medical advice/attention.

Copper Nitrate (Cat. No. 3912-1VL)

Copper nitrate, 625 mg/vial (2 vials). Danger. May intensify fire, oxidizer. Harmful if swallowed. Causes skin irritation. Causes serious eye damage. Very toxic to aquatic life.

TRIZMA® Buffer Concentrate (Cat. No. 3913-50ML)

TRIZMA®-HCl buffer, 1.0 mol/L. Warning. Causes skin irritation. Causes serious eye irritation.

Glutaraldehyde Solution (Cat. No. 3802-75ML)

Glutaraldehyde, 4%, and borate buffer, 67 mmol/L, pH 7.6. (2 x 75 mL)

Hematoxylin Solution, Gill No. 3 (Cat. No. GHS3-100ML)

Hematoxylin (certified), C.I. 75290, 6 g/L, sodium iodate, 0.6 g/L, aluminum sulfate, 52.8 g/L, and stabilizer.

Gill Modified EA Solution (Cat. No. 3915-100ML)

Fast Green FCF (certified), C.I. 42053, 0.017% (w/v), Eosin Y (certified), C.I. 45380, 0.4% (w/v), alcohol, 73% (v/v), absolute methyl alcohol, 25% (v/v), glacial acetic acid, 2% (v/v) and phosphotungstic acid, 0.4% (w/v).

Scott's Tap Water Substitute Concentrate (Cat. No. S5134-100ML)

Contains magnesium sulfate • 7H₂O, 200 g/L, sodium bicarbonate, 20 g/L, and preservative

Special Materials Required but Not Provided

- Absolute Ethanol
- Ethanol, 95%
- Hydrogen Peroxide, 1%
- Acetone, ACS Reagent
- Xylene or xylene substitute

Storage and Stability

Store Diaminobenzidine and Glutaraldehyde Solution in refrigerator (2–8°C). Reagent labels bear expiration date. Store other kit reagents at room temperature (18–26°C). Protect Hematoxylin Solution from light. Reagent labels for Copper Nitrate and Gill Modified EA Solution bear expiration date.

Store Copper Nitrate Solution and Scott's Tap Water Substitute Working Solution at room temperature (18–26°C). Store TRIZMA® Working Solution in refrigerator (2–8°C). If turbidity develops, discard Copper Nitrate Solution, TRIZMA® Buffer Concentrate and Working Solution, Scott's Tap Water Substitute Concentrate and Working Solution.

Store Glutaraldehyde-Acetone Fixative tightly stoppered in refrigerator (2–8°C). Stable provided pH is in range 7.2 to 8.0.

Deterioration

Discard Hematoxylin Solution, Gill No. 3, if solution turns brown (air oxidation) or purple (loss of acidity).

Preparation

Copper Nitrate Solution is prepared by dissolving contents of 1 vial Copper Nitrate in 250 mL deionized water.

TRIZMA® Working Solution (pH 7.6 ± 0.3) is prepared by diluting 1 volume of TRIZMA® Buffer Concentrate with 9 volumes deionized water.

Glutaraldehyde-Acetone Fixative Solution is prepared by adding 25 mL reagent grade acetone to 75 mL Glutaraldehyde Solution.

Scott's Tap Water Substitute Working solution is prepared by diluting 1 volume Scott's Tap Water Substitute Concentrate with 9 volumes deionized water.

Filter Hematoxylin Solution, Gill No. 3, before use.

Precautions

These IVDs are intended for in vitro diagnostic use in a clinical laboratory environment. These IVDs are for professional use by qualified personnel only. Sigma-Aldrich IVDs may be operated by laboratory personnel who are trained to handle human specimens that can be infectious, use microscopes and other laboratory equipment and have color perception and visual acuity to distinguish colors and other objects under a microscope.

Normal precautions exercised in handling laboratory reagents should be followed. Dispose of waste observing all local, state, provincial or national regulations.

Procedure

Specimen Collection

No known test method can offer complete assurance that blood samples or tissue will not transmit infection. Therefore, all blood derivatives or tissue specimens should be considered potentially infectious.

Bone marrow aspirate, peripheral blood or buffy coat preparations may be used. Either heparin or EDTA serves as a satisfactory anticoagulant.³ Unfixed samples may be stored at room temperature (18–26°C) protected from light for several months without loss of activity.⁴

Note

It is recommended that blood films prepared from healthy donors be processed along with patient samples to monitor system performance.

Although myeloperoxidase is generally considered a marker for cells of myelocytic lineage, it is imperative to recognize that monoytoid cells may also display weak peroxidase activity.

Procedure

1. Dissolve contents of 1 Diaminobenzidine vial in 50 mL TRIZMA® Working Solution.
2. Set up a series of staining jars containing the following:
 - Glutaraldehyde-Acetone Fixative Solution
 - Diaminobenzidine Solution prepared in Step 1
 - Copper Nitrate Solution
 - Hematoxylin Solution, Gill No. 3.
 - Scott's Tap Water Substitute Working Solution
 - Gill's Modified EA Solution
 - Two dishes with 95% ethanol
 - Two dishes with absolute ethanol
 - Three dishes with xylene. See NOTE following Step 18
3. Immediately prior to fixing slides, add 0.5 mL 1% hydrogen peroxide to Diamino-benzidine solution prepared in Step 1. Mix well.
4. Fix slides at 4–8°C in Glutaraldehyde-Acetone Fixative Solution for 1 minute.
5. Rinse briefly (30 seconds) in deionized water.
6. Incubate for 45 seconds in Diaminobenzidine/Peroxide solution prepared in Step 3.
7. Rinse briefly (30 seconds) in deionized water.
8. Immerse in Copper Nitrate Solution for 2 minutes with gentle agitation.
9. Rinse briefly (30 seconds) in deionized water.
10. Immerse in Hematoxylin Solution, Gill No. 3 for 8 seconds (4 dips).
11. Rinse in 2 changes of deionized water for 5 seconds with agitation.
12. Immerse for 12 seconds (6 dips) in Scott's Tap Water Substitute Working Solution.
13. Rinse in 2 changes of deionized water for 5 seconds.
14. Immerse in Gill Modified EA Solution for 1 minute.
15. Rinse in 2 changes of 95% ethanol for 3 seconds each.
16. Rinse in 2 changes absolute ethanol for 3 seconds each.
17. Rinse in 3 changes xylene for 3 seconds each.
18. Mount in permanent mounting media and examine microscopically. Color will fade if slides are not mounted.

Note

The procedure may be shortened by eliminating Steps 14–17. Eosinophil peroxidase will appear dark red-brown instead of red-orange. If Gill Modified EA Solution is used, rinsing in alcohol and xylene must be included.

Performance Characteristics

Neutrophils and their precursors show grey-black intracellular granulation. Monocytes stain less intensely. Eosinophils stain red-orange while basophils stain blue. Lymphocytes do not show peroxidase activity.

Blood films prepared from normal donors were stained for myeloperoxidase according to this procedure and by a benzidine method.¹ Neutrophils showed brown-black granulation with this procedure and blue granulation with the benzidine procedure. In both cases, monocytes stained less intensely and lymphocytes did not show myeloperoxidase activity.

If observed results vary from expected results, please contact Sigma-Aldrich Technical Service for assistance.

Analytical Performance Characteristics

The analytical performance results for the given tests conducted on all target structures, confirm 100% sensitivity, specificity, and repeatability.

Cat. No	Product Description	Target	Intra-assay Specificity	Intra-assay Sensitivity	Inter-assay Specificity	Inter-assay Sensitivity
3911	Diaminobenzidine	Sites of Activity	3 of 3	3 of 3	3 of 3	3 of 3
3912	Copper Nitrate	Sites of Activity	3 of 3	3 of 3	3 of 3	3 of 3
3913	TRIZMA® Buffer Concentrate	Osmolality	3 of 3	3 of 3	3 of 3	3 of 3
		pH	3 of 3	3 of 3	3 of 3	3 of 3
3802	Glutaraldehyde Solution	Sites of Activity	3 of 3	3 of 3	3 of 3	3 of 3
GHS3	Hematoxylin Solution, Gill No. 3	Nuclei	3 of 3	3 of 3	3 of 3	3 of 3
3915	Gill Modified EA Solution	Cytoplasm	3 of 3	3 of 3	3 of 3	3 of 3
S5134	Scott's Tap Water Substitute Concentrate	Osmolality	3 of 3	3 of 3	3 of 3	3 of 3
		pH	3 of 3	3 of 3	3 of 3	3 of 3

Warnings and Hazards

Refer to Safety Data Sheet and product labeling for any updated risk, hazard or safety information.

391A:

- H225 Highly flammable liquid and vapour.
- H272 May intensify fire, oxidizer.
- H290 May be corrosive to metals.
- H302 + H312 + H332 Harmful if swallowed, in contact with skin or if inhaled.
- H315 Causes skin irritation.
- H317 May cause an allergic skin reaction.
- H318 Causes serious eye damage.
- H334 May cause allergy or asthma symptoms or breathing difficulties if inhaled.
- H335 May cause respiratory irritation.
- H341 Suspected of causing genetic defects.
- H350 May cause cancer.
- H360 May damage fertility. May damage the unborn child.
- H370 Causes damage to organs.
- H373 May cause damage to organs through prolonged or repeated exposure if swallowed.
- H410 Very toxic to aquatic life with long lasting effects.
- P201 Obtain special instructions before use.
- P210 Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.
- P220 Keep away from clothing and other combustible materials.
- P260 Do not breathe dust/ fume/ gas/ mist/ vapours/ spray.
- P273 Avoid release to the environment.
- P280 Wear protective gloves/ protective clothing/ eye protection/ face protection.
- P304 + P340 + P312 IF INHALED: Remove person to fresh air and keep comfortable for breathing. Call a POISON CENTER/doctor if you feel unwell.
- P305 + P351 + P338 + P310 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Immediately call a POISON CENTER/ doctor.
- P308 + P311 IF exposed or concerned: Call a POISON CENTER/doctor.
- P370 + P378 In case of fire: Use dry sand, dry chemical or alcohol-resistant foam to extinguish.

If during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.

Symbol Definitions

Symbols as defined in EN ISO 15223-1:2021

	Manufacturer		Catalogue Number
	Consult Instructions for Use		Batch Code
	Authorized Representative in the European Community/ European Union		European Union Declaration of Conformity (defined in IVDR 2017/746)
	Use-by Date		In vitro diagnostic medical device
	Temperature Limit		Caution
	Date of Manufacture		Importer
	Indicates the Authorised Representative in Switzerland		

References

- Graham RC, Karnovsky MJ: Glomerular permeability. Ultrastructural cytochemical studies using peroxidases as protein tracers. J Exp Med 124:1123, 1966
- Graham RC, Karnovsky MJ: The early stages of absorption of injected horseradish peroxidase in the proximal tubules of mouse kidney: Ultrastructural cytochemistry by a new technique. J Histochem Cytochem 14:291, 1966
- Hanker JS, Laszlo J, Moore JO: The light microscopic demonstration of hydroperoxidase-positive phi bodies and rods in leukocytes in acute myeloid leukemia. Histochemistry 58:241, 1978
- Hanker JS, Ambrose WW, James CJ, et al: Facilitated light microscopic cytochemical diagnosis of acute myelogenous leukemia. Cancer Res 39:1635, 1979
- Hanker JS, Chandross RJ, Weatherly NF et al: Medusa cells: The morphology and cytochemistry of common amoeboid variants of eosinophils. Histochem J 12:701, 1980
- Hanker JS, Giammara BL: Haematology: Some Advances and Insights from Histochemical and Molecular Biological Studies of Haemoproteins, IN Histochemistry: The Widening Horizons. P.J. Stoward and JM Polak, Editors, John Wiley & Sons Ltd., Sussex, 1991
- Gill GW, Frost JK, Miller KA: A new formula for a half-oxidized hematoxylin solution that neither overstains nor requires differentiation. Acta Cytol 18:300, 1974
- Marmont AM, Damasio E, Zucker-Franklin D: Neutrophils. IN Atlas of Blood Cells - Function and Pathology, Vol 1, D Zucker-Franklin, MF Greaves, CE Grossi, AM Marmont, Editors, Lea and Febiger, Philadelphia, 1981, pp 149-242

Contact Information

To place an order, please visit our web site at [SigmaAldrich.com](https://www.sigmaaldrich.com). For Technical Service, please visit the tech service page on our web site at [SigmaAldrich.com/techservice](https://www.sigmaaldrich.com/techservice).

Revision History

Rev. 3.0	2014
Rev. 4.0	2016
Rev. 5.0	2023
Transferred to new template with current branding. Specified for professional use in intended use and precautions. Revised intended use to align with IVDR guidelines. Updated Material Safety Data Sheet to Safety Data Sheet. Updated contact information. Removed instruction to follow CLSI for specimen collection. Removed EN 980 and changed to EN ISO 15223-1:2021 for symbols. Added adverse event contact information. Added Warnings and Hazards. Added CH-REP information	

Sigma-Aldrich, Inc.,
3050 Spruce Street,
St. Louis, MO 63103 USA
an affiliate of Merck KGaA,
Darmstadt, Germany
+1(314) 771-5765

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Merck KGaA,
64271
Darmstadt,
Germany

MDSS CH GmbH
Laurenzenvorstadt 61
5000 Aarau
Switzerland

The Initial M and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

© 2023 Merck KGaA, Darmstadt, Germany and/or its affiliates. All rights reserved.

The life science business of Merck KGaA, Darmstadt, Germany operates as MilliporeSigma in the U.S. and Canada.

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Referenzen

- Graham RC, Karnovsky MJ: Glomerular permeability. Ultrastructural cytochemical studies using peroxidases as protein tracers. J Exp Med 124:1123, 1966
- Graham RC, Karnovsky MJ: The early stages of absorption of injected horseradish peroxidase in the proximal tubules of mouse kidney: Ultrastructural cytochemistry by a new technique. J Histochn Cytochem 14:291, 1966
- Hanker JS, Laszlo J, Moore JO: The light microscopic demonstration of hydroperoxidase-positive phi bodies and rods in leukocytes in acute myeloid leukemia. Histochemistry 58:241, 1978
- Hanker JS, Ambrose WW, James CJ, et al: Facilitated light microscopic cytochemical diagnosis of acute myelogenous leukemia. Cancer Res 39:1635, 1979
- Hanker JS, Chandross RJ, Weatherly NF et al: Medusa cells: The morphology and cytochemistry of common amoeboid variants of eosinophils. Histochem J 12:701, 1980
- Hanker JS, Giammara BL: Haematology: Some Advances and Insights from Histochemical and Molecular Biological Studies of Haemoproteins, IN Histochemistry: The Widening Horizons. P.J. Stoward and JM Polak, Editors, John Wiley & Sons Ltd., Sussex, 1991
- Gill GW, Frost JK, Miller KA: A new formula for a half-oxidized hematoxylin solution that neither overstains nor requires differentiation. Acta Cytol 18:300, 1974
- Marmont AM, Damasio E, Zucker-Franklin D: Neutrophils. IN Atlas of Blood Cells - Function and Pathology, Vol 1, D Zucker-Franklin, MF Greaves, CE Grossi, AM Marmont, Editors, Lea and Febiger, Philadelphia, 1981, pp 149-242

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Bibliografia

- Graham RC, Karnovsky MJ: Glomerular permeability. Ultrastructural cytochemical studies using peroxidases as protein tracers. J Exp Med 124:1123, 1966.
- Graham RC, Karnovsky MJ: The early stages of absorption of injected horseradish peroxidase in the proximal tubules of mouse kidney: Ultrastructural cytochemistry by a new technique. J Histochem Cytochem 14:291, 1966.
- Hanker JS, Laszlo J, Moore JO: The light microscopic demonstration of hydroperoxidase-positive phi bodies and rods in leukocytes in acute myeloid leukemia. Histochemistry 58:241, 1978.
- Hanker JS, Ambrose WW, James CJ, et al: Facilitated light microscopic cytochemical diagnosis of acute myelogenous leukemia. Cancer Res 39:1635, 1979.
- Hanker JS, Chandross RJ, Weatherly NF et al: Medusa cells: The morphology and cytochemistry of common amoeboid variants of eosinophils. Histochem J 12:701, 1980.
- Hanker JS, Giammara BL: Haematology: Some Advances and Insights from Histochemical and Molecular Biological Studies of Haemoproteins, IN Histochemistry: The Widening Horizons. P.J. Stoward and JM Polak, Editors, John Wiley & Sons Ltd., Sussex, 1991.
- Gill GW, Frost JK, Miller KA: A new formula for a half-oxidized hematoxylin solution that neither overstains nor requires differentiation. Acta Cytol 18:300, 1974.
- Marmont AM, Damasio E, Zucker-Franklin D: Neutrophils. IN Atlas of Blood Cells - Function and Pathology, Vol 1, D Zucker-Franklin, MF Greaves, CE Grossi, AM Marmont, Editors, Lea and Febiger, Philadelphia, 1981, pp 149-242.

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Referencias

- Graham RC, Karnovsky MJ: Glomerular permeability. Ultrastructural cytochemical studies using peroxidases as protein tracers. J Exp Med 124:1123, 1966
- Graham RC, Karnovsky MJ: The early stages of absorption of injected horseradish peroxidase in the proximal tubules of mouse kidney: Ultrastructural cytochemistry by a new technique. J Histochem Cytochem 14:291, 1966.
- Hanker JS, Laszlo J, Moore JO: The light microscopic demonstration of hydroperoxidase-positive phi bodies and rods in leukocytes in acute myeloid leukemia. Histochemistry 58:241, 1978
- Hanker JS, Ambrose WW, James CJ, et al: Facilitated light microscopic cytochemical diagnosis of acute myelogenous leukemia. Cancer Res 39:1635, 1979
- Hanker JS, Chandross RJ, Weatherly NF et al: Medusa cells: The morphology and cytochemistry of common amoeboid variants of eosinophils. Histochem J 12:701, 1980
- Hanker JS, Giammara BL: Haematology: Some Advances and Insights from Histochemical and Molecular Biological Studies of Haemoproteins, IN Histochemistry: The Widening Horizons. P.J. Stoward y JM Polak, Editores, John Wiley & Sons Ltd., Sussex, 1991
- Gill GW, Frost JK, Miller KA: A new formula for a half-oxidized hematoxylin solution that neither overstains nor requires differentiation. Acta Cytol 18:300, 1974
- Marmont AM, Damasio E, Zucker-Franklin D: Neutrophils. IN Atlas of Blood Cells - Function and Pathology, Vol 1, D Zucker-Franklin, MF Greaves, CE Grossi, AM Marmont, Editors, Lea and Febiger, Philadelphia, 1981, pp 149-242

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Referencer

- Graham RC, Karnovsky MJ: Glomerular permeability. Ultrastructural cytochemical studies using peroxidases as protein tracers. J Exp Med 124:1123, 1966
- Graham RC, Karnovsky MJ: The early stages of absorption of injected horseradish peroxidase in the proximal tubules of mouse kidney: Ultrastructural cytochemistry by a new technique. J Histochem Cytochem 14:291, 1966
- Hanker JS, Laszlo J, Moore JO: The light microscopic demonstration of hydroperoxidase-positive phi bodies and rods in leukocytes in acute myeloid leukemia. Histochemistry 58:241, 1978
- Hanker JS, Ambrose WW, James CJ, et al: Facilitated light microscopic cytochemical diagnosis of acute myelogenous leukemia. Cancer Res 39:1635, 1979
- Hanker JS, Chandross RJ, Weatherly NF et al: Medusa cells: The morphology and cytochemistry of common amoeboid variants of eosinophils. Histochem J 12:701, 1980
- Hanker JS, Giammara BL: Haematology: Some Advances and Insights from Histochemical and Molecular Biological Studies of Haemoproteins, IN Histochemistry: The Widening Horizons. PJ Stoward and JM Polak, Editors, John Wiley & Sons Ltd., Sussex, 1991
- Gill GW, Frost JK, Miller KA: A new formula for a half-oxidized hematoxylin solution that neither overstains nor requires differentiation. Acta Cytol 18:300, 1974
- Marmont AM, Damasio E, Zucker-Franklin D: Neutrophils. IN Atlas of Blood Cells - Function and Pathology, Vol 1, D Zucker-Franklin, MF Greaves, CE Grossi, AM Marmont, Editors, Lea and Febiger, Philadelphia, 1981, pp 149-242

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Referenser

- Graham RC, Karnovsky MJ: Glomerular permeability. Ultrastructural cytochemical studies using peroxidases as protein tracers. J Exp Med 124:1123, 1966
- Graham RC, Karnovsky MJ: The early stages of absorption of injected horseradish peroxidase in the proximal tubules of mouse kidney: Ultrastructural cytochemistry by a new technique. J Histochn Cytochem 14:291, 1966
- Hanker JS, Laszlo J, Moore JO: The light microscopic demonstration of hydroperoxidase-positive phi bodies and rods in leukocytes in acute myeloid leukemia. Histochemistry 58:241, 1978
- Hanker JS, Ambrose WW, James CJ, et al: Facilitated light microscopic cytochemical diagnosis of acute myelogenous leukemia. Cancer Res 39:1635, 1979
- Hanker JS, Chandross RJ, Weatherly NF et al: Medusa cells: The morphology and cytochemistry of common amoeboid variants of eosinophils. Histochem J 12:701, 1980
- Hanker JS, Giammara BL: Haematology: Some Advances and Insights from Histocytochemical and Molecular Biological Studies of Haemoproteins, IN Histochemistry: The Widening Horizons. PJ Stoward and JM Polak, Editors, John Wiley & Sons Ltd., Sussex, 1991
- Gill GW, Frost JK, Miller KA: A new formula for a half-oxidized hematoxylin solution that neither overstains nor requires differentiation. Acta Cytol 18:300, 1974
- Marmont AM, Damasio E, Zucker-Franklin D: Neutrophils. IN Atlas of Blood Cells - Function and Pathology, Vol 1, D Zucker-Franklin, MF Greaves, CE Grossi, AM Marmont, Redaktörer, Lea and Febiger, Philadelphia, 1981, sid. 149–242

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Referências

- Graham RC, Karnovsky MJ: Glomerular permeability. Ultrastructural cytochemical studies using peroxidases as protein tracers. J Exp Med 124:1123, 1966
- Graham RC, Karnovsky MJ: The early stages of absorption of injected horseradish peroxidase in the proximal tubules of mouse kidney: Ultrastructural cytochemistry by a new technique. J Histochem Cytochem 14:291, 1966
- Hanker JS, Laszlo J, Moore JO: The light microscopic demonstration of hydroperoxidase-positive phi bodies and rods in leukocytes in acute myeloid leukemia. Histochemistry 58:241, 1978
- Hanker JS, Ambrose WW, James CJ, et al: Facilitated light microscopic cytochemical diagnosis of acute myelogenous leukemia. Cancer Res 39:1635, 1979
- Hanker JS, Chandross RJ, Weatherly NF et al: Medusa cells: The morphology and cytochemistry of common amoeboid variants of eosinophils. Histochem J 12:701, 1980
- Hanker JS, Giammara BL: Haematology: Some Advances and Insights from Histochemical and Molecular Biological Studies of Haemoproteins, IN Histochemistry: The Widening Horizons. P.J. Stoward e JM Polak, Editores, John Wiley & Sons Ltd., Sussex, 1991
- Gill GW, Frost JK, Miller KA: A new formula for a half-oxidized hematoxylin solution that neither overstains nor requires differentiation. Ata Cytol 18:300, 1974
- Marmont AM, Damasio E, Zucker-Franklin D: Neutrophils. IN Atlas of Blood Cells - Function and Pathology, Vol 1, D Zucker-Franklin, MF Greaves, CE Grossi, AM Marmont, Editors, Lea and Febiger, Philadelphia, 1981, pp 149-242

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- Graham RC, Karnovsky MJ: Glomerular permeability. Ultrastructural cytochemical studies using peroxidases as protein tracers. J Exp Med 124:1123, 1966
- Graham RC, Karnovsky MJ: The early stages of absorption of injected horseradish peroxidase in the proximal tubules of mouse kidney: Ultrastructural cytochemistry by a new technique. J Histochem Cytochem 14:291, 1966
- Hanker JS, Laszlo J, Moore JO: The light microscopic demonstration of hydroperoxidase-positive phi bodies and rods in leukocytes in acute myeloid leukemia. Histochemistry 58:241, 1978
- Hanker JS, Ambrose WW, James CJ, et al: Facilitated light microscopic cytochemical diagnosis of acute myelogenous leukemia. Cancer Res 39:1635, 1979
- Hanker JS, Chandross RJ, Weatherly NF et al: Medusa cells: The morphology and cytochemistry of common amoeboid variants of eosinophils. Histochem J 12:701, 1980
- Hanker JS, Giammara BL: Haematology: Some Advances and Insights from Histochemistry and Molecular Biological Studies of Haemoproteins, IN Histochemistry: The Widening Horizons. P.J. Stoward and JM Polak, Editors, John Wiley & Sons Ltd., Sussex, 1991
- Gill GW, Frost JK, Miller KA: A new formula for a half-oxidized hematoxylin solution that neither overstains nor requires differentiation. Acta Cytol 18:300, 1974
- Marmont AM, Damasio E, Zucker-Franklin D: Neutrophils. IN Atlas of Blood Cells - Function and Pathology, Vol 1, D Zucker-Franklin, MF Greaves, CE Grossi, AM Marmont, Editors, Lea and Febiger, Philadelphia, 1981, pp 149-242

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Reference

- Graham RC, Karnovsky MJ: Glomerular permeability. Ultrastructural cytochemical studies using peroxidases as protein tracers. J Exp Med 124:1123, 1966
- Graham RC, Karnovsky MJ: The early stages of absorption of injected horseradish peroxidase in the proximal tubules of mouse kidney: Ultrastructural cytochemistry by a new technique. J Histochn Cytochem 14:291, 1966
- Hanker JS, Laszlo J, Moore JO: The light microscopic demonstration of hydroperoxidase-positive phi bodies and rods in leukocytes in acute myeloid leukemia. Histochemistry 58:241, 1978
- Hanker JS, Ambrose WW, James CJ, et al: Facilitated light microscopic cytochemical diagnosis of acute myelogenous leukemia. Cancer Res 39:1635, 1979
- Hanker JS, Chandross RJ, Weatherly NF et al: Medusa cells: The morphology and cytochemistry of common amoeboid variants of eosinophils. Histochem J 12:701, 1980
- Hanker JS, Giammara BL: Haematology: Some Advances and Insights from Histocytochemical and Molecular Biological Studies of Haemoproteins, IN Histochemistry: The Widening Horizons. P.J. Stoward and JM Polak, Editors, John Wiley & Sons Ltd., Sussex, 1991
- Gill GW, Frost JK, Miller KA: A new formula for a half-oxidized hematoxylin solution that neither overstains nor requires differentiation. Acta Cytol 18:300, 1974
- Marmont AM, Damasio E, Zucker-Franklin D: Neutrophils. IN Atlas of Blood Cells - Function and Pathology, Vol 1, D Zucker-Franklin, MF Greaves, CE Grossi, AM Marmont, Editors, Lea and Febiger, Philadelphia, 1981, pp 149-242

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Referenties

- Graham RC, Karnovsky MJ: Glomerular permeability. Ultrastructural cytochemical studies using peroxidases as protein tracers. J Exp Med 124:1123, 1966
- Graham RC, Karnovsky MJ: The early stages of absorption of injected horseradish peroxidase in the proximal tubules of mouse kidney: Ultrastructural cytochemistry by a new technique. J Histochem Cytochem 14:291, 1966
- Hanker JS, Laszlo J, Moore JO: The light microscopic demonstration of hydroperoxidase-positive phi bodies and in leukocytes in acute myeloid leukemia. Histochemistry 58:241, 1978
- Hanker JS, Ambrose WW, James CJ, et al: Facilitated light microscopic cytochemical diagnosis of acute myelogenous leukemia. Cancer Res 39:1635, 1979
- Hanker JS, Chandross RJ, Weatherly NF et al: Medusa cells: The morphology and cytochemistry of common amoeboid variants of eosinophils. Histochem J 12:701, 1980
- Hanker JS, Giammara BL: Haematology: Some Advances and Insights from Histocytochemical and Molecular Biological Studies of Haemoproteins, IN Histochemistry: The Widening Horizons. P.J. Stoward and JM Polak, Editors, John Wiley & Sons Ltd., Sussex, 1991
- Gill GW, Frost JK, Miller KA: A new formula for a half-oxidized hematoxylin solution that neither overstains nor requires differentiation. Acta Cytol 18:300, 1974
- Marmont AM, Damasio E, Zucker-Franklin D: Neutrophils. IN Atlas of Blood Cells - Function and Pathology, Vol 1, D Zucker-Franklin, MF Greaves, CE Grossi, AM Marmont, Editors, Lea and Febiger, Philadelphia, 1981, pp 149-242

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MDSS GmbH