3050 Spruce Street, St. Louis, MO 63103 USA
Tel: (800) 521-8956 (314) 771-5765 Fax: (800) 325-5052 (314) 771-5757
email: techservice@sial.com sigma-aldrich.com

Product Information

Histopaque[®]-1077 Hybri-Max™ BioReagent

Catalog Number H8889

Product Description

Histopaque-1077 is a solution of polysucrose and sodium diatrizoate, adjusted to a density of 1.077 +/- 0.001 g/ml. This medium facilitates rapid recovery of viable lymphocytes from small volumes of whole blood. This technique is suitable for use in cell-mediated lympholysis and for human lymphocyte antigen (HLA) typing. It may also be employed as the initial isolation step prior to enumeration of T-, B-, and 'null' lymphocytes.

Performance Characteristics

Anticoagulated venous blood is layered onto Histopaque-1077. During centrifugation, erythrocytes and granulocytes are aggregated by polysucrose and rapidly sediment. Lymphocytes and other mononuclear cells remain at the plasma-Histopaque-1077 interface. Erythrocyte contamination is negligible. Most extraneous platelets are removed by low speed centrifugation during the washing steps.

Product Information

Density: 1.076-1.078 g/ml Endotoxin: ≤ 3 EU/ml Solution pH: 8.7-9.1

Splenocyte separation: Clear interface with 90% cell

viability

Sterility: Sterile by USP XXIII

Reconstitution and Use

- To 3 ml whole blood, collected in heparin or EDTA, add 5 ml PBS without calcium and magnesium and mix well by inversion.
- To a 15 ml conical centrifuge tube, add 3 ml of Histopaque-1077 and bring to room temperature.
- Carefully layer 8 ml of the blood-saline mixture onto the Histopaque-1077. Centrifuge at 400 x g for EXACTLY 30 MINUTES at room temperature. Lower temps may cause clumping and poor recovery.

- 4. After centrifugation, use a Pasteur pipette to aspirate the upper layer to within 0.5 cm of the opaque interface containing the mononuclear cells. Discard the upper layer.
- 5. With a Pasteur pipette, carefully transfer the opaque interface to a clean, conical centrifuge tube.
- 6. Add to this tube 10 ml PBS and mix by inversion.
- 7. Centrifuge at 250 x g for **EXACTLY** 10 MINUTES.
- 8. Aspirate the supernatant and discard.
- 9. Resuspend lymphocyte pellet with 5 ml PBS and mix by gentle trituration with Pasteur pipette.
- 10. Centrifuge at 250 x g for **EXACTLY** 10 MINUTES.
- 11. Repeat steps 8,9,and 10, discard supernatant and resuspend pellet in 0.5 ml PBS.

Storage

Stored at 2-8 °C protected from light , Histopaque-1077 has a shelf-life of 2 years.

References

- Lightbody, J., Use of the cell-mediated lympholysis test in transplantation immunity, in Manual of Clinical Immunology, Rose, N.R. and Friedman, H. (Eds.) pp. 851-857 (American Society for Microbiology, Washington, D.C., 1976).
- Amos, D.B. and Pool, P., HLA typing Manual of Clinical Immunology, Rose, N.R. and Friedman, H., (Eds.) pp. 797-804 (American Society for Microbiology, Washington, D.C., 1976).
- Winchester, R.J., and Ross, G., Methods for enumerating lymphocyte populations in Manual of Clinical Immunology, Rose, N.R., and Friedman, H. (Eds.) pp. 64-76 (American Society for Microbiology, Washington, D.C., 1976).

Histopaque is a registered trademark of Sigma-Aldrich Co. LLC

Hybri-Max is a trademark of Sigma-Aldrich Co. LLC

DF,PHC 04/14-1