Immobilon™-Ny+ Blotting Membrane

Charged nylon blotting membrane for superior results from nucleic acid protocols



For superior results in genomics applications including Southern blotting, northern blotting, and colony/plaque lifts, use Immobilon Ny+ charged nylon blotting membrane.

Exceptional Signal and Reprobing Capabilities

Immobilon-Ny+ is a positively charged nylon membrane optimized for reliable and reproducible transfer, immobilization, hybridization, and subsequent reprobing. The density and uniformity of the positively charged surface produces maximum sensitivity with minimal background. The mem-

brane is reinforced for enhanced durability, which is especially important in reprobing experiments.

- Maximum sensitivity with minimal background
- Exceptional reprobing characteristics
- Performs well in both chemiluminescent and radioactive detection system

Performance Data

Excellent Reprobing Characteristics




Figure 1. Signal intensity on charged nylon after round 12 of reprobing: Nylon membranes were cross-linked at 254 nm and underwent multiple reprobes including, in sequence, two hybridizations with 32 P labeled probe, four mock hybridizations, a seventh hybridization with 32P-labeled probe followed by five mock hybridizations and a thirteenth hybridization with ³²P-labeled probe. The membranes were stripped in 0.4M NaOH at 45 °C for 30 minutes between each cycle of hybridization. The mock hybridizations were the same as the true hybridizations except that no 32P-labeled probe was added. Throughout the reprobing process, Immobilon-Ny+ performed better with twice the signal than both Competitor A and Competitor B's charged nylon membrane.

High Hybridization Signal




Figure 2. Direct comparison of three charged nylon membranes in their own optimized protocols. Each membrane was cross-linked at its recommended UV wavelength. Using these cross-linking conditions, Immobilon-Ny+ had twice the hybridization signal of both Competitor A and Competitor B.



Technical Literature

For more information including optimization tips and protocol notes, refer to the following (available online or through your local Millipore office):

- TN057 Comparison of Immobilon-Ny+ Transfer Membrane and 2 Competitive Charged Nylon **Membranes**
- TN054 Optimization of DNA Fixation on Immobilon-Ny+ Transfer Membrane Using UV Light
- -TN055 Transfer of DNA to Immobilon-Ny+ Transfer Membrane
- TN056 Comparison of 4 Stripping Protocols for DNA Probes on Immobilon-Ny+
- TN071 Chemiluminescent Detection
- -TN110 Optimization of Northern Blotting on Immobilion-Ny+
- TN1500ENUS Colony Lifts Using Immobilon-Ny+
- FF250EN00 Immobilion Blotting Membranes Brochure

Ordering Information

Immobilon-Ny+ Blotting Membranes are available in discs, sheets, and rolls.

Description	Qty/Pkg	Catalogue No.
82 mm disc	50	INYC 082 50
85 mm disc	50	INYC 085 50
132 mm disc	50	INYC 132 50
137 mm disc	50	INYC 137 50
6 x 11 sheet	50	INYC 611 50
9 x 12 sheet	50	INYC 091 50
10 x 10 cm sheet	10	INYC 101 00
20 x 20 cm sheet	10	INYC 202 00
22.2 x 22.2 cm sheet	10	INYC 222 2C
	100	INYC 262 60
26 x 26 cm sheet	10	INYC 262 60
$30 \text{ cm} \times 3.3 \text{ m roll}$	1	INYC 000 10

To Place an Order or Receive Technical Assistance

For additional information call your nearest Millipore office: In the U.S. and Canada, call toll-free 1-800-MILLIPORE (1-800-645-5476)In the U.S., Canada and Puerto Rico, fax orders to 1-800-MILLIFX (1-800-645-5439)

Internet: www.millipore.com E-mail: tech service@millipore.com

Millipore Worldwide

AUSTRALIA

Tel. 1 800 222 111 or (02) 9888 8999 Fax (02) 9878 0788

Tel. (01) 877-8926 Fax (01) 877-1654

BALTIC COUNTRIES

Tel. +358 9 804 5110 Fax +358 9 256 5660

BEIGIUM AND **IUXEMBOURG**

Tel. +32 2 726 88 40 Fax +32 2 726 98 84

Tel. (011) 5548-7011 Fax (011) 5548-7923

CANADA

Tel. 1-800-645-5476 Fax 1-800-645-5439

CHINA, PEOPLE'S REPUBLIC OF

(8610) 8519 1250 (8610) 6518 1058 Fax (8610) 8519 1255 Guangzhou:

Tel. (8620) 8755 4049 Fax (8620) 8752 0172

Hong Kong: Tel. (852) 2803 9111 Fax (852) 2513 0313 Shanghai:

Tel. (8621) 5306 9100 Fax (8621) 5306 0838 CZECH REPUBLIC

Tel. 02-2051 3841 Fax 02-2051 4298

DENMARK

Tel. 70 10 00 23 Fax 70 10 13 14

EASTERN EUROPE, C.I.S., AFRICA, MIDDLE EAST AND GULF

Life Sciences Division: Tel. +33 3.88.38.9536 Fax +33 3.88.38.9539 BioPharmaceutical Division:

Tel. +43 1 877-8926 Fax +43 1 877-1654 Laboratory Water Division: Tel. +33¹.30.12.70.00 Fax +33 1.30.12.71.80

FINI AND

Tel. (09) 804 5110 Fax (09) 256 5660

Tel. 01 30 12 7000 Fax 01 30 12 7180

GERMANY

Tel. (06196) 494-0 Fax (06196) 43901

HUNGARY Tel. 1-205 9784 Fax 1-205 9792

Tel. (91) 80-839 46 57 Fax (91) 80-839 63 45 ITALY

Vimodrone (Milano): Tel. (02) 25.07.81 Fax (02) 26.50.324 Roma:

Tel. (06) 52.03.600 Fax (06) 52.95.735

JAPAN

Tel. (03) 5442-9711 Fax (03) 5442-

9736 Life Sciences 9737 BioPharm. 9734 Lab Water

KOREA

Tel. (822) 551-0990 Fax (822) 551-0228

MALAYSIA

Tel. 603-757 1322 Fax 603-757 1711

MFXICO

Tel. (525) 576-9688 Fax (525) 576-8706

THE NETHERLANDS Tel 076-5022000

Fax 076-5022436 NORWAY

Tel. 22 67 82 53 Fax 22 66 04 60

POLAND

Tel. 22-669 12 25 22-663 70 31 Fax 22-663 70 33

PUERTO RICO

Tel. (787) 273-8495 Fax (787) 747-6553

SINGAPORE

Tel. (65) 842 1822 Fax (65) 842 4988

SPAIN AND PORTUGAL

Tel. (34) 917 283 960 Fax (34) 917 292 909

SWEDEN

Tel. 08-628 6960 Fax 08-628 6457

SWITZERLAND

Volketswil: Tel. (01) 908-30-60 Fax (01) 908-30-80 Lausanne:

Tel. (021) 641 2550 Fax (021) 641 2551

TAIWAN

Taipei: Tel. (886-2) 700-1742 Fax (886-2) 755-3267 Hsin Chu City: Tel. (886-3) 571-0178

Fax (886-3) 572-9520 U.K. AND IRFLAND

Tel. +44 1923 816375 Fax +44 1923 818297

Tel. (781) 533-6000 Fax (781) 533-3110

IN ALL OTHER COUNTRIES

Millipore Intertech (U.S.A.) Tel. +1 (781) 533-8622 Fax +1 (781) 533-8630

Printed with soy inks on recycled paper.

Copyright 2001 Millipore Corporation or an

MILLIPORE

Corporation or an affiliated company.

Printed in U.S.A. 4/01 01-098

or an affiliated company.

Lit. No. PF1031EN00

reserved

Millipore is a registered trademark of Millipore

Immobilon is a trademark of Millipore Corporation

affiliated company, Bedford, MA 01803. All rights