

Benzoic and Sorbic Acid in Milk Powder


Purospher® STAR RP-18 endcapped

Chromatographic Conditions

Column: Purospher® Star LP RP-18 endcapped (5 µm) 250x4.6 mm 1.56200.0008
 Injection: 10 µL
 Detection: UV@ 227 nm
 Cell: 1ul/10mm
 Flow Rate: 1.2 mL/min
 Mobile Phase: Phosphate buffer solution (pH6.7). Dissolve 2.5g KH₂PO₄ and 1.91g K₂HPO₄ in 1L H₂O.
 Mix Methanol and phosphate buffer solution 10:90 (v/v).
 Temperature: 30 °C
 Diluent: MeOH

Sample: Weigh 3.00 g of milk powder into a 100 mL volumetric flask. Add 10 mL water, 25 mL of 0.1 mol/L NaOH solution, sonicate for 15 minutes. After the solution is cooled down to room temperature, adjust pH value to 8.0 with 0.5 mol/L of H₂SO₄ solution. Add 2 mL of 92 g/L of K₄Fe(CN)₆ solution, 2 mL of 183g/L of ZnAC solution. Shakes vigorously the flask, cool the solution to room temperature. Add methanol to volume. Filtrate the solution with filter paper and 0.45 µm filter membrane respectively before analysis.

Pressure Drop: 192 Bar (2784 psi)



Chromatographic Data

No.	Compound	Retention Time (min)	Theoretical Plates	Tailing Factor
1	Benzoic acid	8.0	8323	1.05
2	Sorbic acid	10.9	8528	1.03

Benzoic and Sorbic Acid in Milk Powder


Purospher® STAR RP-18 endcapped

Repeatability: Analysis of repeated injections of 4 ppm standard solution (n=5).

Standard	Benzoic Acid	Sorbic Acid
1	97.9	141.0
2	97.3	139.5
3	97.3	139.1
4	97.4	139.6
5	97.6	139.8
Average	97.5	139.8
STDEV	0.219	0.655
RSD (%)	0.2	0.5

Linearity – Benzoic Acid: Analysis of repeated injections (n=5) of standard solutions in the calibration range 1–40 ppm (µg/mL).

Concentration (ppm)	Intensity (mV)
1.0	22.911
2.0	46.666
4.0	97.139
6.0	146.50
8.0	197.38
10.0	248.13
20.0	498.24
40.0	995.40



Benzoic and Sorbic Acid in Milk Powder

Purospher® STAR RP-18 endcapped

Repeatability: Analysis of repeated injections of 4 ppm standard solution (n=5).

Standard	Benzoic Acid	Sorbic Acid
1	97.9	141.0
2	97.3	139.5
3	97.3	139.1
4	97.4	139.6
5	97.6	139.8
Average	97.5	139.8
STDEV	0.219	0.655
RSD (%)	0.2	0.5

Linearity – Sorbic Acid: Analysis of repeated injections (n=5) of standard solutions in the calibration range 1–40 ppm (µg/mL).

Concentration (ppm)	Intensity (mV)
1.0	41.180
2.0	67.673
4.0	139.50
6.0	208.53
8.0	283.83
10.0	354.38
20.0	712.40
40.0	1422.6

