3050 Spruce Street, St. Louis, MO 63103 USA Tel: (800) 521-8956 (314) 771-5755 Fax: (800) 325-5052 (314) 771-5757 email: techservice@sial.com sigma-aldrich.com

# **Product Information**

# Thrombin from bovine plasma

Catalog Number **T6200** Storage Temperature –20 °C

CAS RN 9002-04-4 EC 3.4.21.5 Synonym: Factor IIa

# **Product Description**

Thrombin is an endolytic serine protease that selectively cleaves the Arg–Gly bonds of fibrinogen to form fibrin and release fibrinopeptides A and B.<sup>1,2</sup>

The predominant form of thrombin *in vivo* is the zymogen, prothrombin (factor II), which is produced in the liver. The concentration of prothrombin in normal human plasma is 5–10 mg/dL.<sup>3</sup> Prothrombin is a glycoprotein with a glycan content of ~12%.<sup>3</sup>

Prothrombin is cleaved *in vivo* by activated factor X, releasing the activation peptide and cleaving thrombin into light and heavy chains yielding catalytically active  $\alpha$ -thrombin.  $\alpha$ -Thrombin is composed of a light chain (A chain, MW ~6 kDA) and a heavy chain (B chain, MW ~31 kDa). These two chains are joined by one disulfide bond.<sup>4</sup> The B chain of human thrombin consists of a peptide portion (MW 29,485 Da) and a carbohydrate portion (MW 2,334 Da) with N-linked glycosylation at three Asn residues.<sup>5,6</sup> Bovine thrombin contains 1.7% glucosamine, 1.8% sialic acid, 0.61% galactose, and 0.95% mannose.<sup>7</sup>

Thrombin also contains  $\gamma$ -carboxyglutamyl residues. These modified glutamyl residues are the result of carboxylation by a microsomal enzyme, vitamin K-dependent carboxylase.  $\gamma$ -Carboxyglutamyl residues are necessary for the Ca<sup>2+</sup>-dependent interaction with a negatively charged phospholipid surface, which is essential for the conversion of prothrombin to thrombin.<sup>4</sup> Prothrombin is activated *in vivo* on the surface of a phospholipid membrane that binds the amino terminus of prothrombin along with factors Va and Xa. The activation process starts slowly because factor V is activated to factor Va by the initial, small amount of thrombin. Optimal cleavage sites for thrombin:<sup>2</sup>

- 1. A-B-Pro-Arg-||-X-Y where A and B are hydrophobic amino acids and X and Y are nonacidic amino acids
- 2. Gly-Arg-||-Gly

Thrombin from any mammalian species will clot the fibrinogen of any other mammalian species.<sup>8</sup>

Thrombin cleavage of fibrinogen occurs only at Arg residues; however, the cleavage site is not specific, resulting in 2 products. The primary cleavage product, fibrinopeptide A, is cleaved from fibrinogen after amino acid 16 and sometimes after amino acid 19, while a secondary cleavage product, fibrinopeptide B is produced by cleavage at amino acid 14.<sup>9</sup>

Thrombin does not require divalent metal ions or cofactors for activity. However, Na<sup>+</sup>-dependent allosteric activation of thrombin has been shown to play a role in defining the primary specificity of thrombin to cleave after Arg residues.<sup>10</sup> Thrombomodulin serves as a cofactor for thrombin during the activation of protein C.<sup>11</sup>

Under certain storage conditions, autolytic digestion of  $\alpha$ -thrombin results in formation of  $\beta$  and  $\gamma$ -thrombins, which lack fibrinolytic activity, but retain some activity against synthetic peptide substrates and protein substrates other than fibrinogen.<sup>12</sup> This thrombin preparation is predominantly  $\alpha$ -thrombin.

Thrombin (human and bovine) will catalyze the hydrolysis of several peptide *p*-nitroanilides, tosyl-Arg-nitrobenzyl ester, and thiobenzyl ester synthetic substrates.<sup>13</sup>

Catalytic pH range:<sup>14</sup> 5–10, optimal pH:<sup>14</sup> 8.3 thrombin precipitates ≤pH 5

Human isozymes pl range: 6.35–7.6. Bovine pl range:<sup>15</sup> 7.05–7.1

E<sup>1%</sup><sub>280</sub> = 18.3 (human)<sup>16</sup>

 $E_{280}^{1\%}$  = 19.5 (bovine)<sup>17</sup>

This product is lyophilized from a solution containing sodium chloride and Tris-HCl, pH 7.0. It is activated with bovine lung thromboplastin.

Specific Activity: 40-300 NIH units/mg protein  $(E_{280}^{1\%} = 19.5)$ 

Unit Definition: Activity is expressed in NIH units obtained by direct comparison to a NIH Thrombin Reference Standard, Lot K. The NIH assay procedure uses 0.2 ml of diluted plasma (1:1 with saline) as a substrate and 0.1 ml of albumin solution based on a modification of the method of Biggs.<sup>18</sup> Only clotting times in the range of 15-25 seconds are used for determining thrombin activity. Optimal clotting temperature is 37 °C.

Thrombin concentrations in the literature are typically reported in terms of different units of activity.<sup>18,19</sup> Several conventions are used to express thrombin activity in the literature:

> 1 IOWA unit = 0.83 NIH unit 1 WHO unit = 1 NIH unit 1 NIH unit =  $0.324 \pm 0.073 \,\mu g$ 1 NIH unit = 1 USP unit

# **Precautions and Disclaimer**

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

# **Preparation Instructions**

The product is soluble in water (10 mg/ml), yielding a clear solution.

## Storage/Stability

Stock solutions can be prepared at a concentration of 100 units/ml in a 0.1% (w/v) BSA solution. Stock solutions remain active for one week at 0-5 °C. Solutions are most stable at pH 6.5, as a pH >7 will greatly reduce thrombin activity. Since thrombin solutions adsorb to glass, it is recommended to aliquot the solutions in plastic tubes and store at -20 °C for long-term storage.

Store the lyophilized powder at -20 °C. The product retains activity for at least 5 years.

#### **Related Products**

Synthetic Substrates:13 N-Benzoyl-Phe-Val-Arg-p-nitroanilide (Catalog Number 13042) N-Benzoyl-Phe-Val-Arg-p-nitroanilide hydrochloride (Catalog Number B7632) N-Benzoyl-Phe-Val-Arg 4-methoxy-B-naphthylamide (Catalog Number B1260) Boc-β-benzyl-Asp-Pro-Arg-7-amido-4-methylcoumarin (Catalog Number B4028) Boc-Val-Pro-Arg-7-amido-4-methylcoumarin (Catalog Number B9385) Sar-Pro-Arg p-nitroanilide (Catalog Number S9009) Thrombin generation chromogenic substrate (Catalog Number T3068) N-p-Tosyl-Gly-Pro-Arg 7-amido-4-methylcoumarin (Catalog Number T0273) N-(p-Tosyl)-Gly-Pro-Arg p-nitroanilide (Catalog Number T1637) Inhibitors:20-22 Diisopropylfluorophosphate (Catalog Number D0879) Phenylmethylsulfonyl fluoride (Catalog Number P7626) AEBSF (Catalog Number A8456) Hirudin (Catalog Number H7016) Proflavine (Catalog Number P2508) Antithrombin III (Catalog Number A2221)  $\alpha_1$ -antitrypsin (Catalog Number A9024)  $\alpha_1$ -antiplasmin (Catalog Number A8849) Gabexate mesylate (Catalog Number G2417) Antipain (Catalog Number A6191)  $N_{\alpha}$ -Tosyl-L-lysine chloromethyl ketone hydrochloride (Catalog Number T7254)

## References

- 1. Enzyme Nomenclature: EC 3.4.21.5
- 2. Chang, J.Y., Eur. J. Biochem., 151, 217-224 (1985).
- The Plasma Proteins, 2nd ed., 2, Putnam, F. W., 3. ed: Table 2. See also: The Enzyme Explorer: Plasma and Blood Protein Resource
- Expasy/SwissProt: P00743 4.
- 5. Qian, W.J., et al., J. Proteome Res., 4, 2070-2080 (2005).
- 6. Nilsson, B., et al., Arch. Biochem. Biophys., 224, 127-133 (1983).
- Boyer, P.D., The Enzymes, Academic Press (New 7. York), 3rd ed., Vol. III, p. 277-321 (1971). *The Plasma Proteins*, 2<sup>nd</sup> ed., **2**, Putnam, F. W., ed,
- 8. p. 148.
- 9. Machovich, R., The Thrombin, 1, 63-66 (1984).
- 10. Prasad, S., J. Biol. Chem., 279, 10103-10108 (2004).
- 11. Kisiel, W., Human plasma protein C: isolation, characterization, and mechanism of activation by alpha-thrombin. J. Clin. Invest., 64, 761-769, (1979).
- 12. Boissel, J.P., et al., J. Biol. Chem., 259, 5691-5697 1984).
- 13. Lottenberg, R., et al., Assay of Coagulation Proteases Using Peptide Chromogenic and Fluorogenic Substrates. Meth. Enzymol., 80-C, 341-361 (1981).
- 14. Machovich, R., The Thrombin, 1, 111 (1984).
- 15. Righetti, P.G., and Tudor, G., Isoelectric points and molecular weights of proteins, a new table. Journal of Chromatography, 220, 115-194 (1981).

- 16. Butkowski, R.J. et al., J. Biol. Chem., 252, 4942 (1977).
- 17. Winzor, D.J., and Scheraga, H.A., Arch. Biochem. Biophys., 104, 202-207 (1964).
- 18. Biggs, R., ed., Human Blood Coagulation, Haemostasis and Thrombosis 2<sup>nd</sup> ed., Blackwell Scientific Publications (Philadelphia: 1976), 722.
- 19. The Handbook of Synthetic Substrates, Hemker, H.C., Martinus Nijhoff publisher (1983).
- 20. Lundblad, R.L., et al., Methods Enzymol., 45, 156 (1976).
- 21. Matsuoka, S., et al., JP. J. Pharmacol., 51, 455-463 (1989).
- 22. Wimen, B., Meth. Enzymol., 80, 395-408 (1981).
- 23. Human Blood Coagulation, Haemostasis and Thrombosis, 2nd ed., R. Biggs, ed., p. 722 (1976).
- 24. Chang, Y., Thrombin specificity. Requirement for apolar amino acids adjacent to the thrombin cleavage site of polypeptide substrate. Eur. J. Biochem., 151(2), 217-224 (1985).
- 25. Hakes, D.J., and Dixon, J.E., Anal. Biochem., 202, 293 (1992).
- 26. Gaun, K.L., and Dixon, J.E., Anal. Biochem., 192, 262 (1991).
- 27. De Cristofaro, R., and De Candia, E., J. Thromb. Thrombolysis, 15, 151-163 (2003).
- 28. Sherwood, J.A., Mol. Biochem. Parisitol., 40, 173-181 (1990).
- 29. Berg, D.T., et al., Sceince, 273, 1389-1391 (1996).
- 30. Magnusson, S., The Enzymes, 3rd ed., III, pp. 277-321, Boyer, P.D., ed., Academic Press (1971).

GY,RC,PHC 01/14-1

©2014 Sigma-Aldrich Co. LLC. All rights reserved. SIGMA-ALDRICH is a trademark of Sigma-Aldrich Co. LLC, registered in the US and other countries. Sigma brand products are sold through Sigma-Aldrich, Inc. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see product information on the Sigma-Aldrich website at www.sigmaaldrich.com and/or on the reverse side of the invoice or packing slip.