



3050 Spruce Street  
Saint Louis, Missouri 63103 USA  
Telephone 800-325-5832 • (314) 771-5765  
Fax (314) 286-7828  
email: techserv@sial.com  
sigma-aldrich.com

## Product Information

### Protease from *Streptomyces griseus*

Catalog Number **P6911**  
Storage Temperature  $-20^{\circ}\text{C}$

CAS RN 9036-06-0  
Synonyms: Actinase E; Pronase E

#### Product Description

Protease from *Streptomyces griseus* is a mixture of at least three proteolytic activities including an extracellular serine protease. In general, serine proteases display a wide range of substrate specificities, which are believed to be mediated by an active site composed of one Asp, one His, and a Ser residue in the molecule. This enzyme prefers to hydrolyze peptide bonds on the carboxyl side of glutamic or aspartic acid.

Molecular mass:  $\sim 20$  kDa (monomer)

Protease activity can be completely inactivated by heating above  $80^{\circ}\text{C}$  for 15–20 minutes.

This product is collected from the culture broth of *S. griseus* and purified by successive column procedures. It is highly stable in the pH range of 5.0–9.0, with peak activity at pH 8.8. This protease preparation is a molecular biology reagent with no DNase, RNase, or nickase activity detected. Protease is typically used in nucleic acid isolation procedures in incubations of 0.5–3.0 hours supplemented with 0.2% sodium dodecyl sulfate and 10 mM EDTA. It is compatible with many DNA and RNA isolation buffers and has broad substrate specificity.

Activity:  $\geq 4$  units/mg solid

Unit definition: One unit will hydrolyze casein to produce color equivalent to 1.0 micromole (181  $\mu\text{g}$ ) of tyrosine per minute at pH 7.5 at  $37^{\circ}\text{C}$ .

The protease is incubated for 10 minutes at pH 7.5 at  $37^{\circ}\text{C}$  in a 6 ml reaction volume containing 0.54% casein and 0.041 M potassium phosphate buffer. The reaction is stopped by the addition of 5.0 ml of 0.11 M trichloroacetic acid. Liberated tyrosine equivalents were determined using Folin-Ciocalteu reagent.

#### Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

#### Preparation Instructions

The product is soluble in water (5–20 mg/ml).

#### Storage/Stability

Store the product at  $-20^{\circ}\text{C}$ .

#### References

1. Maniatis, T., *et al.*, Molecular Cloning, A Laboratory Manual, Cold Spring Harbor Laboratory (1982).

JB,KH,MAM 05/07-1

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.