

Product Information

R³ IGF-I, Human

Recombinant analog expressed in *E. coli*

CAS RN 139659-92-0

Synonym: Insulin-like growth factor I

I1146

Product Description

R³ IGF-I is a recombinant analog of insulin-like growth factor containing the complete human IGF-I amino acid sequence with substitution of Arg for Glu³. It was developed as an inexpensive, high quality potent analog of IGF for use as a growth factor supplement for serum-free or low serum cell culture. Recombinantly produced in *E. coli* using a patented expression system, the peptide is harvested from cells in inclusion bodies, which are dissolved and desalted. The active molecule is properly folded under oxidizing conditions and then purified by several liquid chromatography steps. The product is lyophilized from 0.1 M acetic acid and is shown to proliferate Chinese hamster ovary CHO cells at an ED₅₀ of 10 ng/mL.

Molecular mass: 7,676 Da

Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Safety Data Sheet for information regarding hazards and safe handling practices.

Preparation Instructions

Prepare a stock solution of 1 mg/mL by adding 0.1 mL of 10 mM HCl to the vial. Stock solutions of peptide can be stored for at least 3 months at -20 °C or -70 °C. Avoid repeated freeze-thaw cycles.

For concentrations of < 1 mg/mL in buffer, a carrier protein should be added to minimize adsorption of the peptide to plastic or glass surfaces. Bovine serum albumin (BSA) at concentrations between 1-5 mg/mL is recommended as a carrier protein in serum-free conditions. A carrier protein is not necessary in serum-supplemented media. Up to 400 mL of buffer solution containing BSA, or serum-supplemented medium may be added to the stock solution in the vial and mixed. Long-term storage of solutions of < 1 mg/mL is not recommended.

Additional Handling Suggestions

- Do not add the peptide to low protein or protein free media prior to filter sterilization
- Use a low protein binding membrane for filter sterilization
- Filter sterilize solutions at a concentration of 1 mg/mL or greater if no carrier protein is present

Storage and Stability

Store the product at 2 to 8 °C.

Notice

We provide information and advice to our customers on application technologies and regulatory matters to the best of our knowledge and ability, but without obligation or liability. Existing laws and regulations are to be observed in all cases by our customers. This also applies in respect to any rights of third parties. Our information and advice do not relieve our customers of their own responsibility for checking the suitability of our products for the envisaged purpose.

The information in this document is subject to change without notice and should not be construed as a commitment by the manufacturing or selling entity, or an affiliate. We assume no responsibility for any errors that may appear in this document.

Technical Assistance

Visit the tech service page at SigmaAldrich.com/techservice.

Terms and Conditions of Sale

Warranty, use restrictions, and other conditions of sale may be found at SigmaAldrich.com/terms.

Contact Information

For the location of the office nearest you, go to SigmaAldrich.com/offices.

The life science business of Merck operates
as MilliporeSigma in the U.S. and Canada.

Merck and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

© 2010-2024 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.

I1146pis Rev 09/24

