



Product Information

4-Methylumbelliferyl β -D-glucuronide hydrate

Product Number **M 9130**
Storage Temperature -0 °C

Product Description

Molecular Formula: $C_{16}H_{16}O_9$
Molecular Weight: 352.2 (anhydrous)
CAS Number: 6160-80-1
Melting Point: 97-100 °C
Specific Rotation: -105° (0.25% (w/v) in water)
Excitation Maximum: 365 nm¹
Emission maximum: 445 nm¹
Synonym: MUG

β -Glucuronidase (GUS) from *E. coli* has become the reporter enzyme of choice for genetic plant research. 4-Methylumbelliferyl β -D-glucuronide (MUG) is commonly used as a substrate for detecting GUS gene expression in plants.² Since the GUS gene encodes an enzyme not found in plants, this system can be very useful in identifying transformed plants. 4-Methylumbelliferyl β -D-glucuronide (MUG) shows little or no fluorescence. However, when treated with β -glucuronidase, the 4-methylumbelliferone product is fluorescent. The fluorescence at pH 10.3 is approximately 100 times as intense as at pH 7.0 (0.15 M glycine buffer).¹ These fluorescent properties allow MUG to be utilized as a very effective substrate for GUS.^{3,4}

MUG is also used for identifying *E. coli* contamination in drinking water⁵ and for rapid bacterial identification in blood cultures.⁶

Precautions and Disclaimer

For Laboratory Use Only. Not for drug, household or other uses.

Preparation Instructions

This product is soluble in water at a concentration of 0.35 mg/ml. It has also been described as being soluble at 0.1 mM in 0.1 M Sorensen's buffer.⁶

Storage/Stability

MUG at 100 mg/100 ml of H₂O plus 2 drops TRITON[®] X-100, which has been sterile filtered, is stable for 6 months when stored refrigerated.⁷ Solutions of MUG may also be sterilized by autoclaving at 121 °C for 15 minutes.⁸

References

1. J. Org. Chem., **27**, 1074 (1962).
2. Plant Mol. Biol. Rep., **5**, 387-405 (1987).
3. Biochem. J., **61**, 569 (1955).
4. Brot, F.E., et al., Purification and properties of beta-glucuronidase from human placenta. Biochemistry, **17**, 385-391 (1978).
5. Edberg, S.C., et al., National field evaluation of a defined substrate method for the simultaneous detection of total coliforms and *Escherichia coli* from drinking water: comparison with presence-absence techniques. Appl. Environ. Microbiol., **55**, 1003-1008 (1989).
6. Sepulveda, J.L. et al., Rapid presumptive identification of gram-negative rods directly from blood cultures by simple enzymatic tests. J. Clin. Microbiol., **28**, 177-181 (1990).
7. Thompson, J.S., et al., Rapid biochemical test to identify verocytotoxin-positive strains of *Escherichia coli* serotype O157., J. Clin. Microbiol., **28**, 2165-2168 (1990).
8. Feng, P.C. and Hartman, P.A., Fluorogenic assays for immediate confirmation of *Escherichia coli*. Appl. Environ. Microbiol., **43**, 1320-1329 (1982).

TRITON is a registered trademark of Union Carbide.

TMG/JRC 9/02

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.