

3050 Spruce Street
Saint Louis, Missouri 63103 USA
Telephone 800-325-5832 • (314) 771-5765
Fax (314) 286-7828
email: techserv@sial.com
sigma-aldrich.com

# **ProductInformation**

## Wortmannin from Penicillium funiculosum

Catalog Number **W1628** Storage Temperature 2-8 °C

CAS RN: 19545-26-7

Synonym: KY 12420

Molecular formula: C<sub>23</sub>H<sub>24</sub>O<sub>8</sub>

Molecular weight: 428.43 (anhydrous)

Melting point: 240 °C 1

#### **Product Description**

Wortmannin is a highly cell permeable antifungal antibiotic isolated from *Penicillium funiculosum*. It has no antibacterial properties and is specifically active against fungi.<sup>2</sup>

It is a specific and potent inhibitor of myosin light chain kinase <sup>3</sup> and a potent inhibitor of neutrophil activation by inhibiting F-Met-Leu-Phe(FMLP)-stimulated superoxide anion production without affecting intracellular calcium mobilization. <sup>4</sup> It inhibits FMLP-stimulated phospholipase D activation without direct inhibition of the enzyme. <sup>5</sup> It affects various signal transduction cascades by inhibiting phosphatidylinositol-3-kinases (PI3Ks) and blocking the recovery of calcium and potassium channels. This specific activity is important in studies of various receptor-mediated mechanisms. <sup>6-8</sup>

Phosphatidylinositol-3-kinase participates in the signal transduction pathway responsible for histamine secretion following stimulation of high affinity immunoglobulin E receptor (FceRI). Wortmannin blocks these responses through direct interaction with the catalytic subunits (110 kDa) of PI3-kinase.<sup>9</sup>

Wortmannin inhibited the activity of partially purified PI3-kinase from calf thymus at concentrations as low as 1.0 nM and with IC $_{50}$  values of 3.0 nM. Inhibition was irreversible. It inhibited both FceRI-mediated histamine secretion and leukotriene release up to 80% with IC $_{50}$  values of 2.0 and 3.0 nM, respectively.

Additional activities: immunosuppressive activity  $^{10}$ . strong anti-inflammatory activity  $^{11}$ , suppression of cellular responses such as respiratory burst and exocytosis in neutrophils,  $^5$  and catecholamine release in adrenal chromaffin cells $^{12}$  Aggregation and serotonin release in platelets were reported using a final concentration of 1  $\mu$ M of wortmannin in 0.01% DMSO.  $^{13,14}$ 

#### **Precautions and Disclaimer**

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

#### **Preparation instructions**

Sigma routinely tests the solubility at 10 mg/ml in ethyl acetate yielding a clear very faint yellow solution. Also soluble in:

Methanol (5 mg/mL)

Methanol (5 mg/mL)
DMSO (50 mg/mL)

#### Storage/Stability

Store the product desiccated and protected from light at 2-8  $^{\circ}$ C. Under these conditions the product is stable for 3 years.

It is unstable in aqueous solutions of pH 3-8. <sup>1</sup> A solution of 50 mg/ml in DMSO is stable for at least 2 months if stored at -20 °C.

### References:

- 1. Merck Index, 13<sup>th</sup> ed., No. 10110 (2001).
- 2. Korzybski, T., et al., *Antibiotics,* **III**, 1797 (1978).
- 3. Nakanishi, S., et al., *J. Biol. Chem.*, **267**, 2157-63 (1992).
- 4. Okada T., et al., *J. Biol. Chem.*, **269**, 3563-67 (1994).

- 5. Reinhold, S.L., et al., FASEB J., 4, 208-214 (1990).
- 6. Xie, L.H., et al., *J. Physiol.*, **514**, 655-65 (1999).
- 7. Cardenas, M.E., et al., *Trends Biotechnol.*, **16**, 427-33 (1998).
- 8. Stein, R.C., *Endocr. Relat. Cancer*, **8**, 237-248 (2001).
- 9. Yano H., et al., *J. Biol. Chem.*, **268**, 25846-25856 (1993).
- 10. Gunther, R., et al., *Immunopharmacology and Immunotoxicology*, **11**, 559-570 (1989).

- 11. Wiesinger, D., et al., *Experientia*, **30**, 135-136 (1974).
- 12. Ohara-Imaizumi, M., et al., *Biochem. Biophys. Res. Commun.*, **185**, 1016-1021 (1992).
- 13. Hashimoto, Y., et al., *J. Biol. Chem.*, **267**, 17078-17081 (1992).
- 14. Yatomi, Y., et al., *Biochem. J.*, **285**, 745-751 (1992).

ES-S,PHC 11/06-1