

Product Information

MONOCLONAL ANTI-CYCLIN D₃

Clone DCS-22

Mouse Ascites fluid

Product Number **C 7214**

Product Description

Monoclonal Anti-Cyclin D₃ (mouse IgG1 isotype) is derived from the DCS-22 hybridoma produced by the fusion of mouse myeloma cells and splenocytes from an immunized mouse. Recombinant cyclin D₃ protein of human origin was used as the immunogen.¹ The isotype is determined using Sigma ImmunoType™ Kit (Product Code ISO-1) and by a double diffusion immunoassay using Mouse Monoclonal Antibody Isotyping Reagents (Product Code ISO-2).

During the cell cycle of most somatic cells, DNA synthesis (S-phase) and mitosis (M-phase) are separated by two gap phases (G₁ and G₂) of varying duration. Thus, a typical eukaryotic cell sequentially passes through G₁, S, G₂, and M and back into G₁ during a single cycle.³ Regulation of cell cycle progression in eukaryotic cells depends on the expression of cyclin proteins.⁴ These proteins are the regulatory subunits of the cyclin dependent kinases (cdks), which are responsible for the phosphorylation of several cellular targets. Complexes of cyclins and cdks play a key role in cell cycle control. Within the complexes, the cyclin subunit serves a regulatory role, whereas cdks have a catalytic protein kinase activity.⁵ Members of the cyclin family of proteins combine with a cdk subunit to form the active kinase, which initiates G₂ to M and G₁ to S transition. The latter is controlled by cyclins termed G₁ cyclins, which commit the cell to DNA replication. Therefore, the cell cycle can be considered as a cyclin cycle which is controlled by biochemical modifications and formation of complex(es) with cdks.⁶ At least five candidate G₁-phase cyclins, termed cyclins C, D₁, D₂, D₃, and E have been identified in mammalian cells. Each of these cyclins can associate with one or more of the cdk family members. D-type cyclins are induced during the G₁ phase of the mammalian cell cycle in response to a variety of mitogenic growth factors. The three distinct members of the D-type cyclin family are differentially and combinatorially expressed in various cell lineages. Once induced, the D-type cyclins accumulate in complexes with cdks, whose kinase activity may be necessary for driving cells into S phase. The major catalytic partners of the D-type cyclins are cdk4 and cdk6, but at least some D-type cyclins also interact with other cdks, including cdk2 and

cdk5. Cyclin D₁- and D₂-associated cdk4 and/or cdk6 kinase activities have been detected in mid-G₁, prior to the activation of any other known cdk, and they culminate in late G₁ phase. The cyclin D₃-associated cdk4 and/or cdk6 exhibit kinase activities at the G₁/S transition. Cyclins D₁, D₂, and D₃ can be distinguished by their slightly different mobilities on denaturing gels. Under these conditions, the apparent masses are 36, 33-35, and 31-34 kD for cyclins D₁, D₂, and D₃, respectively. Because D-type cyclins probably serve as integrators of growth factor-induced signals with the cell cycle clock, aberrant expression of these proteins might play a role in disrupting the normal timing of events governing G₁ progression and contribute to oncogenesis. A link between tumor formation and inappropriate expression of cyclins has been established.^{7,8} The availability of monoclonal antibody reacting specifically with cyclin D₃ enables the subcellular detection and localization of cyclin D₃ and the measurement of relative differences in cyclin D₃ levels as a function of cell cycle phase.

Reagents

The product is provided as ascites fluid with 0.1% sodium azide as a preservative.

Precautions and Disclaimer

Due to the sodium azide content a material safety sheet (MSDS) for this product has been sent to the attention of the safety officer of your institution. Consult the MSDS for information regarding hazardous and safe handling practices.

Storage/Stability

For continuous use, store at 2-8 °C for up to one month. For extended storage, freeze in working aliquots. Repeated freezing and thawing is **not** recommended. Storage in "frost-free" freezers is **not** recommended. If slight turbidity occurs upon prolonged storage, clarify the solution by centrifugation before use.

Product Profile

Monoclonal Anti-Cyclin D₃ recognizes cyclin D₃ (31-34 kDa) applying ELISA, immunoblotting,^{1,2} immunocytochemistry,¹ immunoprecipitation,¹ and immunohistochemistry¹ (methacarn and formalin-fixed, paraffin-embedded tissue, upon unmasking). The product does not cross-react with other D-type cyclins.¹ Cross-reactivity has been observed with human,^{1,2} monkey, dog, rat, and mouse cyclin D₃.

Monoclonal Anti-Cyclin D₃ may be used for the localization of Cyclin D₃, using various immunochemical assays including immunoblotting,^{1,2} ELISA, immunocytochemistry,¹ immunoprecipitation,¹ and immunohistochemistry.¹

In order to obtain best results, it is recommended that each user determine the optimal working dilution for individual applications by titration assay.

References

1. Bartkova, J., et al., *Int. J. Cancer*, **65**, 323 (1996).
2. Bartkova, J., et al., *Cancer Res.*, **55**, 949 (1995).
3. Freeman, R. S., and Donoghue, D., *Biochemistry*, **30**, 2293 (1991).
4. Pines, J., and Hunter, T., *J. Cell Biol.*, **115**, 1 (1991).
5. Yamashita, M., et al., *Dev. Growth Differ.*, **33**, 617 (1991).
6. Norbury, C., and Nurse, P., *Ann. Rev. Biochem.*, **61**, 441 (1992).
7. Sherr, C. J., *Cell*, **73**, 1059 (1993).
8. Hunter, T., and Pines, J., *Cell*, **79**, 573 (1994).

JWM/KMR 02/02

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.