

## Product Information

# 3,3',5,5'-Tetramethylbenzidine (TMB) Liquid Substrate System for Membranes

Ready-to-use solution

**T0565**

## Product Description

3,3',5,5'-Tetramethylbenzidine (TMB) is a chromogenic substrate suitable for use in colorimetric procedures which utilize horseradish peroxidase (HRP) conjugates.<sup>1-4</sup> These procedures include ELISA and immunoblotting.

This product is supplied as a one-component, ready-to-use HRP substrate that contains TMB in a mildly acidic buffer. It is optimized for use with either nitrocellulose or PVDF Western blotted (transferred and probed) membranes. This solution develops a dark blue, insoluble reaction product in the presence of HRP. Prior to reaction with HRP, the solution is colorless to light yellow. The TMB solution has been optimized to reduce membrane background. Due to the insoluble reaction product, the TMB solution is not compatible with ELISA applications.

This TMB solution has similar sensitivity levels to chemiluminescent detection reagents, as it can detect as little as 0.15 ng. The TMB solution provided is enough for at least 25 mini-gel sized (10 × 10 cm) blots.

Several publications,<sup>5-8</sup> theses,<sup>9</sup> and dissertations<sup>10-16</sup> cite use of this product in their research protocols.

## Precautions and Disclaimer

This product is for R&D use only. Not for drug, household, or other uses. Please consult the Safety Data Sheet for information regarding hazards and safe handling practices.

## Procedure

This procedure assumes that the Western blotting (transfer and probing) steps have already occurred and that HRP-conjugates are bound to the blot.

1. Place the membrane on a clean flat sheet of plastic wrap.
2. Use enough TMB solution to completely cover the membrane's surface. Typically 3 mL is enough to cover a mini-gel sized (10 × 10 cm) membrane.
3. Expose the membrane to the TMB solution at room temperature for 5-15 minutes. Visually monitor the reaction.
4. Remove the substrate when protein bands are visible and the background is still low. High background diminishes the contrast between positive signal and background.
5. Wash the membrane in ultrapure water (such as Cat. No. W4502) for 1 minute.
6. Capture the image of the wet membrane using a camera or scanner.
7. Store the membrane dry and in the dark. If stored correctly, signal should remain on the membrane for a week.

## Storage/Stability

The product as supplied is stable for at least 1 year if stored at 2-8 °C. Keep the substrate out of direct sunlight.

## References

1. Bos, E. *et al.*, *J. Immunoassay*, **2(3-4)**, 187-204 (1981).
2. Wróblewska, B. *et al.*, *Int. J. Food Sci. Tech*, **39(8)**, 839-850 (2004).
3. Doig, N. M. *et al.* *J. Neurosci.*, **30(44)**, 14610-14618 (2010).
4. Szymkiewicz, A., and Chudzik-Kozłowska, J., *Acta Alimentaria*, **43(2)**, 193-291 (2014).
5. Kwok, D. T. K. *et al.*, *Surface Coatings Technol.*, **204(18-19)**, 2892-2897 (2010).
6. Moriarty, T. J. *et al.*, *Mol. Microbiol.*, **86(5)**, 1116-1131 (2012).
7. Bandy, N. J. *et al.*, *J. Bacteriol.*, **196(13)**, 2396-2404 (2014).
8. Rascón-Castelo, E. *et al.*, *Vaccines*, **3(4)**, 973-987 (2015).
9. Maema, Obakeng Goitseone, "Expression feasibility of recombinant enterokinase in *Nicotiana benthamiana*". University of Pretoria, M.Sc. thesis, p. 50 (2014).
10. Xu, Xing, "Strategies for recombinant protein production in maize". Iowa State University, Ph.D. dissertation, pp. 73, 98 (2012).
11. Burkhart, Annette, "The Blood-Brain Barrier *in vitro* Using Primary Culture: Implications for Studies of Therapeutic Gene Expression and Iron Transport". Aalborg University, Ph.D. dissertation, p. 67 (2014).
12. Pruvot, Mathieu, "Eco-epidemiology of production limiting diseases at the wildlife-livestock interface: beef cattle and elk in southwestern Alberta, Canada". University of Calgary, Ph.D. dissertation, p. 293 (2014).
13. Adihkary, Partho Protim, "Role of Progesterone Receptor Membrane Component 1 (PGRMC1) in Cancer Cell Biology". Charles Sturt University, Ph.D. dissertation, p. 82 (2015).
14. Busin, Valentina, "The development of microfluidic paper-based analytical devices for point-of-care diagnosis of sheep scab". Heriot-Watt University, Ph.D. dissertation, p. 89 (2017).
15. Healey, Robert David, "The Supramolecular Chemistry of Natural Sweeteners". University of New South Wales, Ph.D. dissertation, p. 3-49 (2017).
16. Tang, Qingyun, "Engineering of Methyltransferases: Exploring Their Catalytic Promiscuity". Universität Greifswald, Dr. rer. nat. dissertation, p. 31 (2020).

## Notice

We provide information and advice to our customers on application technologies and regulatory matters to the best of our knowledge and ability, but without obligation or liability. Existing laws and regulations are to be observed in all cases by our customers. This also applies in respect to any rights of third parties. Our information and advice do not relieve our customers of their own responsibility for checking the suitability of our products for the envisaged purpose.

The information in this document is subject to change without notice and should not be construed as a commitment by the manufacturing or selling entity, or an affiliate. We assume no responsibility for any errors that may appear in this document.

## Technical Assistance

Visit the tech service page at [SigmaAldrich.com/techservice](https://www.sigmaaldrich.com/techservice).

## Terms and Conditions of Sale

Warranty, use restrictions, and other conditions of sale may be found at [SigmaAldrich.com/terms](https://www.sigmaaldrich.com/terms).

## Contact Information

For the location of the office nearest you, go to [SigmaAldrich.com/offices](https://www.sigmaaldrich.com/offices).

## Troubleshooting Guide

Below are some common problems and their corresponding solutions associated with Western Blotting detection using TMB substrate.

<b>Problem</b>	<b>Cause</b>	<b>Solution</b>
Too much background signal observed on membrane	TMB substrate was left on the membrane too long	Decrease the amount of time that the TMB substrate is on the membrane.
	Too much primary antibody used	Decrease the amount of primary antibody used. Wash with TBST for 5 minutes after the primary antibody incubation.
Nonspecific bands show up on the membrane	Too much secondary antibody used	Decrease the amount of secondary antibody used.
	Too much primary antibody used	Decrease the amount of primary antibody used. Wash the membrane with TBST for 5 minutes after primary antibody incubation.
Signal disappears from membrane	Too much secondary antibody used	Decrease the amount of secondary antibody used
	Membrane not stored correctly	Store the membrane in the dark in ultrapure water.
No signal is observed on the membrane	Signal Degrades over time	Signal will degrade after a week even if membrane is stored in the dark, capture image with a camera or scanner.
	Low amounts of specific protein present	Expose the membrane to TMB substrate for a longer period of time. Include a positive control(s) during analysis.
	Insufficient primary antibody used	Use more primary antibody.
Protein degraded into proteolytic fragments	Insufficient secondary antibody used	Use more secondary antibody.
		Add protease inhibitors to original sample before running a gel.

The life science business of Merck KGaA, Darmstadt, Germany operates as MilliporeSigma in the U.S. and Canada.

MilliporeSigma, and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.  
© 2002-2023 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.

T0565bul Rev 09/23 JDS,ALC,GCY,MAM

