User Guide

Cathivex®-GV Filter Unit (Sterile)

REF SLGV0250S*

CE ₀₄₅₉ **R_x**_{only}

Available languages

English	3	Lietuviškai (Lithuanian).	23
Français (French).	5	Magyar (Hungarian)	25
Italiano (Italian)	7	Český (Czech)	27
Deutsch (German)	9	Polski (Polish)	29
Español (Spanish)	11	繁體中文 (Traditional Chinese).	32
Português (Portuguese)	13	Русский (Russian)	34
Nederlands (Dutch)	15	Norsk (Norwegian).	37
Dansk (Danish)	17	Türkçe (Turkish)	39
Svenska (Swedish).	19	Slovenčina (Slovak)	41
Latviski (Latvian)	21		

*Brazil Product Registration
see page 2.

PR05874w
Rev 06/23
Date of issue: 2023-06

[SigmaAldrich.com](https://www.SigmaAldrich.com)

Millipore®

Definition of Symbols Used

Symbol	Definition	Symbol	Definition
	Sterilized using ethylene oxide		Use-by date
	Non-pyrogenic		Date of manufacture
	Consult instructions for use		Manufacturer
	Do not use if package is damaged and consult instructions for use		Caution
	Do not re-sterilize		CE conformity marking
	Do not re-use		Federal (USA) law restricts this device to sale by or on the order of a licensed healthcare professional.
	Catalogue number		Swiss authorised representative
	Batch code		Unique device identifier
	Temperature limit		

Made in Ireland
Merck Millipore Ltd.
Tullagreen, Carrigtwohill,
Co. Cork, Ireland
T45 KD29
an affiliate of Merck KGaA,
Darmstadt, Germany

MDSS CH GmbH
Laurenzenvorstadt 61
5000 Aarau
Switzerland

Brazilian Registration Holder (BRH)
Merck S/A
Rua Torre Eiffel, nº 100
Lote 3ª, Gleba A
Parque Rincão
Cotia/SP

an affiliate of Merck KGaA,
Darmstadt, Germany

*Brazil Product Registration
Number 81167709004

Cathivex®-GV Filter Unit (Sterile)

SLGV0250S Qty: 50 pack

- 25 mm
- Sterile
- Single use only
- Non-pyrogenic

Intended Purpose

The Cathivex®-GV Filter Unit (Sterile) is a medical device intended for use as an in-line filter to sterilize intravenously administered aqueous solutions in direct patient care and pharmacy admixture applications. Device is intended for medical and/or laboratory professional use only.

Introduction

This filter removes particulates, microbial contamination, and air bubbles in applications where venting and low protein binding membranes are required or desired in sterile filtration of aqueous vitamins, drugs, and other clinical solutions for intravenous infusions. This sterile, single use filter is non-pyrogenic and non-toxic. It consists of a 0.22 µm membrane filter sealed in a polyvinyl chloride (PVC) housing with a female Luer-Lok™ inlet connector and a male Luer-Lok™ outlet connector. A 0.03 µm hydrophobic filter in the dome prevents air-locking by automatically venting air introduced upstream. Clinical benefits include the removal of microorganisms, particles, precipitates, and undissolved powders larger than the membrane's rated pore size from solutions for clinical use.

Chemical Compatibility

These filters are compatible with most aqueous solutions. Compatibility of solution should be tested before use. Visit [SigmaAldrich.com/MedMillex](https://www.sigmaaldrich.com/MedMillex) for more information.

WARNINGS

- Filter solutions with a temperature limit of 45 °C (113 °F).
- For manual use only; do not use on automated systems.
- Not suitable for filtering high viscosity solutions or blood.
- Do not filter solutions containing 5 mg or less of active drug materials unless binding studies have been performed.
- To ensure sterility, do not use this product if the package is damaged.
- Use aseptic technique to maintain sterility removing product from primary package.
- Do not use for direct neuraxial applications which may come in contact with cerebrospinal fluid.
- Do not use the same in-line filter to filter solutions in both directions.
- Sudden changes in pressure may damage filter.
- Sudden loss of pressure could indicate failure of the filter.
- Single use only; do not re-use or resterilize.
- Make sure to wet the GV filter membrane thoroughly before injecting the solution as the hydrophobic vent will not wet out. Improperly wetted filters can become airlocked.
- Change within 24 hours.
- In peripheral administration, the filter should be secured at the venipuncture site. In other intravenous or arterial applications (or in peripheral situations where securing at the venipuncture site is not possible), the filter should be positioned below the patient's midaxillary line to maintain a full liquid level in the filter housing.
- Discard appropriately after single use. See "Disposal" below.

Directions for Use

Use aseptic technique. Remove protective coverings as assembly progresses.

1. Connect primed administration set male adaptor to the filter female Luer-Lok™ connector (inlet) and remove assembly from package.
2. Open the administration set flow control clamp. Prime the filter with solution while tapping to dislodge air bubbles. Close flow control clamp.
3. If using a needle: attach it to the male Luer-Lok™ connector (outlet), clear air by pushing a few drops through to wet filter, and make venipuncture.
4. Regulate flow with administration set flow control clamp.

Specifications and Performance Characteristics

Materials

Filter membrane	Hydrophilic Durapore® polyvinylidene fluoride (PVDF), type GV, pore size 0.22 µm
Vent membrane	Hydrophobic polytetrafluoroethylene (PTFE), pore size 0.03 µm
Housing	Polyvinyl chloride (PVC)

Dimensions

Inlet to outlet	25 mm
Diameter	29 mm
Filtration area	4 cm ²

Temperature limit	45 °C (113 °F)
--------------------------	----------------

Pressure limit at 21 °C	3 bar (45 psi) inlet and differential
--------------------------------	---------------------------------------

Filtration volume	1-100 mL
--------------------------	----------

Hold-up volume	≤ 0.1 mL after air purge (typical volume when vent hole is covered to prevent air flow)
-----------------------	---

Sterilization method	Ethylene oxide gas
-----------------------------	--------------------

Connections	Female Luer-Lok™ inlet; male Luer-slip outlet
--------------------	---

Water flow rate at 21 °C	≥ 4 mL/min at pressure of 50 in. (127 cm) water
---------------------------------	---

Notice

We provide information and advice to our customers on application technologies and regulatory matters to the best of our knowledge and ability, but without obligation or liability. Existing laws and regulations are to be observed in all cases by our customers. This also applies in respect to any rights of third parties. Our information and advice do not relieve our customers of their own responsibility for checking the suitability of our products for the envisaged purpose.

Contact Information

For the location of the office nearest you, go to SigmaAldrich.com/offices.

Disposal

Follow precautions for disposal of items contaminated with potentially infectious or hazardous biological material according to all applicable international, federal, state, and local regulations.

Technical Assistance

Visit the tech service page on our web site at SigmaAldrich.com/techservice.

Any serious incident of this device should be reported to manufacturer and competent authority of country where user is established.

Standard Warranty

The applicable warranty for the products listed in this publication may be found at SigmaAldrich.com/terms.

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Vyrobené v Írsku

Merck, Millipore, Millex, Durapore and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

© 2023 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.

The life science business of Merck operates as MilliporeSigma in the U.S. and Canada.

MERCK