3050 Spruce Street, St. Louis, MO 63103 USA
Tel: (800) 521-8956 (314) 771-5765 Fax: (800) 325-5052 (314) 771-5757
email: techservice@sial.com sigma-aldrich.com

Product Information

Anti-α-Tubulin antibody, Mouse monoclonal clone DM1A, purified from hybridoma cell culture

Product Number T6199

Product Description

Anti- α -Tubulin antibody, Mouse monoclonal (mouse IgG1 isotype) is derived from the hybridoma DM1A produced by the fusion of mouse myeloma cells (NS1) and splenocytes from BALB/c mice immunized with purified chick brain tubulin. The isotype is determined by a double diffusion immunoassay using Mouse Monoclonal Antibody Isotyping Reagents, Product Number ISO2.

Anti- α -Tubulin antibody, Mouse monoclonal recognizes an epitope located at the C-terminal end of the α -tubulin isoform (amino acids 426-430)^{2, 3, 4} in a variety of organisms (e.g., human, bovine, mouse, and chicken). The product may be used in various immunochemical techniques including immunoblotting, ¹ solid-phase RIA, ¹ cell antibody microinjection, ¹ immunocytochemistry, ⁹ immunohistochemistry, ¹⁰ and immunoprecipitation. ⁸

Tubulin is the major building block of microtubules. This intracellular, cylindrical, filamentous structure is present in almost all eukaryotic cells. Microtubules function as structural and mobile elements in mitosis, intracellular transport, flagellar movement, and the cytoskeleton. Tubulin is a heterodimer that consists of α-tubulin and β-tubulin. Both subunits have a molecular weight of approx. 50 kDa and share considerable homology. In addition to α - and β -tubulin, several other tubulins have been identified, bringing the number of distinct tubuling classes to seven. Most of these tubulins have distinct subcellular localization and an emerging diverse set of functions.6 Out of the seven different tubulins four new members of the tubulin family were identified recently, which consist of δ , ξ , η , and ϵ -tubulin. η and ϵ -tubulins were discovered by database searches. Microtubular systems contain at least three α -tubulin isoforms. Two isoforms are coded by two α -tubulin genes, which are both transcribed and code for extremely similar proteins. The third isoform is generated by post-translational modification.⁵ At least three modifications of tubulin subunits have been described: the phosphorylation of β -tubulin from brain, the removal of the carboxyterminal tyrosine form α-tubulin in vertebrate tissues, and the acetylation of the amino group of lysine(s) in α -tubulin.

Monoclonal antibodies recognizing α -tubulin, together with monoclonal antibodies to other tubulin types (β , β -tubulin isotype I +II, β -tubulin isotype III, tyrosine tubulin, and the acetylated form of α -tubulin) provide a specific and useful tool in studying the intracellular distribution of tubulin and the static and dynamic aspects of cytoskeleton.

Reagent

Supplied as a solution in 0.01 M phosphate buffered saline, pH 7.4, containing 1% bovine serum albumin and 15 mM sodium azide.

Antibody concentration: Approx. 1 mg/ml

Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Safety Data Sheet for information regarding hazards and safe handling practices.

Storage/Stability

Store at 2-8 °C for up to one month.

For prolonged storage, freeze in working aliquots at -20 °C. Repeated freezing and thawing, or storage in frost-free freezers, is not recommended. If slight turbidity occurs upon prolonged storage, clarify the solution by centrifugation before use. Working dilutions should be discarded if not used within 12 hours.

Product Profile

Immunoblotting: a working antibody concentration of 0.5-1 μ g/ml is recommended using a total tissue extract from chicken gizzard.

Immunocytochemistry: a working antibody concentration of 0.5-1 μ g/ml is recommended using cultured chicken fibroblasts (CFB).

Note: In order to obtain the best results using various techniques and preparations, we recommend determining optimal working dilutions by titration.

References

- 1. Blose, S.H., et al., J. Cell Biol., 98, 847-858 (1984).
- Breitling, F., and Little, M., J. Molec. Biol., 189, 367-370 (1986).
- 3. Wolff, A., et al., Biol. Cell, 63, 319-326 (1988).
- Serrano, L., et al., Anal. Biochem., 159, 253-259 (1986).
- LeDizet, M., and Piperno, G., J. Cell Biol., 104, 13-22 (1986).
- Dutcher, S.K., et al., Curr. Opin. Cell Biol., 13, 49-54 (2001).

- 7. McKean, P.G., et al., *J. Cell Sci.*, **114**, 2723-2733 (2001).
- 8. Eiserich, J.P., et al., *Proc. Natl. Acad. Sci. USA*, **96**, 6365-6370, (1999).
- 9. Perez, R.G., et al., *J. Neurol.*, **17**, 9407-9414 (1997).
- 10. Machaca, K. and Haun, S., et al., *J. Biol. Chem.*, **49**, 38710-38715 (2000).

DS,PHC 11/15-1