

3050 Spruce Street Saint Louis, Missouri 63103 USA Telephone (800) 325-5832 (314) 771-5765 Fax (314) 286-7828 email: techserv@sial.com sigma-aldrich.com

ProductInformation

FIBROBLAST GROWTH FACTOR-ACIDIC (aFGF) From Bovine Pituitary Glands

Product No. F 5267

Product Description

Fibroblast Growth Factor-Acidic (aFGF) is a potent mitogenic agent for a wide variety of mesoderm-derived cells including Balb/c 3T3 fibroblasts, capillary and endocardial endothelial cells, myoblasts, vascular smooth muscle cells, mesothelial cells, gial, and astroglial cells, and adrenal cortex cells.^{1,2} Isolated from bovine brain by a method involving heparin affinity chromatography,^{3,4} aFGF is a 15.9 kDa protein⁵ having two potential binding domains for heparin.⁶ In cells not growth-inhibited by heparin (such as baby hamster kidney cells), heparin potentiates the actions of aFGF, but in cells inhibited by heparin (such as bovine brainderived capillary endothelial cells), no such potentiation is observed.⁷ The closely related protein Fibroblast Growth Factor-Basic (Product No. F 5392), purified from bovine pituitary glands, acts upon the same cellular receptors as aFGF but with differing specific activities, depending on the cell type.⁸ Both aFGF and bFGF may play important roles in vivo in cell proliferation and differentiation associated with embryogenesis, tissue regeneration, CNS development, wound healing, angiogenesis, and tumor progression.² Although bFGF is found in a variety of organs, aFGF has been found only in brain, hypothalamus and retina. Acidic FGF has several synonyms, including endothelial cell growth factor, heparin-binding growth factor (class I or alpha), retina-derived growth factor, and astroglial growth factor 1.9

The biological activity of Fibroblast Growth Factoracidic is determined in a cell proliferation assay using fetal bovine heart endothelial cells. The EC_{50} is defined as the effective concentration of growth factor that elicits a 50% increase in cell growth in a cell based bioassay.

Reagents

 α FGF is lyophilized from 25 mM sodium phosphate, 50 mM NaCl, pH 7.0 containing Bovine Serum Albumin, 100 μ g/vial as a carrier protein.

Reconstitution and Use

To prepare a stock solution, reconstitute with 1-5 ml sterile tissue culture media containing 0.1-1.0% BSA. This may be diluted immediately before use to the final working concentration of aFGF, generally 0.3-30 ng/ml. Additional filtration is not recommended and may result in product loss due to adsorption onto filter membrane.

Storage/Stability

Store vial at -20 °C frozen.

After reconstitution, store aliquots at 2-8 °C for two weeks or frozen. Prolonged storage of product or repeated freezing and thawing is not recommended.

References

- 1. Gospodarowicz, D. and Moran, J., Proc. Natl. Acad. Sci. U.S.A., **71**, 4648 (1974).
- 2. Gospodarowicz, D. et al., Endocrine Rev., **8**, 95 (1987).
- Gospodarowicz, D. and Moran, J., Proc. Natl. Acad. Sci. U.S.A., 71, 4584 (1974).

JWM/DAA 3/2004

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.