

Product Information

Glutathione-Agarose

Set of 3 pre-packed columns (2.5 mL each),
1:1 suspension in a 0.5 M NaCl + 20% ethanol solution

G3907

Product Description

Affinity chromatography with glutathione-agarose permits rapid, mild, non-denaturing and highly selective purification of glutathione-binding enzymes,¹⁻¹² such as glutathione-S-transferase, glutathione peroxidase, and glyoxalase I. This product consists of pre-packed columns where the resin has glutathione attached through its sulfur to epoxy-activated 4% cross-linked beaded agarose, to give a 12-atom spacer.

Several publications have cited use of this G3907 product in their protocols.¹³⁻¹⁴

Product Profile

- Three columns pre-packed with 2.5 mL each of a saline suspension of glutathione-agarose
- Ligand: 10 to 20 μ moles per mL swollen resin
- Binding capacity: \sim 5-10 mg glutathione S-transferase per mL resin

Storage/Stability

If stored properly at 2-8 °C, the product has a shelf life of two years. If exposed to heat and moisture, the beads may not swell properly or may not be active.

Precautions and Disclaimer

For R&D use only. Not for drug, household, or other uses. Please consult the Safety Data Sheet for information regarding hazards and safe handling practices.

Preparation Instructions

Reagents

Equilibration Buffer: Phosphate buffered saline (PBS), 10 mM phosphate buffer (pH 7.4) and 150 mM NaCl (such as Cat. No. P3813). Add appropriate protease inhibitors when preparing cell lysate.

Triton X-100: Cat. No. T9284.

PBS-T: PBS containing 1% Triton™ X-100.

Elution Buffer: 5 mM to 10 mM reduced glutathione (such as Cat. No. G4251) in 50 mM Tris-HCl, pH 8.0, prepared fresh. Final pH: \sim 7.5.

Cleansing Buffer 1: 0.1 M borate buffer, pH 8-9, with 0.5 M NaCl. Prepare using boric acid (such as Cat. No. B0394). Adjust pH with NaOH.

Cleansing Buffer 2: 0.1 M acetate buffer, pH 4, with 0.5 M NaCl. Prepare acetate buffer using sodium acetate (such as Cat. No. 236500). Adjust pH with acetic acid (such as Cat. No. 695092).

Storage Buffer: 1 M NaCl with 1 mM sodium azide as preservative.

Procedure

1. Equilibrate the resin with several column volumes of Equilibration Buffer. Do not let the column run dry at any time.
2. Prepare cell lysate in appropriate buffer.
 - Tris or phosphate buffers, pH 6.5-9.5, are typical lysis buffers which are compatible with glutathione affinity chromatography.
 - Salt concentrations of up to 1 M do not interfere with binding.
 - Protease inhibitors such as EDTA (Cat. No. E7889) or PMSF (Cat. No. P7626) are often included in the lysis buffer. Serine protease inhibitors included in the lysis buffer will not interfere with subsequent thrombin or factor Xa treatment, as these inhibitors are removed before the proteolysis step.
 - Binding of GST to glutathione-agarose is unaffected by 1% Triton™ X-100, 1% Tween®-20, 1% CTAB, 10 mM DTT, or 0.03% SDS.
3. Add Triton™ X-100 to a final concentration of 1% (v/v).

- Centrifuge 10 minutes at 10,000 × *g* at 4 °C, or filter through a 0.45 µm filter to clear cell lysate.

Note: Use only clarified supernatant. To prevent clogging the column, highly viscous samples containing chromosomal DNA or RNA should be sonicated or treated with nuclease to reduce the viscosity. Cellular debris and particulate matter must be removed by centrifugation or filtration.

- Load the clarified supernatant onto the column under gravity flow.

Note: Depending on the sample and the flow rate, not all the protein may bind. Multiple passes over the column, or closing the loaded column and incubating it on a rotator, may improve the binding efficiency.

- Wash resin four times with PBS-T at 4 °C.
- Elute GST from the resin with Elution Buffer (3 times, 1 mL each). For batch purification, mix each elution step gently for 2 minutes.
- Analyze each fraction by SDS-PAGE.
- Free glutathione may be removed from sample by dialysis against buffer of choice.

Regeneration

- Wash with approximately 5 resin volumes of 0.1 M borate buffer, pH 8-9, containing 0.5 M NaCl (Cleansing Buffer 1).
- Wash with at least 5 column volumes of deionized water.
- Wash with at least 5 column volumes of 0.1 M acetate buffer, pH 4, containing 0.5 M NaCl (Cleansing Buffer 2).
- Wash with approximately 5 resin volumes of deionized water. For long-term storage, store in 1 M NaCl containing a bacteriostatic agent.
- Equilibrate with Equilibration Buffer before use.

References

- Simons, P.C., and Vander Jagt, D.L., *Anal. Biochem.*, **82(2)**, 334-341 (1977).
- Simons, P.C., and Vander Jagt, D.L., *Methods Enzymol.*, **77**, 235-237 (1981).
- Scopes, R., in *Protein Purification, Principles & Practice*. Springer-Verlag (New York, NY), pp. 30, 34 (1982).
- Toribio, F. *et al.*, *J. Chromatogr. B. Biomed. Appl.*, **684(1-2)**, 77-97 (1996).
- Huang, Yi, "Localization of Alpha-2 Adrenergic Receptor Subtypes Using Subtype-Specific Antibodies". University of Arizona, Ph.D. dissertation, pp. 47, 54, 166, 167 (1995).

- Spencer Valero, Lilian Maritza, "Monoclonal Antibodies Binding to Malarial Merozoite Surface Protein 1 Protect *in vivo* Against *Plasmodium yoelii* Infection". The Open University, Ph.D. dissertation, pp. 24-26, 40-41, 60-64 (1997).
- Li, J., and Tsai, M.-D., *Biochemistry*, **41(12)**, 3977-3983 (2002).
- Li, J. *et al.*, *Biochemistry*, **43(14)**, 4394-4399 (2004).
- Li, J. *et al.*, *Biochemistry*, **44(40)**, 13246-13256 (2005).
- Guo, Yi, "Study of the Structure and Function Relationship of Oncoprotein Gankyrin". The Ohio State University, Ph.D. dissertation, pp. 25, 43, 77 (2009).
- Zhang, Wen Xiao (Wendy), "Mechanism of Tocopherol Transfer by human α-Tocopherol Transfer Protein (α-hTTP)". Brock University, Ph.D. dissertation, p. 78 (January 2010).
- Aubry, C. *et al.*, *Org. Biomol. Chem.*, **4(5)**, 787-801 (2006).
- Chan, Y.K., and Gack, M.U., *Nat. Immunol.*, **17(5)**, 523-530 (2016).

Notice

We provide information and advice to our customers on application technologies and regulatory matters to the best of our knowledge and ability, but without obligation or liability. Existing laws and regulations are to be observed in all cases by our customers. This also applies in respect to any rights of third parties. Our information and advice do not relieve our customers of their own responsibility for checking the suitability of our products for the envisaged purpose.

The information in this document is subject to change without notice and should not be construed as a commitment by the manufacturing or selling entity, or an affiliate. We assume no responsibility for any errors that may appear in this document.

Technical Assistance

Visit the tech service page at SigmaAldrich.com/techservice.

Standard Warranty

The applicable warranty for the products listed in this publication may be found at SigmaAldrich.com/terms.

Contact Information

For the location of the office nearest you, go to SigmaAldrich.com/offices

The life science business of Merck KGaA, Darmstadt, Germany operates as MilliporeSigma in the U.S. and Canada.

MilliporeSigma, Millipore, MILLIPLEX and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

© 2022 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.

G3907dat Rev 01/22 AC,GCY