

Quick Start Guide

GenElute™ HP Endotoxin-Free Plasmid Megaprep Kit

NA0600

Reagents to Prepare

- Add 1.25 mL of RNase A Solution to the Resuspension Solution. Store at 2-8 °C.
- Add 300 mL of 95-100% ethanol to Wash Solution 2. After each use, tightly cap the diluted Wash Solution to prevent ethanol evaporation.

Protocol

Vacuum format.

Harvest Bacteria

1. Pellet cells from 200-1000 mL overnight culture by spinning for 10 min at 5,000 x g. Discard supernatant.

Prepare Column

2. Place Column onto vacuum manifold.
3. Add 40 mL of Column Preparation Solution to column, apply vacuum and allow liquid to pass through.

Resuspend and Lyse Bacteria

4. Resuspend cells in 40 mL Resuspension Solution. Pipet to mix.
5. Add 40 mL Lysis Solution. Gently invert 6-8 times to mix. Do not vortex. Allow to clear for 3-5 min.

Prepare Cleared Lysate

6. Add 40 mL Neutralization Solution (chilled prior to use) to the lysed cells and gently invert 6-8 times to mix.
7. Attach the Lysate Filter to a 45 mm neck bottle (≥250 mL).
8. Empty entire bag of Lysate Clearing Agent into lysate and gently invert 6-8 times to mix. Incubate 5 min at room temperature.
9. Add lysate to assembled Lysate Filter and apply vacuum to clear cell debris.
10. Add 50 mL Binding Solution to the bottle containing the cleared lysate and swirl to mix.

Bind DNA to Column

11. Add cleared lysate to prepared column, apply vacuum and allow all lysate to pass through.

Wash to Remove Contaminants

12. Add 40 mL of Wash Solution 1 to column and allow it to pass through.
13. Add 60 mL Wash Solution 2 to column and allow it to pass through.
14. Leave vacuum on for 20-30 min to dry the column.

Elute Purified DNA

15. Attach a VacCap to a 50 mL collection tube and attach the column to the VacCap.
16. For maximum plasmid yield, add 20 mL of Endotoxin-Free Water to column and apply vacuum for 2-3 min.
17. For maximum plasmid concentration, add 5 mL of Endotoxin-Free Water and apply vacuum for 2 mins. Add another 5 mL of Endotoxin-Free Water and apply vacuum for another 5 min.

Notice

We provide information and advice to our customers on application technologies and regulatory matters to the best of our knowledge and ability, but without obligation or liability. Existing laws and regulations are to be observed in all cases by our customers. This also applies in respect to any rights of third parties. Our information and advice do not relieve our customers of their own responsibility for checking the suitability of our products for the envisaged purpose.

The information in this document is subject to change without notice and should not be construed as a commitment by the manufacturing or selling entity, or an affiliate. We assume no responsibility for any errors that may appear in this document.

Contact Information

For the location of the office nearest you, go to [SigmaAldrich.com/Offices](https://www.sigmaaldrich.com/Offices).

Technical Assistance

Visit the tech service page on our web site at [SigmaAldrich.com/TechService](https://www.sigmaaldrich.com/TechService).

Standard Warranty

The applicable warranty for the products listed in this publication may be found at [SigmaAldrich.com/Terms](https://www.sigmaaldrich.com/Terms).