

ProductInformation

PAP
(Peroxidase-Anti-Peroxidase Soluble Complex)
Antibody developed in Mouse

Product No. **P 3039**

Product Description

Antiserum is developed in mouse using Horseradish Peroxidase as the immunogen. Soluble antigen-antibody complexes of horseradish peroxidase anti-peroxidase (PAP) are prepared by a modification of the method reported by Sternberger, et al.¹

Each milliliter of product contains 3.5-9.0 mg PAP soluble complex as determined by the method of Lowry (prior to the addition of BSA).

Reagents

The product is provided as a solution in 8 mM sodium acetate, 15 mM ammonium acetate buffered saline, pH 7.2, containing 0.5% BSA, with 0.01% thimerosal as a preservative.

Product Profile

To ensure functional product reactivity, each lot of PAP preparation is assayed utilizing standard histological staining techniques.² A minimum working dilution of 1:200 was determined using indirect immunoperoxidase techniques on chicken fibroblasts.

By immunoelectrophoresis (IEP), a single arc of precipitation is observed when the product (prior to the addition of BSA) is electrophoresed and then allowed to diffuse against anti-mouse whole serum or anti-mouse IgG. No reaction is observed against free horseradish peroxidase.

Storage/Stability

For continuous use, store at 2-8 °C for up to one month. For extended storage, the solution may be frozen in working aliquots. Repeated freezing and thawing is **not** recommended. Storage in "frost-free" freezers is **not** recommended. If slight turbidity occurs upon prolonged storage, clarify the solution by centrifugation before use.

References

1. Sternberger, L.A., et al., J. Histochem. Cyto-chem., **18**, 315 (1970).
2. Sternberger, L.A., Immunocytochemistry (John Wiley & Sons Inc., New York), p. 122 (1979).

PCS 2/99

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications.

Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply.

Please see reverse side of the invoice or packing slip.