


MultiScreen[®]_{HTS} 384-well Filter Plate with Glass Fiber Filter

protocol note



High throughput method for in-plate receptor-ligand binding assays and radioisotope counting

Introduction

Heterogeneous filter binding assays have been used for the purpose of characterizing receptor-ligand binding. To date, most assays have been performed on 96-well platforms. The new MultiScreen_{HTS} 384-well filter binding format increases the density of samples that can be screened in a single plate, making it possible to significantly reduce the amount of reagents for each assay. MultiScreen_{HTS} filter plates are also designed for binding reactions and analyses to be performed directly in the plate.

General Considerations

This protocol note provides an overview of the basic receptor-ligand binding reaction which can be performed in MultiScreen_{HTS} 384-well filter plates and provides guidance in scaling this assay from the 96-well format. Individual binding assay parameters (such as receptor concentration, ligand concentration, buffer components, etc.) can differ greatly from one system to another and will necessarily need to be optimized for each system.

Reactions

Displacement binding is a method by which a constant amount of radiolabeled ligand is competitively displaced by a non-labeled competitor ligand. This is done with a titration of the non-labeled ligand to estimate the relative binding affinity for the receptor (IC₅₀, K_i) or is used to screen large libraries of compounds for their ability to displace the radioligand.

Materials

1. MultiScreen_{HTS} 384-well Filter Plate with glass fiber filter (MZFB NOW or MZFC NOW)
2. MultiScreen_{HTS} Vacuum Manifold (MSVM HTS 00)
3. Robotic liquid handling system for automated assays (strongly preferred) or multichannel pipette for manual assays
4. Liquid scintillation counting equipped with plate reader for 384-well plates (i.e. MicroBeta Trilux[®] Plate Reader)
5. Liquid scintillant cocktail (i.e. Packard Supermix[™] or PerkinElmer Optifluor[™])
6. Oven capable of maintaining 25 °C to 50 °C
7. Assay reagents:
 - a) Reaction buffer (defined by user) including buffer (i.e. Tris), salts (i.e. NaCl), essential co-factors
 - b) Receptor preparation. May be in the form of purified recombinant receptor, partially purified cellular membranes or tissue homogenates
 - c) Labeled ligand. ³H or ¹²⁵I labeled with high specific activity (typically 50Ci/mmol or higher).
 - d) Unlabeled ligand and/or test compounds for performing displacement binding assays

Displacement Binding

For each well in the MultiScreen_{HTS} 384-well filter plate, 100 µL total; reaction volume:

10 µL	Radioligand dilution at approx. K _d concentration
10 µL	Compound(s) being tested at single concentration or serial dilutions
10 µL	Diluted receptor preparation
70 µL	Incubation buffer
<hr/>	
= 100 µL	total well

Note: the suggested volume is 50 – 100 µL for binding reactions in the 384-well format. When converting from 100 – 200 µL assays performed in 96-well format to 384-well format, keep in mind that adding the same mass of receptor in a smaller volume will result in increasing its effective concentration and can lead to radioligand depletion (i.e. binding affinity appears lower than actual since a higher concentration of ligand is required). This is of particular concern if the final receptor concentration is significant compared to the binding affinities of the radioligands or the compounds being measured. Therefore, the amount of receptor should be reduced according to volume to prevent ligand depletion as long as good radiometric signal is maintained. These parameters should be determined by the user.

Protocol

General Binding Protocol for In-Plate Receptor-Ligand Binding*:

1. If required, add 100 μ L 0.5% PEI to each well of the MultiScreen_{HTS} 384 (Cat # MZFBNOW or MZFCNOW) plates to be tested and allow to incubate at 4 °C, 2 – 24 hours. (Do not allow PEI solution to dry. Use a lid to prevent evaporation.)
2. Wash treated wells 2 times with cold (4 – 8 °C) incubation buffer.
3. As required, thaw or reconstitute receptor preparation and radiolabeled ligand on ice. Prepare or acquire stock solutions of compounds to be screened.
 - a) Dilute radiolabeled ligand in incubation buffer.
 - b) Dilute compounds to be screened in incubation buffer.
 - c) As required, dilute receptor preparation in incubation buffer.
4. In each well, add compounds first then radiolabeled ligand, mix then add diluted receptor preparation (final volume of reaction should be 50 – 100 μ L).
5. Incubate 1 – 2 hours at room temperature.
6. Vacuum filter reaction volume using the vacuum manifold (8 – 15" Hg).
7. Wash each test well 5 – 10 times by adding and then vacuum filtering 100 μ L of incubation buffer.
8. After the last wash, apply vacuum until the glass fiber filter is almost dry.
9. Allow plates to dry at 50 °C for at least 1 hour. Alternatively, dry plates overnight at 37 °C.
10. When using tape, place clear sealing tape over bottom of plate.
11. Add 5 – 10 μ L of scintillation cocktail (e.g., Optifluor) to each well.
12. Place clear sealing tape over top of plate[†].
13. Let the plate stand for at least 2 hours (best signal may be obtained with overnight for incubation).
14. Assay plates for radioactivity on a beta scintillation counter (e.g., MicroBeta Trilux Plate Reader [coincidence mode recommended] or equivalent: see Radioisotope Counting).

* If out-of-plate incubation is desired, transfer incubations from separate solid-bottom plate to prewashed filter plate wells (steps 1 and 2). Quickly continue from step 6.

† If using a MicroBeta Trilux 32 plate stacker, use every other slot for a total of 16 plates. If 32 plates are desired, tape bottom of cassette instead of bottom of filter plate.

Radioisotope Counting

Radioactivity (scintillation counting) can be measured directly in the Multiscreen_{HTS} 384-well filter plate. The design allows the plate to be analyzed from the top (top counting) or from both sides simultaneously (coincidence counting). For best signal to noise ratios, it is recommended to count in the coincidence mode. Low energy isotopes such as ¹⁴C and ³H require the use of scintillation fluid. For best results, it is recommended to use 5 – 10 µL of scintillation cocktail per well when counting a 384-well filter plate. For optimal results and consistency, the scintillation cocktail should be allowed to mix with the sample for at least an hour before counting. In some instances, longer incubations (up to 24 hours) may result in even higher counting efficiencies. The type and amount of scintillation cocktail, and how much

time should elapse prior to counting to produce the best and most consistent results, should be determined by the user.

Typical counting efficiencies for ³H labeled compounds are listed in Table 1.

counting efficiency value before fitting data. The same amount of bound radioligand in the 96-well filter plate will give more CPM than an equivalent amount in the 384-well filter plate. However, the determined specific

$$\text{Counting Efficiency} = \frac{\text{CPM measured in counting instrument}}{\text{Disintegration per minute (DPM)}}$$

Counting efficiency values for 96-well filter plates are typically greater than for 384-well plates. This is especially true when the underdrain of the 96-well filter plate is removed. The design of the MultiScreen_{HTS} 384-well filter plate allows for direct coincidence counting through the underdrain without any need to remove it.

It is important when converting from CPM to moles of bound ligand per mass of protein in the assay (i.e. fmol/mg) to use the appropriate

bound amount (i.e. fmol/mg) will be the same if the appropriate counting efficiency value is used.

The signal to noise, however, is typically the same or better in the 384-well format depending on counting mode, scintillant, etc. (Figure 1).

Important Note: Different assays have different counting criteria and levels of signal. For the most accurate results, counting efficiency should be determined by the operator by adding a known DPM of isotope and counting.

Coincidence Counting Efficiencies

³ H Millipore Plate Type	Cat#	# Wells	Packard Optifluor	Wallac Supermix
			% ³ H Efficiency at 50K DPM	% ³ H Efficiency at 50K DPM
MultiScreen _{HTS} 384 FB (MZFB NOW)		384	8.7	13.7
MultiScreen _{HTS} FB (MSFB NOB)		96	29.4/51.5*	38.9/67.1*
MultiScreen _{HTS} 384 FC (MZFC NOW)		384	7.1	6.7
MultiScreen _{HTS} FC (MSFC NOB)		96	36.6/60.3*	45.8/72.5*
¹²⁵ I Plate Type	Cat#	# Wells	Packard Optifluor	Wallac Supermix
			% ¹²⁵ I Efficiency at 50K DPM	% ¹²⁵ I Efficiency at 50K DPM
MultiScreen _{HTS} 384 FB (MZFB NOW)		384	20.4	22.0
MultiScreen _{HTS} FB (MSFB NOB)		96	58.1	80.0
MultiScreen _{HTS} 384 FC (MZFC NOW)		384	18.4	18.7
MultiScreen _{HTS} FC (MSFC NOB)		96	63.9	82.2

Table 1. Typical radiometric counting efficiencies values for 50,000 DPM ³H and ¹²⁵I labeled compounds as determined in the Wallac Microbeta Trilux Counting Instrument with 10 µL (384-well) or 50 µL (96-well) Packard Optifluor or Wallac Supermix (scintillation cocktail). Values presented are the average of 20 efficiency determinations.

* ³H values for 96-well filter plate are listed with underdrain on/with underdrain removed.




Figure 1. The design of the MultiScreen_{HTS} 384-well filter plate allows for radioactivity detection equivalent to the 96-well MultiScreen_{HTS} filter plate. Signal to noise was determined for the MultiScreen_{HTS} 384-well and the MultiScreen_{HTS} 96-well filter plates with FB glass fiber with ³H labeled ligand. Counting was performed in the Wallac MicroBeta TriLux Liquid scintillation counter in default coincidence counting mode with Optifluor (PerkinElmer) liquid scintillant cocktail. CPM were determined as a function of initial DPM loading. The CPM values were divided by the average CPM of the remaining wells that did not contain radioactive material.

Data Analysis

Displacement Binding

Radioligand-binding inhibition is determined with a constant radioligand concentration usually at or around the K_d value for the receptor. Serial dilutions of an unlabelled competitor ligand were normalized to a control binding experiment without unlabelled ligand (% Control) (see Figure 2). Relative affinity values (IC₅₀) are determined by fitting binding inhibition values by non-linear regression using Prizm[®] data software (www.graphpad.com). For screening of compounds, a single concentration is used as a first screen to identify those compounds which prevent radioligand binding, followed up by the quantitative displacement binding described above to determine IC₅₀ values.




Figure 2. Radiolabeled ligand-binding inhibition is determined with a constant radiolabeled ligand concentration and serial dilutions of unlabelled competitor ligand. Data are normalized to the specific binding of radiolabeled ligand observed in the absence of competitor (% Control).

Summary

- The use of filter plates for in-plate receptor-ligand binding assays eliminates the need for a separate incubation plate thereby decreasing the number of manipulations and reducing radioactive solid waste.
- Optimizing 384-well plates for coincidence counting gives highly quantitative data that can be collected at very low levels of radioactivity. The data remains linear over several orders of magnitude of signal.
- Receptor-ligand binding assays in 384-well MultiScreen_{HTS} filter plates demonstrate the same robustness and reliability as 96-well assays but use half the reagents in half the assay volume.

To Place an Order or Receive Technical Assistance

For additional information call your nearest Millipore office.

In the U.S. and Canada, call toll-free **1-800-MILLIPORE (1-800-645-5476)**

In the U.S., Canada and Puerto Rico, fax orders to **1-800-MILLIFX (1-800-645-5439)**

Internet: www.millipore.com

Tech Service:

www.millipore.com/techservice

pure e-commerce

Now you can buy Millipore products online @

www.millipore.com/purecommerce

Millipore and MultiScreen are registered trademarks of Millipore Corporation.

Optifluor and MicroBeta are trademarks of PerkinElmer.

Supernix is a registered trademark of Wallac Oy.

Prizm is a registered trademark of Graphpad.

Lit. No. PC1051EN00 Rev. - 12/04 04-098

Printed in the U.S.A.

© 2004 Millipore Corporation, Billerica, MA U.S.A. All rights reserved.

Related Literature

AN086EN00: MultiScreen_{HTS} 384-well Receptor-ligand Binding Application Note

PF1150EN00: MultiScreen_{HTS} Glass Fiber Filter Plates Data Sheet

Millipore Worldwide

AUSTRALIA

Tel. 1 800 222 111
or (02) 9888 8999
Fax (02) 9878 0788

AUSTRIA

Tel. 0820 87 44 64
Fax 0820 87 44 44

BALTIC COUNTRIES

Tel. +358 203 05 645
Fax +358 203 05 644

BELGIUM

Tel. 070 225 645
Fax 070 225 644

BRAZIL

Tel. (011) 5548-7011
Fax (011) 5548-7923

CANADA

Tel. (800) 645-5476
Fax (800) 645-5439

CHINA, PEOPLE'S REPUBLIC OF

Beijing:
Tel. (86-10) 8519 1250
(86-10) 6518 1058
Fax (86-10) 8519 1255
Guangzhou:
Tel. (86-20) 8752 0187
(86-20) 8752 0173
Fax (86-20) 8752 0172
Hong Kong:
Tel. (852) 2803 9111
Fax (852) 2513 0313
Shanghai:
Tel. (86-21) 5306 9100
Fax (86-21) 5306 0838

CZECH REPUBLIC

Tel. 02 2051 3841
Fax 02 2051 4298

DENMARK

Tel. 7010 5645
Fax 7010 5644

EASTERN EUROPE, C.I.S., AFRICA, MIDDLE EAST AND GULF

Tel. +33 3 88 38 9536
Fax +33 3 88 38 9539

FINLAND

Tel. 0203 05 645
Fax 0203 05 644

FRANCE

Tel. 0825 045 645
Fax 0825 045 644

GERMANY

Tel. 01 805 045 645
Fax 01 805 045 644

HUNGARY

Tel. 01 209 04 33
Fax 01 209 02 95

INDIA

Tel. 080 839 46 57
Fax 080 839 63 45

IRELAND

Tel. +44 1923 816375
Fax +44 1923 818297

ITALY

Tel. 848 845 645
Fax 848 845 644

JAPAN

Tel. (03) 5442-9711
Fax (03) 5442-9736

KOREA

Tel. (822) 551-0230
Fax (822) 551-0228

LUXEMBOURG

Tel. +32 2726 88 40
Fax +32 2726 98 84

MALAYSIA

Tel. 03-7957-1322
Fax 03-7957-1711

MEXICO

Tel. (55) 5576 9688
Fax (55) 5576 8706

THE NETHERLANDS

Tel. 0900 7 645 645
Fax 0900 7 645 644

NORWAY

Tel. 810 62 645
Fax 810 62 644

POLAND

Tel. 22 669 12 25
Fax 22 663 70 33

PORTUGAL

Tel. +34 91 728 39 60
Fax +34 91 729 29 09

PUERTO RICO

Tel. (787) 273-8495
Fax (787) 747-6553

SINGAPORE

Tel. 6842-1822
Fax 6842-4988

SPAIN

Tel. 901 516 645
Fax 902 011 644

SWEDEN

Tel. 0771 200 645
Fax 0771 200 644

SWITZERLAND

Tel. 0848 645 645
Fax 0848 645 644

TAIWAN

Tel. 886-2-2792-9333
Fax 886-2-2792-6555

U.K.

Tel. 0870 900 46 45
Fax 0870 900 46 44

U.S.A.

Tel. (800) 645-5476
Fax (800) 645-5439

IN ALL OTHER COUNTRIES

Millipore International
Tel. +1 781 533 8622
Fax +1 781 533 8630

MILLIPORE