

3050 Spruce Street Saint Louis, Missouri 63103 USA Telephone 800-325-5832 • (314) 771-5765 Fax (314) 286-7828 email: techserv@sial.com sigma-aldrich.com

ProductInformation

Lectin-Sepharose[®] 6MB from *Triticum vulgaris*

Product Number L 6257 Storage Temperature 2-8 °C

Product Description

At low pH (below pH 3), this lectin is a monomer (17 kDa by sedimentation velocity). However, it is a dimer (35 kDa by sedimentation velocity) at neutral to slightly acidic pH.^{1,2} By SDS-PAGE analysis, the monomers migrate as 18 kDa proteins.³

The absorption maximum (λ_{max}) for the native dimer is 272 nm with a molar extinction coefficient (E^M) of 1.09 x 10⁵. The pl varies by lectin isoform (isolectins I, IIa, III - pl = 8.7 +/- 0.3 and isolectin IIb - pl = 7.7 +/- 0.3).⁴

Lectins are proteins or glycoproteins of non-immune origin that agglutinate cells and/or precipitate complex carbohydrates. Lectins are capable of binding glycoproteins even in presence of various detergents.⁵ The agglutination activity of these highly specific carbohydrate-binding molecules is usually inhibited by a simple monosaccharide, but for some lectins, di, tri, and even polysaccharides are required.

Lectins are isolated from a wide variety of natural sources, including seeds, plant roots and bark, fungi, bacteria, seaweed and sponges, mollusks, fish eggs, body fluids of invertebrates and lower vertebrates, and from mammalian cell membranes. The precise physiological role of lectins in nature is still unknown, but they have proved to be very valuable in a wide variety of applications *in vitro*, including:

- 1. blood grouping and erythrocyte polyagglutination studies.
- 2. mitogenic stimulation of lymphocytes.
- 3. lymphocyte subpopulation studies.
- 4. fractionation of cells and other particles.
- 5. histochemical studies of normal and pathological conditions.

Sigma offers a range of lectins suitable for the above applications. Most Sigma lectins are highly purified by affinity chromatography, but some are offered as purified or partially purified lectins, suitable for specific applications. Many of the lectins are available conjugated to (conjugation does not alter the specificity of the lectin):

- 1. fluorochromes (for detection by fluorimetry).
- 2. enzymes (for enzyme-linked assays).
- 3. insoluble matrices (for use as affinity media).

Please refer to the table for general information on the most common lectins.

The inhibition of agglutination activity by di-N-acetylglucosamine $(GlcNAc)_2$ on this wheat germ lectin is reported to be aproximately 600 times greater than that of N-acetylglucosamine (GlcNAc). Tri-N-acetylglucosamine (GlcNAc)₃ is reported to be about 3000 times more inhibitory than GlcNAc.⁶

This product is immobilized on Sepharose[®] 6MB. The use of this immobilized lectin in purifying membrane glycoproteins has been published.⁷ This product has also been used to bind the oocyte receptors responsible for mediating the effects of insulin and IGF-1. The solubilized *Xenopus* oocyte receptor was applied to the lectin, and elution was performed with N-acetylglucosamine.⁸

This resin should be regenerated using 0.5 M NaCl containing Mg^{2+} , Mn^{2+} , Ca^{2+} , and Zn^{2+} (1 mM each). The resin should be incubated for 30 minutes, washed with fresh regeneration solution, and then be re-equilibrated with running buffer. If the resin is to be stored, the solution should contain all of the above metal ions plus a bacteriostat.

Precautions and Disclaimer

For Laboratory Use Only. Not for drug, household or other uses.

Preparation Instructions

This Sepharose conjugate is a suspension is in 0.9% NaCl and 20% ethanol. It should be centrifuged for 30 seconds at 1,000 x *g* to pellet. The supernatant should then be discarded and replaced by binding buffer dictated by the experiment.

Storage/Stability Aggregation is thought to occur in the presence of high concentrations of 2-mercaptoethanol.

Lectin			Mitogeni		
	MW (kDa)	Subunits	Blood Group	Sugar	Activity
Abrus precatorius			-		+
Agglutinin	134	4		gal	
Abrin A (toxin)	60	2		gal	
Abrin B (toxin)	63.8	2(αβ)		gal	
Agarius bisporus	58.5	-	-	β-gal(1→3)galNAc	
Anguilla anguilla	40	2	Н	α-L-Fuc	
Arachis hypogaea	120	4	Т	β-gal(1→3)galNAc	
Artocarpus integrifolia	42	4	Т	α-gal→OMe	+
Bandeiraea simplicifolia					
BS-I	114	4	А, В	α -gal, α -galNAc	
BS-I-A ₄	114	4	А	α-galNAc	
BS-I-B ₄	114	4	В	α-gal	
BS-II	113	4	acq, B, Tk, T	glcNAc	
Bauhinia purpurea	195	4	_	β-gal(1→3)galNAc	+
Caragana arborescens	60; 120 ^a	2/4	_	galNAc	
Cicer arietinum	44	2	-	fetuin	
Codium fragile	60	4	-	galNAc	
Concanavalin A	102	4	-	α-man, α-glc	+
Succinyl-Concanavalin A	51	2	-	α-man, α-glc	+ ^b
Cytisus scoparius	_	-	-	galNAc, gal	
Datura stramonium	86	2(αβ)	_	(glcNAc) ₂	
Dolichos biflorus	140	4	A ₁	α-galNAc	
Erythrina corallodendron	60	2	_	β-gal(1→4)glcNAc	+
Erythrina cristagalli	56.8	2(αβ)	_	β-gal(1→4)glcNAc	
Euonymus europaeus	166	4(αβ)	В, Н	α-gal(1→3)gal	+
Galanthus nivalis	52	4	(h)	non-reduc. α-man	
Glycine max	110	4	_	galNAc	+ ^c
lelix aspersa	79	_	А	galNAc	
lelix pomatia	79	6	A	galNAc	
athyrus odoratus	40-43	4(αβ)	_	α-man	+
ens culinaris	49	2	_	α-man	+
imulus polyphemus	400	_ 18	_	NeuNAc	
Bacterial agglutinin	_	_	_	galNAc, glcNAc	
ycopersicon esculentum	71	_	_	$(g cNAc)_3$	
Jaackia amurensis	130	2(αβ)	0	sialic acid	+
laclura pomifera	40-43	2(αβ) 2(αβ)	_	α -gal, α -galNAc	
/omordica charantia	115-129	2(αβ) 4(αβ)	_	gal, galNAc	
Nomoralca charantia Naja mocambique mocambiqu			_	yai, yaintAu	
		—	-	-	
Naja naja kaouthia	-	2	_ (b)	- « D mon	
Varcissus pseudonarcissus	26 T		(h)	α-D-man	

			Mitogenic		
Lectin	MW (kDa)	Subunits	Blood Group	Sugar	Activity
Perseau americana	-	-	-	-	
Phaseolus coccineus	112	4	-	-	
Phaseolus limensis	247(II)	8	А	galNAc	+
	124(III)	4			
Phaseolus vulgaris					
PHA-E	128	4	-	oligosaccharide	+
PHA-L	128	4	_	oligosaccharide	+
PHA-P				-	
PHA-M					
Phytolacca americana	32	_	-	(glcNAc)₃	+
Pisum sativum	49	4(αβ)	-	α-man	+
Pseudomonas aeruginosa PA-I	13-13.7	_	-	gal	+ ^c
Psophocarpus tetragonolobus	35	1	_	galNAc, gal	
Ptilota plumosa	65; 170	_	В	α-gal	
Ricinus communis				C	
Toxin, RCA ₆₀	60	2	-	galNAc, β-gal	
Toxin, RCA ₁₂₀	120	4	-	β-gal	
Sambucus nigra	140	4(αβ)	-	αNeuNAC(2→6)gal	+ ^c
				galNAc	
Solanum tuberosum	50; 100 ^a 1, 2		_	(glcNAc)₃	
Sophora japonica	133	4	А, В	β-galNAc	
Tetragonolobus purpureas	120(A)	4	Н	α-L-fuc	
	58(BA)	2	Н	α-L-fuc	
	117(C)	4	Н	α-L-fuc	
Triticum vulgaris	36	2	_	(glcNAc) ₂ , NeuNAc	+
Ulex europaeus					
UEAI	68	_	Н	α-L-fuc	
UEA II	68	_	_	(glcNAc) ₂	
Vicia faba	50	4(αβ)	-	man, glc	+
Vicia sativa	40	4(αβ)	-	glc, man	+
Vicia villosa	139	4	$A_{1+}T_n$	galNAc	
A ₄	134	4	A ₁	galNAc	
B ₄	143	4	T _n	galNAc	
Vigna radiata	160	4	-	α-gal	
Viscum album	115	4(αβ)	-	β-gal	
Wisteria floribunda	68	2	-	galNAc	

^a Concentration-dependent molecular weight

^b Non-agglutinating and mitogenic
^c Mitogenic for neuraminidase-treated lymphocytes

References

- 1. Nagata, Y., and Burger, M. M., Wheat germ agglutinin. Molecular characteristics and specificity for sugar binding. J. Biol. Chem., 249, 3116-3122 (1974).
- Harris, E. L., and Angal, S., Protein Purification 2. Methods: A Practical Approach, Oxford University Press (New York, NY: 1989), p. 270.
- Rice, R. H., and Etzler, M.E., Subunit structure of wheat 3. germ agglutinin. Biochem Biophys Res Commun., 59, 414 (1974).
- Kronis, K. A., and Carver, J. P., Specificity of isolectins 4. of wheat germ agglutinin for sialyloligosaccharides: a

360-MHz proton nuclear magnetic resonance binding study. Biochemistry, 21, 3050-3057 (1982).

5. Rueben, L., et al., Activities of lectins and their immobilized derivatives in detergent solutions. Implications on the use of lectin affinity chromatography for the purification of membrane glycoproteins. Biochemistry, 16, 1787-1794 (1977).

- Allen, A. K., et al., The purification, composition and specificity of wheat-germ agglutinin. Biochem. J., 131, 155-162 (1973).
- Lotan, R., et al., Activities of lectins and their immobilized derivatives in detergent solutions. Implications on the use of lectin affinity chromatography for the purification of membrane glycoproteins. Biochem., 16, 1787-1794 (1977).
- Janicott, M., et al., The insulin-like growth factor 1 (IGF-1) receptor is responsible for mediating the effects of insulin, IGF-1, and IGF-2 in Xenopus laevis oocytes. J. Biol. Chem., 266, 9382 (1991).

Sepharose is a registered trademark of Pharmacia Biotech.

IRB/MWM/JRC/NSB/SAG 4/03

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.